
DOC022.97.80057

# H-Series Meters

10/2013, Edition 2

English..... 3  
Français..... 20

## Table of contents

[Additional information](#) on page 3

[Specifications](#) on page 3

[General information](#) on page 5

[Installation](#) on page 9

[User interface and navigation](#) on page 11

[Startup](#) on page 13

[Operation](#) on page 14

[Maintenance](#) on page 15

[Troubleshooting](#) on page 17

## Additional information

Additional information is available on the manufacturer's website.

## Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	Handheld meters: 9 x 20 x 5 cm (3.5 x 8 x 2 in.) Benchtop meters: 20 x 13 x 8 cm (5 x 8 x 3 in.)
Weight	Handheld meters: 1300 g (3.0 lb.) Benchtop meters: 900 g (2.0 lb.)
Battery enclosure	Water resistant
Battery requirements	4-ANSI 15 A or IEC-LR6 (AA Alkaline)
Power consumption	Backlight on and Bluetooth™ active: 1 W Backlight on and Bluetooth™ inactive: 50 mW
Power source	Internal power source: 4 AA alkaline or rechargeable nickel metal hydride (NiMH) batteries; battery life: > 200 hours External power source: 100 to 240 VAC, 50/60 Hz input; 4.5 to 7.0 VDC; 100 mA (benchtop meters only)
Storage temperature	-20 to +40 °C (4 to 140 °F)
Operating temperature	5 to 40 °C (41 to 104 °F)
Operating humidity	Relative humidity: 50% maximum at 25 °C (77 °F), non-condensing
Input connectors	Handheld meters: 8-pin ISFET, BNC with phono jack temperature, 12-pin conductivity (H170G only) Benchtop meters: 8-pin ISFET BNC with phono jack temperature, 12-pin conductivity (H270G only), 2 mm reference, USB and external AC
USB adapter	Peripheral
Data memory (internal)	Up to 999 measurement results at user selectable intervals from 1 to 1999 seconds
Data storage	Automatic in store mode; user enables data logging mode; data is user stored, recalled and deleted
Data export	USB connection to PC; transfer the data log or as data is read
Languages	English
Temperature correction	Off, automatic and manual (parameter dependent)
Measurement	Continuous measurement
Protection rating	Handheld meters: IP67 Benchtop meters: IP42
Certifications	CE

Specification	Details
<b>pH</b>	
Calibration	Up to five points: 1.68, 4.01, 6.86, 7.00, 9.18, 10.01, 12.45
Accuracy	±0.01 pH
Resolution	0.01 pH
Range	-2.00 to 19.99
<b>mV</b>	
Calibration	None
Accuracy	±1 mV
Resolution	Autoranging, 0.1 and 1
Range	Autoranging, ±199.9 mV to ±1999 mV
<b>Temperature</b>	
Calibration	None
Accuracy	±0.5 °C
Resolution	0.1 °C (0.1 °F)
Range	-5 to 105 °C (23 to 221 °F)
<b>ISE</b>	
Calibration	Up to five points
Accuracy	Probe dependent
Resolution	0.1 ppm–0.1 ppt
Range	Autoranging, -0.0 ppm to 1999 ppt
<b>Conductivity</b>	
Calibration	Up to five points
Accuracy	±1% full scale or ±1 digit
Resolution	0.01 µS, 0.1 µS, 1 µS, 0.01 mS, 0.1 mS
Range	Autoranging: 0.00 to 19.99 µS, 20.0 to 199.9 µS, 200 to 1999 µS, 2.00 to 19.99 mS, 20.0 to 199.9 mS
<b>TDS</b>	
Calibration	Up to five points
Accuracy	±1% full scale or ±1 digit
Resolution	0.01 ppm, 0.1 ppm, 1 ppm, 0.01 ppt, 0.1 ppt, 1 ppt, 0.1 mg/L, 1 mg/L, 0.01 gal/L, 0.1 gal/L
Range	Autoranging, ppm: 0.00 to 9.99 ppm, 10.0 to 99.9 pm, 100 to 999 ppm, 1.00 to 9.99 ppt, 10.0 to 99.9 ppt, 100 to 200 ppt mg/L: 0.00 to 199.9 mg/L, 200 to 1999 mg/L, 2.00 to 19.99 gal/L, 20 to 50 gal/L
<b>Salinity</b>	
Calibration	None (derived from conductivity)
Accuracy	±0.1 ppt (-2 to +35 °C or 28.4 to 95 °F)
Resolution	0.1 ppt, 1%

Specification	Details
Range	0 to 42, ppt 0 to 4.2%
<b>Dissolved oxygen</b>	
Calibration	One or two points, user-selectable to any value
Accuracy	±1.5% full scale
Resolution	0.1%, 0.01 ppm or mg/L
Range	0.0% to 199.9% saturation, 0 to 19.99 ppm or mg/L Salinity correction: automatic with conductivity probe Barometric pressure compensation: automatic
<b>Barometric pressure</b>	
Calibration	Factory calibration
Accuracy	±1.5 hPa (10 to 40 °C or 50 to 104 °F)
Resolution	1 mm Hg or 1 hPa 0.01 in Hg±
Range	225 to 900 mm Hg or 300 to 1200 hPa (8.86 to 35.43 in. Hg)

## General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

## Safety information

### NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

### Use of hazard information

#### ▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

#### ▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

#### ▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. <i>Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.</i>

## Certification

### Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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### FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

## Product overview

### NOTICE

Always disconnect power to the meter when electrodes are changed. Only use the meter as instructed in this manual or the meter performance can decrease.


The H-series handheld and benchtop meters operate with glass sensor electrodes with BNC connectors or non-glass probes with ISFET (ion sensitive field effect transistor) silicon chip sensors. The meters use one pH electrode (a BNC pH electrode or an ISFET pH probe) at a time. When the meter power is set to on, the meter automatically identifies the type of electrode that is attached.

The H-Series meters are available in eight models:

- Waterproof handheld meters with Bluetooth™ technology. Refer to [Figure 1](#):
  - H160G—pH and ORP
  - H170G—pH, ORP, conductivity, TDS and salinity
- Benchtop meters with a USB output. Refer to [Figure 2](#):
  - H260G—pH and ORP
  - H270G—pH, ORP, conductivity, TDS and salinity
  - H280G—pH, ORP, conductivity, TDS, salinity and dissolved oxygen (DO)
- Benchtop meters with Bluetooth™ technology and a USB output. Refer to [Figure 2](#):
  - H260GB—pH and ORP
  - H270GB—pH, ORP, conductivity, TDS and salinity
  - H280GB—pH, ORP, conductivity, TDS, salinity and DO


**Note:** The difference between the benchtop meters is that the GB benchtop meters have Bluetooth™ technology and a USB output, and the G benchtop meters only have a USB output. Unless noted, when the benchtop series is documented in this manual, the benchtop meter includes all of the H-series benchtop meters (both the G and the GB benchtop meters).

**Figure 1 Handheld meter**

1 Rubber dust caps	5 ISFET pH probe connector (8-pin)
2 Conductivity probe connection (12-pin, H170G only)	6 LCD display
3 BNC probe connector	7 Power button
4 3.5 mm phono jack for glass pH electrode, ISE, ORP or DO temperature sensors	8 Battery cover

**Figure 2 Benchtop meter**


1 Power button	6 ISFET pH probe connector
2 Conductivity probe connector (12-pin, H270G and H280G only)	7 AC power connector
3 USB connector	8 External reference connector
4 3.5 phono jack for glass pH electrode, ORP, ISE or DO temperature sensors	9 LCD display
5 BNC connector for glass pH electrode, ISE, ORP or DO (H280G only) probes	10 Battery cover

## Product components

Make sure that all components have been received. Refer to [Figure 3](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

**Figure 3 Handheld and benchtop meter components**

1 Waterproof handheld meter	4 Benchtop meter
2 SmartLogger II software	5 USB cable
3 AA Alkaline batteries (4x)	6 AC-DC power supply kit (power supply and three plugs: US, EU, UK)

## Installation

### ▲ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

## Electrical installation

### Connect to AC power

### ▲ DANGER

Electrocution hazard. If this equipment is used outdoors or in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used for connecting the equipment to its main power source.

## ⚠ WARNING


Fire hazard. Use only the power supply that is specified for this instrument.

The benchtop meters use AC power with an AC power adapter kit. Refer to [Product components](#) on page 8. The kit includes an AC-DC power supply, USB/DC adapter and AC power cord. Refer to [Figure 4](#) for AC power connections.

**Note:** Always set power to off before any power connections are made.

**Figure 4 AC power connection**

## Install the batteries

## ⚠ WARNING

Explosion hazard. Incorrect battery installation can cause the release of explosive gases. Be sure that the batteries are of the same approved chemical type and are inserted in the correct orientation. Do not mix new and used batteries.

## NOTICE

Do not tighten the screws too much or instrument damage can occur.

## NOTICE

Only do this procedure if the power to the meter is set to off or disconnected. Do not complete this operation with probes attached to the meter. Remove all probes from the meter or instrument damage can occur.


The meter uses AA alkaline or rechargeable NiMH batteries. To conserve the battery life, refer to [Configure the meter](#) on page 13 to configure automatic shutdown of the meter.

Refer to [Figure 5](#) to install the batteries.

### Items to collect:

- Phillips screwdriver (for handheld meters)
- Flathead screwdriver (for benchtop meters)
- AA Alkaline batteries (4x)

**Figure 5 Battery installation**


## User interface and navigation

### Keypad description

Figure 6 shows the handheld and benchtop meter keypads. Table 1 gives the function of the keys on the keypad and the meters that use that function.



**Figure 6 Keypad description**

**Table 1 Keypad functions**

Key	Action	H160G H260G H260GB	H170G H270G H270GB	H280G H280GB
	ON/OFF: Set the meter power to on or to off.	x	x	x
	pH mode	x	x	x

**Table 1 Keypad functions (continued)**

Key	Action	H160G H260G H260GB	H170G H270G H270GB	H280G H280GB
	mV mode	x	x	x
	ISE mode	x	x	x
	TDS mode		x	x
	Conductivity mode		x	x
	Salinity mode		x	x
	Calibration mode	x	x	x
	DO mode			x
	Setup mode	x	x	x
	Barometric pressure mode			x
	Time and date display	x	x	x
	Store a reading	x	x	x
	Recall a stored reading.	x	x	x
	Delete a stored reading.	x	x	x
	Scroll through values, setup screens and options.	x	x	x
	Set the backlight to on. After 2 minutes without a keystroke, the backlight is set to off.	x	x	x
	Start/stop the data log	x	x	x
	Select an option, setting or value.	x	x	x
	Cancel an option, setting or value.	x	x	x

## Display description

Figure 7 shows the measurement modes and values, data connection and storage options, battery status, temperature values, stabilization lock and connection statuses shown on the display.

Figure 7 Display

1 Data log indicator	10 Data storage memory location
2 Temperature and data values	11 Automatic temperature compensation (ATC) value (pH, conductivity, TDS or salinity)
3 Measured value	12 Calibration mode indicator
4 Stabilization lock	13 Temperature unit
5 Hold indicator	14 Low battery indicator
6 ISFET probe indicator	15 Bluetooth™ connection indicator
7 Measurement mode	16 PC data transfer icon
8 Measurement units	17 Setup mode
9 Storage options	

## Startup

### Set the power to on

Push the power button to set the power to on or off. Make sure that the power supply (AC power or battery power) is correctly installed.

### Configure the meter

Do the operations in order:

1. Set the power to on.
2. Push **SETUP** to put the meter in setup mode.
3. Use the arrow keys to select an option, then push **ENTER**:

Option	Function	Description
<b>CLr ALL CAL</b>	Clear calibration points	Erases all the calibration points. Make sure to calibrate the meter. <sup>1</sup>
<b>int</b>	Data log interval	Keeps up to 999 data points in intervals from 1 to 1999 seconds (default = 10) in the data log.
<b>year</b>	Year	Use the number keys to enter the correct year.
<b>date</b>	Date format	Sets the date to mm/dd/yy or dd/mm/yy format.
<b>m/d date</b>	Month and day	Use the number keys to set the correct month and day.
<b>time</b>	Time	Sets the correct time in a 24-hour format. Use the number keys to set the correct time.
<b>oFF</b>	Automatic shutdown	Sets the shutdown parameters from 1 minute to 000 minutes (continuous power). The meter beeps 1 minute before shutdown. Make sure to power cycle for the automatic shutdown to occur. Automatic shutdown is disabled during: data logging, Bluetooth™ transfers and when connected to the USB port.
<b>Snd</b>	Sound options	Sets sound alerts to on or off. Three sounds tell the user about different functions: <ul style="list-style-type: none"> <li>• <b>One beep:</b> Stabilization occurred with the stabilization lock on.</li> <li>• <b>Two beeps:</b> An error occurred. The error code shows on the display. Refer to <a href="#">Troubleshooting</a> on page 17.</li> <li>• <b>Three beeps:</b> Measurement stability in calibration mode, regardless of the stabilization lock setting.</li> </ul>
<b>°C °F</b>	Temperature display	Sets the temperature unit from °C to °F.
<b>READY</b>	Stabilization lock	Sets the stabilization lock to on or off: <ul style="list-style-type: none"> <li>• The "READY" icon shows when an endpoint occurs and the stabilization lock is set to on.</li> <li>• After stabilization, the display locks the value and the meter ignores slight measurement changes in pH, conductivity or TDS.</li> <li>• Sets the stabilization lock to off during titrations or slight change detection. The display automatically unlocks after a significant measurement change is found.</li> </ul>

<sup>1</sup> This does not erase calibration data for an ISFET probe.

4. Push **ENTER** to keep the changes and go back to setup mode.  
**Note:** To exit and not keep the changes, push **CANCEL**.

## Operation

### Configure the Bluetooth™ wireless connection

Transfer data and control several instrument functions from a remote location to a PC with the Bluetooth™ wireless connection.

**Note:** The arrow icon shows on the display when the meter is in communication with the SmartLogger II software on a PC. The arrow icon flashes when data is transferred. Refer to the SmartLogger II documentation to setup the wireless connection to a PC.

1. Set the power to on. The meter looks for the Bluetooth™ connection with a PC.
2. On the PC, select the option to find or add a Bluetooth™ device. Next, the PC prompts if the user wants to pair with a found device.
3. When the PC prompts for a pass key or PIN, enter the PIN (default = **1234**). The PC shows if the pair is successful and the meter shows the Bluetooth™ icon.

4. If the pair is not successful, do steps 1 through 4 again.

**Note:** If there is more than one meter in range, each meter is identified with the model number and the serial number (e.g., H170G LP SN1 2755).

## Send data to a PC

Transmit data in storage to a PC and see real-time measurement values from a remote location on a PC. Refer to the SmartLogger II PC software guide for operation instructions for the PC to USB connection.

1. Set the power to off.
2. Connect the USB cable to the USB port of the PC.
3. Use the USB drivers to make a USB connection to the PC.
4. Set the power to on. Refer to the SmartLogger II documentation for data transfer information.

## Maintenance

### Replace the batteries

The display shows "bAt" when the batteries are too low to give a reliable measurement. The low battery icon (refer to [Display description](#) on page 13) shows when there is approximately 25 hours of battery power. Measurement errors are possible when the batteries are low. Refer to [Install the batteries](#) on page 10 to replace the batteries.

**Note:** The date and time must be set again when the batteries are removed or are fully discharged.

### Clean the instrument

Clean the exterior of the instrument with a moist cloth and a mild soap solution and then wipe the instrument dry.

### ISFET pH probe maintenance

#### ▲ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

#### ▲ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

The expected life of an ISFET probe is approximately 18 months. The reference electrode has a KCl gel that is diluted over time. The reference is sealed and is non-refillable. Replace the probe when it becomes difficult to calibrate.

- Keep the probe dry with the protective shroud on when not in use.
- Clean oil, fat, food particles, starch, protein or other materials from the probe tip after use.
- Never use sharp metal objects (e.g., a needle, a pin, etc.) to clean the sensor surface.
- Remove the probe from environments with static electricity. Electrostatic discharge (ESD) can permanently damage the probe.
- Remove the probe from environments that will damage the sensor, such as hydrofluoric acid or abrasive samples.

- Remove the probe from environments that will damage the epoxy materials used in the probe tip (e.g., keep the probe away from acetone, toluene, methylene chloride, xylene and other strong organic solvents).
- Do not use the probe in temperatures more than 60 °C (140 °F). Thermal cycling can decrease the life of the probe.
- For semi-solids use, carefully twist the probe to make sufficient contact with the sample to the sensor.
- When semi-solids are tested, make sure that solid objects (i.e., bone or gristle) do not scratch the sensor.
- Cool samples to room temperature to maximize probe life.
- Always use new buffers and new rinse solutions.

### Prepare the probe for storage

#### NOTICE

Do not use the probe for long-term pH measurement applications.

**Note:** *No electrode storage solution is necessary.*

A new probe usually has visible reference gel at the tip of the probe, as well as in the interior of the rubber dust cap. To remove the gel, carefully clean with a soft-bristled toothbrush and mild soapy water (a few drops of dish soap in a warm cup of water). The gel can continue to show for two to five days. Do not put the rubber dust cap on the probe until all of the reference gel is removed. Do not complete the steps that follow until all reference gel is removed from the probe and the interior and exterior of the rubber dust cap.

### Prepare the probe for storage

1. Put the new probe (or the probe that was in extended storage) in pH 7 buffer for at least 5 minutes.
2. Stir the probe in pH 7 buffer solution to dislodge air bubbles.
3. Use fresh deionized water to rinse the probe.
4. Dry the probe with a lint-free cloth.
5. Put the rubber dust cap back on.
6. Keep the probe in dry storage when not in use.

### Calibrate the ISFET probe

Refer to [Prepare the probe for storage](#) on page 16 before calibration.

#### NOTICE

Do not use the probe for applications that cycle between hot and room temperature samples.

**Note:** *Remove other active or non-active measurement devices during pH or conductivity measurement. Other devices, even AC power, can cause interference.*

**Note:** *Only do calibrations away from the sun. Direct sunlight can cause unstable readings or difficulty in calibration.*

1. Connect the probe to the meter.
2. Set the meter power to on.
3. Clean the probe with new deionized water and dry with a lint-free cloth.
4. Put the probe in the pH 7 buffer.
5. Clean the probe with deionized water and dry with a lint-free cloth.
6. Put the probe in the second buffer (pH 4 or pH 10).
7. Read the results.
8. If the result is not correct, the probe is not correctly hydrated.  
Soak the probe for another 5 minutes in pH 7 buffer, then do the calibration again.

## Clean the ISFET probe

### NOTICE

Do not use sharp metal objects (a needle, a pin, etc.) to clean the sensor. This can scratch the sensor and cause permanent damage to the probe.

For use in dairy, cheese or meat applications, soak the probe in Pepsin Cleaning Solution for 15 minutes before the probe is cleaned.

Regularly clean a non-glass probe:

1. Remove the rubber dust cap from the probe, then rinse with new deionized water.
2. Use a soft-bristle toothbrush and a mild detergent (a few drops of dish soap in a warm cup of water) to carefully clean the probe.
3. Rinse with new deionized water to remove all debris from the sensor surface.
4. Dry the probe with a lint-free cloth.
5. Calibrate the probe again. Refer to [Calibrate the ISFET probe](#) on page 16, then [Prepare the probe for storage](#) on page 16.

## Repair the ISFET probes

### NOTICE

Do not use the probe to find out if the buffer is above 60 °C (140 °F). If the probe is suddenly put into very hot liquid, the probe can be permanently damaged.

Monitor how long the probes are in dry storage. If the probes are in storage for an extended period of time, the KCl gel at the reference junction can crystallize.

1. Heat pH 7 buffer to approximately 45 to 60 °C (115 to 140 °F).
2. Soak the probe for 2 minutes.
3. Put the probe in room temperature pH 7.00 buffer and let cool.

## BNC electrode maintenance

Prevent unstable readings:

- Keep the probe in an electrode storage solution.
- Start measurements with at least a 2-point calibration. Update often with 1, 2 or 3-point calibrations.
- Use new buffers and new deionized rinse solution.
- Use buffer solutions with pH values no larger than 3 pH units apart. Ideally, the buffers bracket the anticipated pH values of the unmeasured samples.
- Use deionized water to rinse residual buffer and sample solutions from the probe after calibration and measurement.
- Calibrate at the same temperature as the sample solution. Although the meter has an ATC, get the best results when the calibration buffers and the sample are the same temperature.
- Keep the connectors clean and dry. Dirty or damp connectors can cause unstable readings.

## Troubleshooting

Problem	Possible cause	Solution
No display	Automatic shutdown set the power to off.	Set the power to on again.
	There is no power.	Replace the batteries.

Problem	Possible cause	Solution
Unstable reading	The probe is dirty.	Clean the probe.
	The probe or the meter connectors are dirty.	Clean the probe contacts on the cable connector. Clean the meter with methanol on a cotton swab. Let dry completely. Connect the probe to the meter again.
	No flow in the reference junction	Clean the warm buffer.
	ISFET probe is not correctly hydrated.	Soak the probe for at least 5 minutes in a pH 7.00 buffer.
	Interference from other devices	Remove other devices from the solution.
	Interference from direct sunlight	Use protection for the probe from direct sunlight.
	The probe is in a very low ionic strength solution.	A stable reading is not possible.
	The pH or temperature of the solution changes.	A stable reading is not possible until pH and temperature are constant.
The meter continually shows -2.00 or 19.99 with an electrode attached. The "ISFET" icon does not show when ISFET probe is attached.	The ISFET probe is not sensed by the meter or the probe.	Set the meter power to off. If the "ISFET" icon does not show, replace the ISFET probe.
	Out of calibration	<a href="#">Calibrate the ISFET probe</a> on page 16.
	The probe is not in the solution.	Put the probe in liquid. Carefully shake the probe to make sure that air bubbles are not caught on the sensor surface.
	No probe is connected to the meter.	Set the meter power to off. Connect the pH probe. Set the meter power to on.
	The probe is dirty.	Clean the probe.
	The probe is damaged.	Replace the probe. Contact technical support.
Screen flashes during calibration.	The probe sensor surfaces are dirty or it is necessary for the probe to be conditioned again.	Clean and condition the probe again.
	No flow in the reference junction	Clean the warm buffer.
	The buffers are contaminated or expired.	Calibrate with new buffers again.
	Interference from other devices in the solution	Remove all devices from the solution.
	The battery is low.	Replace the batteries if the battery icon shows low battery power.
	Interference from direct sunlight	Use protection for the probe from direct sunlight.
	The probe is too old.	Replace the probe.

## Error codes

Table 2 shows the codes that can occur for various reasons. Error codes show instrument malfunction or user error.

**Table 2 Error code descriptions**

Error Code	Description	Solution
E02	The ISFET pH probe is damaged.	Replace the probe. Contact technical support for probe replacement information.
E03	Clean the probe.	Clean the probe. If the error continues, replace the probe.
E04	Glass pH probe slope error. The slope is smaller than 85% or larger than 102% of 59.16 mV per pH unit.	Clean the probe. If the error continues, replace the probe.
E06	ISFET pH probe slope error.	Clean the probe. Soak in 40 °C (113 °F) pH 7.00 buffer for 2 minutes. Calibrate again. If the error continues, replace the probe.
E07	Clean the probe.	Clean the probe. Make sure that no air bubbles on the sensor surface or foreign objects or materials on the sensor.
E08	Too long to calibrate. Signal is not stable during calibration.	Clean the probe. Disconnect stirrers and other AC power sources. Make sure that the calibration solution temperature is constant. Replace the probe if the error continues.
E13	The ISFET pH probe temperature sensor is damaged.	Replace the probe. Contact technical support for probe replacement information.
E14	The battery is extremely low.	Replace the batteries immediately. Damage to the accuracy and function of the meter could have occurred.
E15	Replace the battery immediately. The accuracy and function of the meter could be compromised.	Contact technical support for service information.
E20	The conductivity temperature sensor is damaged.	Replace the probe. Contact technical support for probe replacement information.
E25	The slope is smaller than 60% or larger than 140% of the nominal.	Replace DO probe membrane and fill solution. Replace the probe if the error continues. Nominal = 0.1822 ppm/mV.
E26	The mV reading is more than $\pm 10$ mV from nominal.	Replace DO probe membrane and fill solution. Replace the probe if the error continues. Nominal = 0.0 mV at 0% saturation; 45mV at 100%.
E27	DO probe temperature error	Attach the temperature sensor to the 3.5 mm phono jack. DO readings are highly temperature dependent and a temperature sensor must be attached.
E28	Barometric pressure sensor error	Contact technical support for repair information.
E30	ISE electrode calibration error	The slopes are not the same sign or are not within 25% of each other. Calibrate again in the correct solutions.
E40	Unrecognized host command	Command from a host PC is not found. Use only valid commands.
E42	Invalid input	The value entered during the setup is invalid. Enter a different value.
E44	No probe is installed	No probe is installed for the applicable parameter. Set the meter power to off. Install the correct probe. Set the meter power to on.

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)