
M295Ag
Silver Electrode

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The design principles of the M295Ag Silver Electrode

Silver Electrode

Operating Instructions

The M295Ag Silver Electrode can be used directly or after coating for the determination of silver halides either by direct measurements or by potentiometric titration.

It can be used together with a standard REF601 ($\text{Hg}/\text{Hg}_2\text{SO}_4$) reference electrode or with a double junction reference electrode eg. REF251 or REF451.

Preparing for measurement

1. Remove the protection cap from the electrode.
2. Before starting a measurement, rinse the electrode with distilled water.
3. Check that the electrode body and sensing element are not damaged.

Cleaning and maintenance

Silver rod

The sensing part of the electrode can be polished with finishing abrasive bands, part no. B90X550 or with silver polish.

In order to have reproducible and accurate results with the electrode, a few treatments are required. These treatments will depend on the last measurement performed:

Cleaning

The electrode must be carefully cleaned before use, but in many cases a thorough rinsing with distilled water will suffice. Owing to the high chemical resistance of the electrode, most cleaning chemicals and solvents may be used (the electrode tip and stem consist of pure silver and PTFE respectively).

- To remove most electrode deposits which cannot be removed with pure water: flush with dilute nitric acid and/or ammonium hydroxide followed by distilled water.
- In case of severe contamination, clean the electrodes with silver polish and then flush with acetone and distilled water.

Storage

Replace the protection cap, and store the electrode in its box at room temperature.

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