ISE Ca²⁺ 96 60C

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


www.hach.com

© Hach Company/Hach Lange GmbH, 2013. All rights reserved. Printed in U.S.A.

08/2013, Edition 1

ISE Ca²⁺ 96 60C EN

Additional information is available on the manufacturer's website.

Caution
 ⚠ *Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.*

Introduction
 This probe is a combination calcium selective electrode with a non-refillable double junction reference. Use the probe with a meter or titrator for calcium determination in aqueous solutions.

Technical information

Range	0.02 mg/L to 4 g/L (5 × 10 ⁻⁷ M to 1 M) Ca ²⁺
pH range	3.5 to 11
Linear range	5 mg/L to 4 g/L (10 ⁻⁴ M to 1 M) Ca ²⁺
Slope	27 mV/decade, 90 to 110% at 25 °C (77 °F)
Ionic strength adjuster (ISA)	3 M KCl
Operating temperature	0 to 50 °C (32 to 122 °F)
Interferences	Ba ²⁺ , Al ³⁺ , Sr ²⁺

Preparation
 Refer to the instrument manual for detailed procedures.

1 Connect the probe
 Before initial use, soak the probe for 30 minutes in a 10 mg/L Ca²⁺ solution.

2 Calibration
 Add 1 mL ISA per 100 mL of standard. Use a minimum of two standards in an increasing order of concentration.

3 Measurement
 Add 1 mL ISA per 100 mL of sample.

Maintenance and storage
 To clean: rinse with distilled water before each measurement.
 Storage: clean and dry the probe. Install the sensor protection cap for storage.

Troubleshooting
 Make sure that the ionic strength and temperature of the samples and standards are the same.

Make sure that the stir speed is the same for all the standards and samples.
 Remove deposits from the membrane.
 Replace the probe if damaged.

Warranty
 6 months for manufacturing faults only. The warranty does not cover normal wear or improper use.

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