Radiometer Analytical SAS

Certified Conductivity Standards Leitfähigkeitssäfte mit Zertifikat Etalons de conductivité certifiés		
Type	Value	Part no.
KCl 1D	111.3 mS/cm ±0.5% (25°C)	S51M001
KCl 0.1D	12.85 mS/cm ±0.35% (25°C)	S51M002
KCl 0.01D	1408 µS/cm ±0.5% (25°C)	S51M003
NaCl 0.05%	1015 µS/cm ±0.5% (25°C)	S51M004
NaCl 25 µS/cm	25.0 µS/cm ±5% (25°C)	S51M013

2-pole Conductivity Cells

D31M102 • Printed by Radiometer Analytical SAS • France • 2005-09E

Conductivity Solutions Leitfähigkeitslösungen Solutions de conductivité		
Type	Value	Part no.
KS910	0.1M KCl (12.88 mS/cm ±1% (25°C))	C20C250
KS220	0.01M KCl (1.413 mS/cm ±1% (25°C))	C20C270
KS330	0.001M KCl (146.9 µS/cm ±2.5% (25°C))	C20C280


Operating Instructions

Mode d'Emploi


RADIOMETER ANALYTICAL SAS
72 rue d'Alsace, 69627 Villeurbanne Cedex, France
E-mail: radiometer@analytical.com Web: www.radiometer-analytical.com
Tel. : +33 (0)4 78 03 38 38 - Fax : +33 (0)4 78 68 88 12

Radiometer
analytical
A Hach Company Brand

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Contents

Introduction	2
Operation	3
Cell constant	3
Calibration procedure	3
Sample measurement procedure	4
Cleaning and storage	4
Measuring tips!	5
Theory	6
Conductivity	6
Resistivity	7
The cell constant	7
Temperature effects	8
Cable resistance	10
Cable capacitance	10
Troubleshooting	11
Applications	12
Accessories	13
Adapters	13
Certified conductivity standards	13
Additional calibration solutions	13
Specifications	14

Introduction

Accessoires

Conductivity is measured in a wide range of industries and gives an indication of total ionic concentration of the sample solution.

Radiometer Analytical offers a selection of 2-pole conductivity cells for a wide variety of applications.

2-pole cells have a traditional design based on two plates of platinum. They are ideal for routine measurement of conductivity and for use with a sample changer due to easy rinsing.

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English

Cell constant "K"

Each cell is calibrated during manufacturing and the "actual" cell constant value "K" is displayed on the calibration certificate supplied with the conductivity cell.

The "K" value is expressed in cm^{-1} . Cell constant intervals for each cell are given on page 14.

Note: the CDM92, CDM210 and CDM230 Conductivity Meters express "K" in cm^{-1} . However, some older meters use "K" values expressed in " $\text{cm}^{\prime\prime}$ ".

The following method is used in our laboratories to determine the cell constant:

1. The conductivity cells are placed in deionised water and left for 1 to 2 hours.
2. The cell constant "K" is determined using a CDM210 or CDM230 Conductivity Meter and a 0.1 D or 0.01 D KCl standard solution.

Calibration procedure

1. Place the conductivity cell in deionised water and leave for 1 to 2 hours.
 2. Rinse the conductivity cell thoroughly with standard solution, to eliminate traces of deionised water.
 3. Mount the cell in a suitable cell holder. The SAM7 Sample Stand is ideal for this.
- Note:** the conductivity cell should be placed at approximately 1 cm from the bottom and the walls of the beaker.
4. Connect the conductivity cell to the conductivity meter.
 5. Carry out the calibration in a standard solution with a concentration as close as possible to that of the sample. The following solutions may be used.

For platinised cells: use 0.1 M or 0.1D KCl standard at 25°C.

For non-platinised cells: use 0.01 M or 0.01D KCl standard at 25°C.

The cell constant may, however, change due to contamination and should be regularly checked using the method described above.

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Sample measurement procedure
The sample measurement conditions must be as close as possible to the calibration conditions, e.g. temperature of the sample and position of the conductivity cell.

Before performing measurements, rinse the conductivity cell thoroughly with the sample.

Note: check that the conductivity cell is immersed in the sample; refer to page 14 for details on the minimum immersion depth. The conductivity cell should ideally be placed at approximately 1 cm from the bottom and the walls of the beaker.

Recommended conductivity cells by application are given on page 12.

Cleaning and Storage

The 2-pole conductivity cells require very little maintenance to be kept in optimal condition.

WARNING: do not wipe the plates of a platinised conductivity cell with paper tissue, as this risks removing the platinised layer. It is recommended to dip the conductivity cell several times in a beaker of deionised water.

- Between measurements, rinse the cell with deionised water.
- If the cell has been exposed to a solvent immiscible with water, clean the cell with a solvent miscible with water e.g. ethanol or acetone, and rinse carefully with water.
- Note:** do not use acetone on epoxy body conductivity cells.
- In case of deposits on the cell, dip the cell in RENOVO•N Normal Cleaning Solution for approximately 1 hour.
- IMPORTANT:** rinse thoroughly with distilled water after cleaning.
 - For short and long-term storage, it is recommended to rinse the cell and place it in deionised water.
- Note:** it is possible to store the conductivity cell in dry conditions. However, before use it is necessary to re-hydrate the cell. To do so, place the conductivity cell in deionised water and leave for 8 hours.

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Measuring tips!

- In order to obtain accurate conductivity measurements, you must be sure of the cell constant value. It is therefore recommended to control the cell constant periodically. In order to do this, perform a daily conductivity measurement in a standard. If the result obtained is in accordance with the theoretical value, continue your measurements, if not, perform a new calibration to re-determine the cell constant.

- The calibration conditions must be as close as possible to the working conditions you use when performing measurements in a sample.
- Conductivity measurement is temperature dependant (if the temperature increases so will the conductivity value). When performing conductivity measurements, it is advisable to choose a cell with a built-in temperature sensor or use a separate temperature sensor.
- For most accurate measurements, it is recommended to thermostat calibration solutions and samples.
- Normally conductivity results are reported after correction to a reference temperature. Radiometer Analytical's conductivity meters do this automatically.

Theory

Conductivity

Conductivity is the ability of a solution to pass an electric current between two electrodes. The current in a solution is carried by charged ions, therefore, the more ions present in a solution, the higher the conductivity. A solution will resist current flow if a small amount of ions is present.

In practice, an alternating current at a suitable frequency is used to limit polarisation effects.


Fig. 1: Ionic migration in solution

How conductivity is measured?

In a traditional 2-pole cell, a voltage (V) is applied to two flat platinum plates immersed in a solution, and the resulting current (I) is measured. The resistance of the solution (R) can be calculated using Ohm's law ($V = R \times I$).

$$R = V/I$$

The conductivity of a solution, κ , is calculated on the basis of the conductance, (G), and is given by the expression:

$$\kappa = G \bullet K \text{ (S/cm)}$$

κ = conductivity (S/cm)
 G = conductance (S), where, $G = 1/R$
 K = cell constant (cm⁻¹)

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

Once the conductivity of a solution is known the resistivity, ρ , can, be calculated as follows:

$$\rho = \frac{1}{K} (\Omega \cdot \text{cm})$$

The cell constant

2-pole conductivity cell

Fig. 2: Traditional 2-pole cell

The cell constant is based on the geometry of the cell and can in theory be calculated from the dimensions of the cell. As the cell dimensions change, the cell constant varies as a ratio of l to A . See figure 2 above.

$$K = \frac{1}{A} (\text{cm}^{-1})$$

K = cell constant (cm^{-1})
A = surface area of the electrodes (cm^2)
l = distance between the electrodes (cm)

In the theoretical cell using 1 cm squares of platinum, 1 cm apart, the cell constant (I/A) is 1.0 and the conductance reading in siemens (S) is numerically equal to the conductivity in S/cm . For low conductivity solutions, the electrodes can be placed closer together (reducing I) and/or increasing the surface to give cell constants of 0.1 cm^{-1} . This increases the conductance between the platinum plates and makes it easier for the electronic circuit to obtain a result.

However, in reality it is not possible to calculate the cell constant as the surface area cannot be measured due to the roughness of the plates and because the current is not confined to the space between the plates. This is why we calibrate with traceable conductivity standards of known conductivity. For the calibration procedure; refer to page 3.

Note: It is possible to notice changes in cell constant values over time and it is therefore highly recommended to check the cell constant periodically. It is quite normal that for the same conductivity cell you will find a slight variation in the conductivity of the solution used for calibration. This is especially true for non-platinised cells, due to the polarisation effect.

Temperature effects

The conductivity of a solution will increase with temperature. In order to compare measurements taken at different temperatures, the conductivities are corrected to a reference temperature using the temperature coefficient.

For example: ultrapure water has a temperature coefficient of $5.2\%/\text{°C}$ whereas ionic salts are at the $2\%/\text{°C}$ level. Acids, alkalis and concentrated salt solutions are at a somewhat lower level, typically $1.5\%/\text{°C}$.

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
Determination of the temperature coefficient (θ)

By measuring the conductivity of a sample at temperature T_1 close to T_{ref} and another temperature T_2 , you can calculate the temperature coefficient by using the following equation:

$$\theta = \frac{(\kappa(T_1) - \kappa(T_2)) \cdot 100}{(T_1 - T_2) \cdot \kappa(T_1)} \quad \text{Unit: \%/\text{°C}}$$

T_2 should be selected as a typical sample temperature and should be approximately 10°C different from T_1 .

Note: do not forget to disable the "Temperature correction" function on your conductivity meter.

When temperature correction is selected on the meter, the conductivity at the reference temperature is displayed. Note that the correction is accurate only within a limited temperature range around T_1 and T_2 . On the figure below:

$$T_1 = 26^\circ\text{C}, T_2 = 14^\circ\text{C} \text{ and } T_{\text{ref}} = 25^\circ\text{C}$$


Fig. 1: Transport des ions en solution

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K = constante de cellule (cm^{-1})

Cable resistance

A cable of a given length has a resistance. The influence of the cable resistance on measured conductance (G_m) is as follows:

$$G_m = \frac{G_s}{1 + (R_c \cdot G_s)}$$

G_m = measured conductance (S)

G_s = conductance of the solution (S)

R_c = cable resistance (Ω)

The cable resistance influences measurements performed using a 2-pole conductivity cell because the cable resistance will be in series with the resistance of the solution.

The cable resistance is mainly important for measurements of high conductance.

Note: enter the cable resistance value when programming the conductivity meter.

Cable capacitance

A cable of a given length has a given capacity. The cable capacitance influences the measurements of very low conductance. It can be necessary to enter a cable capacitance in order to obtain accurate conductivity measurements at low concentrations (below 4 $\mu\text{S}/\text{cm}$).

Note: most cables used with our 2-pole conductivity cells have low capacitance. Consequently, it is not necessary to enter this value, for standard length cables (approximately 1.30 m), as any value below 350 pF will not affect measurements. However, it is necessary to enter a cable capacitance for cables with lengths greater than 1.30 m or when using conductivity cells with built-in temperature sensor.

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Troubleshooting

English

It is quite normal that for the same conductivity cell you find that the cell constant can vary slightly.

If your cell constant deviates a little from the previous determined value, clean your conductivity cell, see page 4 and re-determine the cell constant.

Please note, It is recommended to control the cell constant periodically in order to ensure the reliability of conductivity or resistivity measurements.

WARNING: do not wipe the plates of a platinised conductivity cell with paper tissue, as this risks removing the platinised layer.

For low conductivity measurements. Work in a stable environment avoiding contamination from the air e.g. CO₂. It is also advisable to continuously renew the solution, e.g. use a flow system.

Applications

Application/Type	CDC241-9	CDC267-9	CDC641T	CDC741T	CDC745-9	CDC749	XE100
Various aqueous and non-aqueous media			✓				✓
Use with sample changer					✓		
Built-in temperature sensor				✓	✓		
Continuous measurements	✓					✓	
Microsamples						✓	
Titration	✓						
Epoxy body					✓		
Viscous media	✓						
Highly resistive media		✓					
Field use					✓		
Use in glass tubes						✓	

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Adapters

Plug adapter UHF-F/MAB6-M	A94P001
Plug adapter MAB6-F/UHF-M + PYE	A94P002
Plug adapter MAB6-F/Perena + Banana	A94P003

Certified conductivity standards

Each 500 ml bottle of Radiometer Analytical certified conductivity standards is placed in an airtight tin to protect them from evaporation and microbiological growth, allowing them to keep their values for a full two years from the date of production. Each bottle comes with its own certificate of conformity and traceability following ISO Guide 31⁽¹⁾.


Additional calibration solutions

Radiometer Analytical also produces a range of molar KCl calibration solutions, delivered in 500 ml bottles.

⁽¹⁾ ISO Guide 31 provides guidelines for certified reference materials

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* Interval of production

Note: some older meters use "K" values expressed in "cm".

** For screw cap electrodes, refer to the table below.

Cables	CL117	CL119	CL136	CL336	CL319
Cable resistance (Ω)	0.180	0.155	0.145	0.350	0.360
Cable capacitance (pF)	170	130	70	200	390