EN

ISENA381 probe Introduction

The probe is a digital, combination ion selective probe with a double junction, refillable junction reference and a built-in temperature sensor.

Warning
 ⚠ *Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.*

User manual
 A user manual with all instructions is on the manufacturer's website.

Collect the following items:
 Compatible meter; beakers; Na⁺ standard solutions; rinse bottle; ISA (0.4 g) powder pillows; prepared ISA rinse solution; lint-free cloth.

- 1 Connect the probe**
- 2 Calibration**
Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.
- 3 Measurement (direct method)**
- 4 Storage**
 ⚠ *Solution: 25 mL of 100 mg/L Sodium standard plus one ISA pillow (0.4 g).*

Technical information

Range	0.023 mg/L (1x10 ⁻⁶ M) to 23,000 mg/L (1 M) Na ⁺
Slope	-59 mV/pNa (90 to 110% at 25 °C per Nernstian theoretical value)
Operating temperature	0 to 50 °C (32 to 122 °F)
Minimum sample volume	25 mL
Filling solution	0.02 M NH ₄ Cl

Warranty
 1 year on the probe. This warranty covers manufacturing defects, but not improper use or wear.

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