

Brackish And Seawater Methods List for Automated Ion Analyzers

Flow Injection Analysis



19 January 2018

Introduction


Lachat Instruments has many methods for the analysis of nutrients and other ions in seawater and brackish waters. A key feature of these methods is that Lachat's brackish chemistries **do not require matrix matching**.

→The methods can be used with samples over the range of 0-35 ppt salinity!

→Standards are prepared in DI water as well.

→All Lachat methods are designed to prevent the "memory effect" (i.e. carryover), which can be problematic for segmented flow analyzers.


The data shown below is from a high range nitrate method, which covers the range of 5-500 mg N/L. The high standard was injected in duplicate, followed by 4 replicate blank injections:



Carryover Study: 500 mg N/L followed by 4 blanks; Carryover passed

Filename: 5-500 Supp run 3 .omn Date: 11 June 2009

Below is a calibration Diagram, covering the range of 1.25-100 µg N/L as ammonia!



File Name: Supp 1 10 27 08.OMN

Acq. Date: 27 October 2008

The pages that follow contain information on Lachat's brackish water methods. These cover a wide range of analytes and analyte levels. Ranges from sub single ppb to 10's or even 100's of ppm are available.



Form:

The method either determines this form of the analyte or converts the analyte to this form for determination.

00	Form given by previous three numbers	01	Phosphate (PO_4^{3-})
02	Calcium (Ca^{2+})	03	Potassium (K^+)
04	Nitrate (NO_3^-)	05	Nitrite (NO_2^-)
06	Ammonium (NH_4^+), Ammonia (NH_3)	07	Chloride (Cl^-)
08	Boric Acid (H_3BO_3)	09	Iodide (I^-)
10	Sulfate (SO_4^{2-})	11	Sulfite (SO_3^{2-})
12	Fluoride (F^-)	13	Chromium (VI) (Cr)
18	Total Iron ($\text{Fe}^{2+} + \text{Fe}^{3+}$)	19	Iron (II) (Fe^{2+})
21	Bromide (Br^-)	23	Molybdenum (VI) (Mo)
24	Hydronium (H_3O^+ , H^+)	25	Hydroxide (OH^-)
26	Magnesium (Mg^{2+})	27	Silicate (SiO_2)
29	Sulfide (S^{2-})	31	Calcium carbonate (CaCO_3)
32	Sodium cation (Na^+)	33	Aluminum (inorganic) (Al)
34	Aluminum (organic) (Al)	35	Chlorate (ClO_3^-)
36	Hypochlorite (OCl^-)	40	Perchlorate
42	Sulfur dioxide		

Chemistry:

Some analytes have more than one chemistry.

Example:

Ammonia	10-107-06-1	phenolate, phenate
	10-107-06-2	salicylate
	10-107-06-5	gas diffusion

Concentration:

Each range of concentrations for an analyte is given by a single letter. See the methods list for the ranges. Some methods cover more than one range.

Heaters:

Standard heater: Standard heaters have a 175 cm section of 0.032" i.d. (0.8mm) and a 650 cm section of 0.032" i.d. tubing

Non-standard heater: Has a different type and/or length of tubing than that listed above. (Controller and heater block are the same; only the tubing is different).

In the list of methods that follows:

Designation in the methods list means the method is EPA accepted as equivalent for NPDWR, NPDES, or both

^ Designation in the methods list means the method is equivalent for NPDES reporting under the MUR

NEW!
MEMI

^R **Designation means the manifold is compliant with RoHS-2.** (Verify with your sales person that the method can be sold into the EU).

Flow Injection Analysis

Ammonia

10-107-06-5-J ^{RA}	0.01-1.0 0.1-20	0.002 0.02	mg N/L as NH ₃	Waters	Gas Diffusion method. Salicylate/DCIC. May be used for TKN as well as brackish/saline samples. 660 nm. Requires a standard heater.	16-Jan-15
31-107-06-1-B ^{RA}	5 – 600 0.36-42.86	0.7	µg N/L as NH ₃ µM N/L as NH₃	Brackish / Seawaters	Alkaline phenol-based method; Can be used for determination of samples w/ 0 to 35 ppt salinity; 630 nm. Requires a standard heater. NPDES Equivalent (350.1)	18-Sep-08
31-107-06-1-F ^A	0.005 – 2.0 0.36-142.86	0.002	mg N/L as NH ₃ µM N/L as NH₃	Brackish / Seawaters	Alkaline phenol-based method; 630 nm. Can be used for determination of samples w/ 0 to 35 ppt salinity; Requires a standard heater. NPDES Equivalent (350.1)	12-Nov-07
31-107-06-1-G ^{RA}	1.25 – 100 0.089-7.143	0.41	µg N/L as NH ₃ µM N/L as NH₃	Brackish / Seawaters	Alkaline phenol-based method; DCIC. 630 nm 2-cm detector method; QC8500 only; Can used for determination of samples w/ 0 to 35 ppt salinity; Requires a non-standard heater. NPDES Equivalent (350.1)	26-Jan-10
31-107-06-1-H ^{RA}	0.25 – 30.0 0.018-2.143	0.025	mg N/L as NH ₃ mM N/L as NH₃	Brackish / Seawaters	Alkaline phenol-based method; 630 nm. high range method; Can used for determination of samples w/ 0 to 35 ppt salinity; <u>Ultra-High Throughput method</u> (>120 samples/hr) Requires a standard heater.	31-Oct-08
31-107-06-1-I	5-500 0.3571-35.71	0.47	µg N/L as NH ₃ µM N/L as NH₃	Brackish / Seawaters	Alkaline phenol-based method, citrate/tartrate buffer. 630 nm, Requires a standard heater.	21-Feb-12
31-107-06-1-Q ^A	0.005-2.0 0.36-142.86	0.0022	mg N/L as NH ₃ µM N/L as NH₃	Brackish / Seawaters	Alkaline phenol-based method, citrate buffer. 630 nm, Requires a standard heater. NPDES Equivalent to 349.0	17-Aug-10

Bromide

30-135-21-1-B ^R	5.0 – 60.0 0.0625-0.751	0.22	mg Br ⁻ /L mM Br⁻/L	Brackish / Seawaters	Phenol Red method. 590 nm. Follows Standard Methods (4500-Br-D) 590 nm.	3-Sep-03
----------------------------	-----------------------------------	------	---	-------------------------	---	----------

Method No	Range	MDL	Units	Matrix	Comments	Rev Date
Chromium						
31-124-13-1-A ^R	2– 200 0.038-3.85	0.66 0.0127	µg Cr/L as Cr(VI) µM Cr/L as Cr(VI)	Brackish / Seawaters	Hexavalent chromium in seawater/brackish waters. Diphenylcarbazide; 540 nm.	24-Aug-09
Iron						
31-126-18-1-A ^R	0.5-30 mg Fe/L	0.024	mg Fe/L	Brackish / Seawaters	Total soluble iron as Fe (II and III); TPTZ indicator. 600 nm. Inert sample probe required.	15-Sep-03
31-126-18-1-B ^R	0.05 – 0.500 0.895-8.95	0.004	mg Fe/L µg Fe/L	Brackish / Seawaters	Total soluble iron as Fe (II and III); TPTZ indicator. 600 nm. Inert sample probe required.	15-Sep-03
31-126-19-1-A ^R	0.50 – 30.0 0.00895-0.537	0.23	mg Fe/L mM Fe/L	Brackish / Seawaters	Total soluble iron as Fe (II); TPTZ indicator 600 nm. Inert sample probe required.	26-Nov-08
Kjeldahl Nitrogen (TKN)						
10-107-06-5-J	0.1-5.0 0.25-20	0.02 0.05	mg N/L	Waters	Kjeldahl Digests, Salicylate/DCIC 660 nm. copper catalyst. Gas diffusion method. Sea/brackish water. Can also be used for Ammonia	26-Sept-12
10-107-06-6-D [^]	0.5 – 20	0.25	mg N/L	Waters	Kjeldahl digests; copper catalyst ; inline distillation method; NPDES Equivalent (351.2); samples w/ particulates not suitable. 660 nm. <u>Can be used with brackish/seawater digests.</u> Requires an in-line module and a standard heater or two heated channels (with one heater non-standard).	31-Jul-09
Nitrate + Nitrite						
30-107-04-1-C ^{R^}	0.05 – 2.0 3.57-142.86	0.0029	mg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 520 nm multi-range method; NPDES Equivalent (353.2).	20-Nov-08
31-107-04-1-A ^{R^}	0.1 – 10 0.0071-0.713	0.0049	mg N/L mM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 520 nm. NPDES Equivalent (353.2)	2-May-08
	17.5 – 70 1.25-5.0	0.126	µg N/L µM N/L			

Method No	Range	MDL	Units	Matrix	Comments	Rev Date
31-107-04-1-C ^{RA}	0.07 – 0.70 5.0-50.0	0.00168	mg N/L µM N/L	Brackish/ Seawaters	Sulfanilamide/NED Cd reduction method; 520 nm.NPDES Equivalent (353.2)	2-May-08
31-107-04-1-D ^{RA}	0.5 – 14 0.036-1.0	0.2	µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 540 nm NPDES Equivalent (353.2). Requires a standard heater.	2-May-08
31-107-04-1-E ^{EA}	5 – 400 0.36-28.57	0.51	µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 540 nm NPDES Equivalent (353.2).	19-Aug-03
31-107-04-1-F ^{RA}	0.25 – 14 0.018-1.0	0.042	µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED. Cd reduction method; 540 nm. 2-cm detector method; QC8500 only ; NPDES Equivalent (353.2). Requires a standard heater.	8-Jul-08
31-107-04-1-G ^{RA}	0.25 – 10 0.018-0.714 0.01– 1.0 0.714-71.43	0.05 0.002	mg N/L mM N/L mg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 540 nm <u>Ultra High Throughput method</u> . (>120 samples/hr); multi-range; NPDES Equivalent (353.2)	24-Apr-08
31-107-04-1-H ^{RA}	0.25 – 30 0.18-2.143	0.025	mg N/L mM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 540 nm can also use w/ non-saline matrix; NPDES Equivalent (353.2)	28-Oct-08
31-107-04-1-I ^R	5 – 500 0.357-35.71	0.025	mg N/L mM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 540 nm. Ultra high level, inline dialysis method	12-Jul-09
31-107-04-1-J ^{RA}	1-100 0.071-7.143	0.2	µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cd reduction method; 520 nm. NPDES Equivalent (353.2)	30-Jun-10
31-107-04-1-K ^R	1-100 0.071-7.143	0.2	µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED, Cd reduction method. Imidazole buffer . 540 nm. NPDES Equivalent (353.2)	20-May-11
31-107-04-5-A ^R	0.01 – 5.0 1.43-357.14	0.009	mg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Nitrate Reductase method ; 540 nm. Enzymatic reagents must be purchased from NECi.	11-Feb-09
31-107-04-6-A ^R	0.05 – 5.0 0.0036-0.357 0.2-20.0 0.0143-1.43	0.006	mg N/L mM N/L mg N/L mM N/L	Brackish / Seawaters	UV Nitrate Reduction ; Sulfanilamide/NED 540 nm. Multi-range method PATENTED. Requires an in-line module with UV lamp.	19-Jun-09

Method No	Range	MDL	Units	Matrix	Comments	Rev Date
Nitrite						
<i>See also IC section</i>						
31-107-05-1-A ^{RA}	17.5 – 70 1.25-5.0	0.01	µg N/L as NO ₂ ⁻ µM N/L as NO₂⁻	Brackish / Seawaters	Nitrite only; Sulfanilamide/NED 540 nm. NPDES Equivalent (353.2)	13-May-08
31-107-05-1-B ^{RA}	0.1 – 15 0.007-1.07	0.01	mg N/L as NO ₂ ⁻ mM N/L as NO₂⁻	Brackish / Seawaters	Nitrite only; Sulfanilamide/NED 540 nm. NPDES Equivalent (353.2)	29-Oct-08
Nitrogen - Total Nitrogen						
31-107-04-3-A ^R	25 – 1000 1.79-71.43	4.90	µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cadmium reduction , 540 nm. Total N; alkaline persulfate inline digestion method; samples w/ particulates not suitable.	3-Feb-10
31-107-04-4-B ^R	0.02– 5.00 1.43-357.14 1.0– 40.0 0.071-2.86	0.0068 0.111	mg N/L µM N/L mg N/L mM N/L	Brackish / Seawaters	Sulfanilamide/NED Cadmium reduction , 520 nm. Total N; manual persulfate digestion w/ Cd reduction; low-flow method; total phosphorus can be measured from same digest (31-115-01-4-B) ; multi-range method.	16-Jun-08 21-Feb-12
31-107-04-4-C ^R	0.02– 5.00 1.43-357.14 1.0– 40.0 0.071-2.86 NO₂ + NO₃-N: 2.5-500 NO₂- N: 1-125	0.0068 0.111 0.44 0.2	mg N/L µM N/L mg N/L mM N/L µg N/L µM N/L µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cadmium reduction , 520 nm. Total N; manual dual persulfate digestion Imidazole buffer ; low-flow method; total phosphorus can be measured from same digest (31-115-01-4-B); Can also use for NO₂ + NO₃ and NO₂. (Support for NO ₂ + NO ₃ and NO ₂ included).	

Method No	Range	MDL	Units	Matrix	Comments	Rev Date
Orthophosphate						
31-115-01-1-G ^{RA}	62 – 310 2.0-10.0		µg P/L as PO ₄ ²⁻ µM P/L as PO₄²⁻	Brackish / Seawaters	Orthophosphate; molybdate based method; 880 nm. NPDES Equivalent (365.5). Requires a standard heater	13-May-08
31-115-01-1-H ^{RA}	5 – 400 0.16-12.9	1.0	µg P/L as PO ₄ ²⁻ µM P/L as PO₄²⁻	Brackish / Seawaters	Orthophosphate; molybdate based method; 880 nm. NPDES Equivalent (365.5). Requires a standard heater	13-May-08
31-115-01-1-I ^{RA}	1-100 0.032-3.23	0.25	µg P/L as PO ₄ ²⁻ µM P/L as PO₄²⁻	Brackish / Seawaters	Orthophosphate; molybdate based method; 880 nm. NPDES Equivalent (365.5) Requires a standard heater	13-May-08
31-115-01-1-J ^{RA}	0.01 – 2.0 0.323-64.52	0.002	mg P/L as PO ₄ ²⁻ µM P/L as PO₄²⁻	Brackish / Seawaters	Orthophosphate; molybdate based method; 880 nm. NPDES Equivalent (365.5)	30-Nov-07
	0.5-20 0.016-0.645	0.01	mg P/L as PO ₄ ²⁻ mM P/L as PO₄²⁻		Requires a standard heater	
31-115-01-1-W ^{RA}	0.25 – 20 0.008-0.645	0.007	µg P/L as PO ₄ ²⁻ µM P/L as PO₄²⁻	Brackish / Seawaters	Orthophosphate; molybdate based method; 880nm <u>2-cm detector method</u> ; QC8500 only; for samples with very low or no silicate ; NPDES Equivalent (365.5). Requires a non-standard Heater PN 58112 allows replicates from a single sample tube.	22-Feb-08
31-115-01-1-Y ^{R**A}	0.5 – 100 0.016-3.23	0.164	µg P/L as PO ₄ ²⁻ µM P/L as PO₄²⁻	Brackish / Seawaters	Orthophosphate; molybdate based method; 880nm 2-cm detector method; QC8500 only ; for samples with high silicate; NPDES Eq. (365.5) Requires a non-standard heater PN 58112 allows replicates from a single sample tube.	29-Feb-08
Phosphorus, Total (Acidic Persulfate)						
31-115-01-3-D ^R	0.050 – 1.0 1.63-32.36	0.002	mg P/L µM P/L	Brackish / Seawaters	Total P; molybdate method; 880 nm. inline persulfate digestion; samples w/ particulates not suitable. glass calibration vials. Requires an in-line sample prep module.	5-Jul-07


Method No	Range	MDL	Units	Matrix	Comments	Rev Date
31-115-01-3-F ^R	2-100 0.065-3.23	0.59	µg P/L µM P/L	Brackish / Seawaters	Total P; molybdate based method; 880nm, inline persulfate digestion; samples w/ particulates not suitable Requires glass standard and sample vials an in-line sample prep module, and non-standard heater.	13-Oct-08
31-115-01-4-A ^{RA}	12.5 – 400 0.40-12.9	1.66	µg P/L µM P/L	Brackish / Seawaters	Total P; molybdate based method; 880 nm. manual persulfate digestion; Requires a standard heater and autoclave for the digestion. NPDES Equivalent (365.3)	17-Sep-03
31-115-01-4-B ^R	0.005 – 1.0 0.16-32.23	0.0038	mg P/L µM P/L	Brackish / Seawaters	Total P; molybdate based method; 880nm, manual persulfate digestion; low-flow method; dual-Digest . Total N Can be measured simultaneously using 31-107-04-4-B; multi-range method Can also analyze particulate phosphorus and orthophosphorus with this method. Requires a standard heater and block digester for the digestion.	12-Dec-09
	0.25 – 10 0.008-0.323	0.0358	mg P/L mM P/L			
OP:	5 – 1000 0.16-32.23	0.7	µg P/L			
	0.25 – 10 0.008-0.323	0.013	mg P/L			
pP:	0.1-5.0	0.015	mg P/L			
Silicate						
31-114-27-1-A ^{RA}	1,202-6,009 20 – 100	0.2	µg SiO ₂ /L µM SiO₂/L	Brackish / Seawaters	Molybdate reactive method; Stannous chloride reductant. 820 nm. Requires a standard heater. NPDES Equivalent (USGS I-2700-85)	17-Sep-03
31-114-27-1-B ^{RA}	75.0-300.45 1.25 – 5.0	0.01	µg SiO ₂ /L µM SiO₂/L	Brackish/ Seawaters	Molybdate reactive method; Stannous chloride reductant. 820 nm. Requires a standard heater and plastic sample vials. NPDES Equivalent (USGS I-2700-85)	17-Sep-03
31-114-27-1-D ^{RA}	10 – 1700 0.166-28.29	1.43	µg SiO ₂ /L µM SiO₂/L	Brackish/ Seawaters	Molybdate reactive method; Stannous chloride reductant. 820 nm. Requires a standard heater. NPDES Equivalent (USGS I-2700-85)	17-Sep-03

Method No	Range	MDL	Units	Matrix	Comments	Rev Date
31-114-27-1-E ^{RA}	2.5 – 100 0.042-1.66	0.606	µg SiO ₂ /L µM SiO₂/L	Brackish / Seawaters	Molybdate reactive method; 820 nm, 2cm detector method; QC8500 only; Requires a standard heater and plastic sample and standard vials. NPDES Equivalent (USGS I-2700-85)	28-Feb-08
31-114-27-1-F ^{RA}	0.5 – 30 0.0083-0.499	0.05 0.00083	mg SiO ₂ /L mM SiO₂/L	Brackish / Seawaters	Molybdate reactive method; Stannous chloride reductant 820 nm, NPDES Equivalent. (USGS I-2700-85) . Requires a standard heater.	23-Oct-08
31-114-27-2-A ^R	60.09-6009 1-100	0.6 0.1	µg SiO ₂ /L µM SiO₂/L	Brackish / Seawaters	Molybdate reactive method; Ascorbic acid reductant 820 nm. NPDES Equivalent (366.0) Requires a standard heater and plastic sample and standard vials.	23-Nov-10


Urea

31-206-00-1-A ^R	10 – 400 0.714-28.57	2.9	µg N/L as Urea µM N/L as Urea	Brackish / Seawaters	Diacetyl monoxime/thiosemicarbazide method. 530 nm. Requires non-standard heater and 60 position sample racks.	16-Sep-03
31-206-00-1-B ^R	0.025 – 5.00 1.79-357.14	0.004	mg N/L as Urea µM N/L as Urea	Brackish / Seawaters	Multi-range method. Diacetyl monoxime/thiosemicarbazide. 530 nm. Requires non-standard heater and 60 position sample racks.	7-Dec-07
	0.2 – 20 0.0143-1.429	0.026	mg N/L mM N/L as Urea			


Selected Parts. Please note this is NOT a complete listing




50902 or 50902L
Tee Fitting




50045L
Reagent Line Weight




50046L
50036



50065 Flow cell flare, 13 cm




50067 Flow cell flare, 35 cm




31933 Flow cell, 10mm


Nipples




50015 Large




50015A Small




50014L Reducing




50019 Stainless Steel




Interference Filters
89### ## indicates the wavelength
400 nm



24927 Switching Valve
50962 Flare Kit for Nitrate
50963 Flare Kit for Sulfate
24012 Connector, Amber
24015 Washer




50237A Cadmium Column




50254 Cadmium Column Maintenance kit

Union Fitting




50913




50901

Small Pump Tube Adapter




50906




85287 for QC8500

Large Pump Tube Adapter




85293 for QC8500




50907


Pump Tube Adapter




50905




85269 for QC8500




50009
Tube Connector, Amber




31077 for QC8000 Valve
Valve Connector, White




85245
One-piece fitting
for QC8500 valve




28051L
Rheodyne Nut & Ferrule




28058 Rheodyne Ferrule




28057
Upchurch Nut & Ferrule




50013
O-ring, Viton




50024L
Large Clip, Gray




50023L
Small Clip, Gray




50913 Teflon
Union Fitting
(ultem nuts)




50007 or 85258 for QC8500
Large Collar, Black




50006 or 85257 for QC8500
Small Collar, White



50060
O-ring Probe



50100L
Degassing Tube



50031
Scalpel



28081
Blades for PEEK
Tubing Cutter




50021
Transmission Tubing, PVC, 0.030" ID



50029
Transmission Tubing, PVC, 0.060" ID




50091 Microloop Tubing, Teflon, 16 cm
50092 Microloop Tubing, Teflon, 12.5 cm



28080 PEEK Tubing Cutter



Teflon Tubing
50927 Manifold, 0.022" ID, Green
50928 Manifold, 0.032" ID
41300L Tefzel Tubing, 0.040" ID x 0.060" OD
30928 Zeus Tubing, 0.032" ID



Coils

Size	Coil Support Only	Wrapped, Teflon 0.022" ID	Wrapped, Teflon 0.032" ID
1" or 4.5 cm	50016L	50981	50916
2" or 7cm	50018L	50982	50918
2.5" or 8 cm	50017	50983	50917
4" or 12cm	50020	50984	50920
8" or 22 cm	50022L	N/A	50922



Pump Tubing

534XX	PVC
544XX	Duraprene
494XX	Silicone
654XX	Acidflex

XX is the number that specifies the pump tube color:

05	Orange-Yellow, 0.020" ID	13	Blue-Blue, 0.065" ID
06	Orange-White, 0.025" ID	14	Green-Green, 0.073" ID
07	Black-Black, 0.030" ID	15	Purple-Purple, 0.081" ID
08	Orange-Orange, 0.035" ID	16	Purple-Black, 0.090" ID
09	White-White, 0.040" ID	17	Purple-Orange, 0.100" ID
10	Red-Red, 0.045" ID	18	Purple-White, 0.110" ID
11	Gray-Gray, 0.051" ID	19	Yellow-Blue, 0.060" ID
12	Yellow-Yellow, 0.056" ID		




Alternating Coils 12 cm

Coil Support Only	Wrapped, Teflon 0.022" ID	Wrapped, Teflon 0.032" ID
50039	50985	50921



Bottles
28193L Glass Bottle, 1000 mL
35102 Glass Bottle, 2000 mL
43915 Glass Bottle, 100 mL (includes cap)



50012
O-ring Remover



28101L
End Plug



50062
Drill Bit

Lachat Instruments Brand Loveland, Colorado USA

United States: 800-247-7613 tel | 970-669-3050 tel | 970-461-3915 fax | sales@lachatinstruments.com | techhelp@hach.com

Outside United States: Contact the Lachat office or distributor serving you. | sales@lachatinstruments.com | intltech@hach.com

www.lachatinstruments.com

Printed in U.S.A. ©Hach Company, 2016. All rights reserved.

LACHAT
INSTRUMENTS

