DOC022.97.90756

QP 1680-TOC/TNb

11/2021, Edition 1

Basic User Manual
Manuel d'utilisation de base
Manual básico del usuario
Manual Básico do Usuário

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Section 1 Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x H x D)	38 x 44 x 70 mm (15.0 x 17.4 x 27.6 in.)
Weight	29 to 32 kg (64.0 to 70.6 lb)
Pollution degree	2
Overvoltage category	II
Parameters	TOC, TNb
Oxidation method	TOC: Non-dispersive infrared detection (NDIR); TNb: Chemiluminescence
Temperature	Furnace temperature 1050 °C (1922 °F) maximum
Measuring range	0 to 30,000 mg/L
Detection limit	TC, TIC, NPOC: 50 µg/L; TNb: 20 µg/L
Repeatability	Until 10 mg/L: < 5%; More than 10 mg/L: < 2%
Sample volume	10 to 1000 µL
Power requirements	Analyzer: 100-240 VAC ±10%, 50/60 Hz; 16 A (with safety ground); 750 W maximum PC: 100-240 VAC ±10%, 50/60 Hz; 16 A (with safety ground); 90 W maximum Monitor: 100-240 VAC ±10%, 50/60 Hz; 16 A (with safety ground); 100 W maximum
Operating temperature	20 to 30 °C (68 to 86 °F)
Storage temperature	-20 to 60 °C (-4 to 140 °F)
Relative humidity	20 to 80%, non-condensing
Altitude	2000 m (6562 ft) maximum
Environmental conditions	Indoor use
Gas supply	Type: Oxygen or synthetic air; Quality: 99,998% minimum (4.8); Pressure: 3 to 10 bar (43.5 to 145 PSI)
Gas connection	1/8-inch Swagelok male connector for 1/8-inch OD tubing ¹

¹ Supplied with the analyzer.

Specification	Details
Gas consumption	QP 1680-TOC: 150 mL/min ² QP 1680-TOC/TNb: 200 mL/min ² QP 1680-TNb: 150 mL/min
Certifications	CE, UKCA, certified to UL safety standards
Warranty	1 year (EU: 2 years)

Section 2 Online user manual

This Basic User Manual contains less information than the User Manual, which is available on the manufacturer's website.

Section 3 General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

3.1 Safety information

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

3.1.1 Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

² The analyzer uses 300 mL/min more to remove NPOC.

3.1.2 Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol indicates that the marked item can be hot and should not be touched without care.
	This symbol identifies the presence of a strong corrosive or other hazardous substance and a risk of chemical harm. Only individuals qualified and trained to work with chemicals should handle chemicals or perform maintenance on chemical delivery systems associated with the equipment.
	This symbol indicates that the object is heavy.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

3.2 Intended use


The QP 1680-TOC/TN_b analyzer is intended for use by individuals who measure water quality parameters in the laboratory.

3.3 Product overview

The QP 1680-TOC/TN_b analyzer measures non-purgeable organic carbon (NPOC) and total bound nitrogen (TN_b) in aqueous samples. The analyzer also measures total carbon (TC), total inorganic carbon (TIC) and calculates total organic carbon (TOC).

The analyzer has an integrated sampler that analyzes carbon and nitrogen in different water matrices. The analyzer connects to a PC with the TE Instruments Software (TEIS2) for operation and automation. Refer to [Figure 1](#) and [Figure 2](#).


Figure 1 Product overview—Front and side views

1 QP 1680-TOC/TN analyzer	9 Safety cover
2 Furnace compartment	10 Stirrer and sampler
3 Front cover	11 LED indicator
4 Filter holder	12 USB port type B
5 Particle filter, 5 µm	13 Power supply connection and power switch
6 IC vessel	14 1/8" swagelock oxygen port
7 Humidifier	15 Waste connection
8 Halogen scrubber	16 Fan

LED color	Description
Green	The analyzer is prepared for use.
Yellow	The analyzer is not prepared for use. Some parameters are not within the measurement range.
Blue	The analyzer is in operation.
Red	The analyzer has no connection to the TEIS2 software.

Figure 2 Product overview—Top view


1 Furnace injection point	5 Reagent container (ultra-pure water)
2 Waste port	6 Reagent container (acid)
3 IC injection point	7 Optional waste container
4 Sampler arm	8 Vial tray, 65 or 96 samples

3.4 Product components

Make sure that all components have been received. Refer to [Figure 3](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 3 Product components

1 Analyzer	8 Sample tray
2 Sample introduction needle, stainless steel	9 USB cable
3 Combustion tube	10 24-mL sample vials, glass (100x)
4 Filter holder	11 1-L containers (3x)
5 Humidifier assembly	12 USB flash drive with TEIS2 software
6 Halogen scrubber	13 PTFE blue-stripe tubing, 1/8-inch OD, 2 m (includes the nut and ferule)
7 IC vessel	14 Power cord

Section 4 Installation

▲ DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

4.1 Installation guidelines

▲ WARNING

Fire hazard. This product is not designed for use with flammable liquids.

▲ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

▲ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

NOTICE

This instrument is sensitive to electromagnetic and electromechanical interference. These interferences may have an effect on the analysis performance of this instrument. Do not put this instrument near equipment that may cause interference.

Install the instrument:

- Indoors in a clean, dry, well-ventilated and temperature-controlled location
- In a location with a flat surface to prevent the unit from movement
- In a location with minimum mechanical vibrations and electronic noise
- In a location with no electromagnetic interferences from equipment such as transmitters, power switching, etc.
- In an environmental enclosure that supplies protection from precipitation and direct sunlight
- In a location where there is sufficient clearance to make plumbing and electrical connections
- In a location where the display is easily seen by the user
- In a location with a compatible waste container for the drain tube

4.2 Prepare the reagents and standards

▲ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

- Wear protective clothing, safety glasses or face protections and rubber gloves.
- Prepare the reagents in a fume hood.
- Use only glass or PTFE labware.
- Make sure that all of the bottles vent to air after installation.
- Make sure that the applicable accident prevention regulations are observed.
- Dispose of substances correctly and in accordance with the applicable regulations.

For NPOC and IC analysis, a solution of strong acid (e.g., HCl 3 mol/L) is added with a capacity to transfer 300 ppm of IC from the sample vial (40 mL) into CO₂. Use 500 mL of acid solution to prepare approximately 750 NPOC samples. Adjust the acid concentration or amount of acid added based on the sample pH or buffering strength to get a sample pH less than 2. Make sure that the sample pH is less than 2 before the analysis.

Table 1 Reagents and standards

Chemical	Specification
Ultra-pure water (UPW)	Electrical conductivity at 25°C: 0.0555 µS/cm maximum; electrical resistivity at 25°C: 18 MΩ·cm minimum <i>Note: UPW is used to prepare standards or dilutions. 2.5 L of UPW is necessary for the analyzer installation.</i>
Phosphoric acid (H ₃ PO ₄) 85% w/w ³	Used to prepare the 9.33% phosphoric acid solution (1 mol/L) <i>Note: The recommended acid solution is phosphoric acid (H₃PO₄). As an alternative, other strong acids can be used.</i>
Phosphoric acid solution (H ₃ PO ₄) (1 mol/L) 9.33% w/w	54.9 mL is used to prepare 500 mL of solution.
Hydrochloric acid (HCl) 36% w/w	Used to prepare the 10.4% hydrochloric acid solution.
Hydrochloric acid solution (HCl) (3 mol/L) 10,4% w/w	144.7 mL is used to prepare 500 mL of solution.
Sulfuric acid (H ₂ SO ₄) 96% w/w	Used to prepare the 13.5% sulfuric acid solution
Sulfuric acid solution (H ₂ SO ₄) (1,5 mol/L) 13,5% w/w	70.3 mL is used to prepare 500 mL of solution.

4.3 Remove the transport locks

▲ WARNING

Personal injury hazard. Instruments or components are heavy. Use assistance to install or move.

NOTICE


Use the bottom of the analyzer to lift the analyzer. Do not use the plastic covers or the sampler arm to lift the analyzer or damage to the analyzer can occur.

Remove the magnetic lock and transportation strips from the analyzer. Refer to [Figure 4](#) and the steps that follow:

1. Remove the white stirrer lock.
2. Remove the white strip.
3. Remove the red warning strip..
4. Discard the red warning strip
5. Install the white strip again.

³ Concentration of the solution in terms of percentage composition. x% w/w solution equals X grams of solute dissolved in 100 grams of solution.

Figure 4 Remove the transport locks

4.4 Connect to carrier gas

▲ DANGER

Fire and explosion hazard. Oil and grease can burn with explosive violence when there is oxygen. Make sure that all parts that are in contact with oxygen do not have grease or oil.

The instrument cannot operate without a carrier gas connection. An oxygen or synthetic air supply must be within 2 m (6.5 ft) of the rear of the analyzer. Use the supplied gas inlet tubing to connect oxygen or synthetic air. Refer to [Specifications](#) on page 3 for the gas supply specifications.

Connect the analyzer to a carrier gas as follows:

1. Open the gas air supply for some seconds to remove unwanted material.
2. Use a wrench to connect the supplied PTFE blue-stripe tubing to the gas supply and the Oxygen port. Refer to [Figure 5](#).
3. Do not open the gas supply. Make sure that the carrier gas connection is leak tight.

Figure 5 Connect to carrier gas

4.5 Connect the waste tubing

The analyzer has a waste connector to discard the water after the analysis. Refer to [Figure 6](#) and the steps that follow.

1. Use 4-mm OD tubing to connect the waste. The waste connector is a push-in connector. Push the tubing firmly into the waste connector.
2. Put the other end of the tubing in a sink or a waste container within 3 m (9.8 ft) of the analyzer.

Figure 6 Connect the waste tubing


4.6 Assemble the combustion tube

Do not start the analyzer without the combustion tube installed. Refer to [Figure 7](#) and the steps that follow.

1. Remove the top cover of the furnace.
2. Put the combustion tube in the furnace.
3. Use the thumbscrews to attach the furnace tube.
4. Do not connect the gas line to the combustion tube.

Figure 7 Assemble the combustion tube

4.7 Assemble the sampler arm

1. Use the supplied 4-mm allen screw to attach the head of the sampler arm to the sampler.
2. Put the sample introduction needle into the guide of the sampler arm.

4.8 Connect the PC and power

⚠ DANGER	
	Electrical shock and fire hazards. Make sure that the supplied cord and non-locking plug meet the applicable country code requirements.
⚠ DANGER	
	Electrocution hazard. Protective Earth Ground (PE) connection is required.

⚠ DANGER

Electrocution and fire hazard. Attach the power cord so that it cannot touch hot surfaces.


A PC with TEIS2 software installed is necessary to operate the analyzer.

1. Put the PC near to the analyzer. Make all of the necessary PC connections.
2. Use the supplied USB cable to connect the PC to the analyzer. Refer to [Figure 8](#).

Note: Make sure that the USB cable is not parallel to a power cable or electrical noise interference on the connection can occur.

3. Make sure that the power switch of the analyzer is set to OFF. Connect the supplied power cord to a grounded electrical outlet.

Figure 8 Connect the PC and power

4.9 Install the software

Note: Do not start the TEIS2 software at this moment. Not loading the correct database will cause the instruments to malfunction.

Use the supplied USB flash drive to install the TEIS2 software on the PC as follows.

1. Connect the USB flash drive to the PC.
2. Go to USB > Software > 3rd part. Install two software programs:
 - a. **AbdeRdr90_en-Us.exe**
 - b. **vcredist_x86.exe** (internet version) or **dotNetFx40_Full_x86_x64.exe** (off-line version)
3. Go to USB > Software > TEIS > V2.x.x.
4. Install **Setup_TEIS2.msi**.
5. Open the TEIS2 software. "No database found" shows.

6. Click OK.
7. Select the database from the USB flash drive. Go to USB > Database > Backup 20XX.XXX > TEIS_DATA.FDB
8. Close the software and restart the PC.


Section 5 Startup

Items to collect: flow meter


The first time the analyzer is used, do the steps that follow:

1. Set the PC to on.
2. Set the power switch of the analyzer to on.
3. Open the gas supply. Make sure that the gas pressure is less than 8 bar (116 PSI).
4. Start the TEIS2 software. Login with service engineer credentials.
5. Wait until all of the devices in the software show as connected.
6. Use the flow meter to examine the gas connections. Then, install the correct item.
 - a. Connect a flow meter to the lower right connection of the humidifier. The flow must be $150 \text{ mL/min} \pm 5 \%$.
 - b. Fill the humidifier with ultra-pure water. Carefully connect the humidifier.

Note: Make sure to push the humidifier in a straight line into the connector or damage to the connector can occur.

- c. Connect the flow meter to the carrier gas tubing in the top of the furnace. The flow must be $150 \text{ mL/min} \pm 5 \%$.
- d. Connect the carrier gas tubing at the top of the furnace. Install the top cover of the furnace.

- e. Connect the flow meter to the middle connection of the IC vessel.


The flow must be 150 mL/min \pm 5 %.

- f. Carefully install the IC vessel. Apply pressure in the middle of the IC vessel.


- g. Connect the flow meter to the bottom connection of the halogen scrubber. The flow must be 150 mL/min \pm 5 %.

- h. Carefully install the halogen scrubber. Apply pressure in the middle of the halogen scrubber.

- i. Connect the flow meter to the bottom connection of the filter. The flow must be 150 mL/min \pm 5 %.

- j. Install the filter holder with the particle filter.

7. In the TEIS2 software, open **Visual devices**.

The **Total flow** box shows the total flow in mL/min. For analyzers without nitrogen detector, the flow must be 150 mL/min \pm 5 %. For analyzers with nitrogen detector, the flow must be 200 mL/min \pm 5 %.

8. Start the method "Test method: Send sampler to IC port to check system pressure". The analyzer sends the Sampler to the IC port to close the IC port inlet.
The **Condenser** box shows the temperature of the condenser. When set at 1 °C, the value must be 1 °C ± 0.5 °C.
The **Oxygen pressure** box shows the internal pressure of the oxygen or synthetic air. The value must be 1800 to 2300 mbar.
9. Remove the containers from the analyzer.
In **Visual devices**, make sure that the values for the weigh unit cups are 0 g ± 10 g.
10. Install the containers filled with the applicable reagents in the weight units. Use the color label to identify the tubing for each container.
 - Weight unit 1: ultra pure water. Connect the two green tubes.
 - Weight unit 2: acid. Connect the red tube.
 - Weight unit 3: waste. Use 4-mm PTFE tubing for the waste. As an alternative, put the tube in a sink.
11. In **Visual devices**, examine the cooling and heating values.
For analyzers with nitrogen detector, the cooling value must be 10 °C ± 5 %. The heater value must be 50 °C ± 5 %.
12. Examine the furnace temperature values.
For analyzers without nitrogen detector, the value must be 680 °C ± 5 %. For analyzers with nitrogen detector, the value must be 720 °C ± 5 %.
13. Install the tray in the analyzer. Make sure that the tray shows in **Visual devices**.
14. Open and close the front cover. Make sure that the status of the front cover shows in **Visual devices**.

5.1 Examine the sampler tubing

Make sure to fix the sampler tubing with the tubing clamps. Make sure that the routing of the sampler tubing is correct. The tubing must let the sampler arm move freely.

Section 6 Maintenance

▲ WARNING

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

▲ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

▲ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

6.1 Maintenance schedule

▲ WARNING

Multiple hazards. Only qualified personnel must repair and do maintenance tasks of compressed gas cylinders and accessories.

Table [Table 2](#) shows the recommended schedule of maintenance tasks. Facility requirements and operating conditions may increase the frequency of some tasks.

Table 2 Maintenance schedule

Task	1 day	1 week	30 days	90 days	365 days
Examine for leaks and spills. Clean as necessary.	X				
Make sure that the tubing in the reagent containers touches the bottom of the container.		X			
Make sure that the waste container has sufficient capacity. Discard waste as necessary.		X			
Examine the level in the humidifier. Fill as necessary.		X			
Examine the injection needle for damage. Replace if damaged.		X			
Examine the movement of the sampler		X			
Clean the furnace injection port.		X			
Replace all UPW in the UPW reagent container.			X		
Replace the particle filter.			X		
Examine all tubing from the reagent containers for damage.			X		
Examine for dust at the air intake points of the analyzer. Replace the fan filter if necessary.			X		
Examine for leaks at the reagent containers. Replace the containers as necessary.			X		
Examine the O-ring of the furnace tube (red, top of the furnace tube)				X	
Examine the O-ring of the furnace injection port (black)				X	
Examine the halogen scrubber. Replace as necessary.				X	
Replace the catalyst.				X	
Annual maintenance ⁴					X

⁴ An annual maintenance inspection by a factory-trained and certified service engineer is recommended. The annual maintenance includes important checks and part replacements to keep the analyzer in good condition.

6.2 Clean spills

▲ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

1. Obey all facility safety protocols for spill control.
2. Discard the waste according to applicable regulations.

6.3 Stop the analyzer

Do the steps that follow to stop the analyzer for short periods (4 hours or more).

1. In the TEIS2 software, go to **Task manager** and select **System method**.
2. Select **Standby method**.
3. Click **Add system method**. The method shows in the **Task manager**.
4. Select **hold** for all other tasks in the task manager. Click **Start**.
5. Wait until the analyzer completes the Standby procedure.

6.4 Shutdown the analyzer

NOTICE

Do not set the analyzer to off until the shutdown procedure is complete and the furnace temperature has decreased to room temperature or damage to the analyzer will occur.

Do the steps that follow to prepare the analyzer to be stopped for a long period (more than 7 days):

1. In the TEIS2 software, go to **Task manager** and select **System method**.
2. Select **Shutdown method** from the list.
3. Click **Add system method**. The method shows in the **Task manager**.
4. Select **hold** for all other tasks in the task manager. Click **Start**.
5. Wait until the analyzer completes the shutdown procedure.
6. Set the analyzer to off.

Section 7 Troubleshooting

▲ WARNING

Burn hazard. The furnace gets extremely hot. Do not touch. The outer surface of the analyzer may become hot. Keep chemicals away from hot surfaces.

Make sure to set the analyzer to off. Let the analyzer temperature decrease to room temperature before a disassembling operation.

Problem	Possible cause	Solution
Leakage or loss of flow	There is a problem in the furnace injection point.	Add a drop of UPW on top of the furnace injection point. If there are air bubbles, remove and fully clean the injection port. As an alternative, replace the injection port.
	There is a problem in IC injection point.	<ul style="list-style-type: none"> • Add a drop of UPW on top of the IC injection point. If there are air bubbles, use a spray bottle to put some UPW into the IC injection point. • If the IC injection point continues to push the UPW out, the condenser has a possible blockage. Set the condenser temperature to 20 °C (68 °F). Wait 30 minutes. Use the spray bottle to put some UPW into the IC injection point. The flow should come back to the standard value. Contact technical support.
	There is a problem in particle filter.	The analysis can form CuCl, which is a tiny powder that can cause a blockage in the filter. The blockage can increase the pressure, which can cause the measured flow to be lower than expected. Use compressed air to clean the particle filter. Replace the particle filter.
	There is a leak after the furnace tube is assembled.	The O-ring that seals the bottom of the furnace tube can expand or shrink during furnace operation. Remove the complete furnace and use a flathead screwdriver to tighten the O-ring.
No air bubbles in the reagent bottle UPW	The system startup does not work.	Make sure that the system starts up correctly. If the startup is unsuccessful, the oxygen valves can be damaged. Contact technical support to replace the oxygen valves.
Area counts too high	There is a problem in the needle or needle sleeve.	<ul style="list-style-type: none"> • Examine the needle and the needle sleeve for damage. • Make sure that the sleeve starts at 21 mm from the bottom of the needle. Replace the needle if necessary.
	There is a problem in the particle filter or a leak in the system. The pH value of the acid is incorrect.	<ul style="list-style-type: none"> • Examine the particle filter. Replace the particle filter if necessary. • Make sure that the pH value of the acid is less than 2. • Compare the current flow of the instrument to the flow of the samples in the calibration line used to calculate the current samples. The flow has an effect on the area measured from a sample. A lower flow gives a higher area count and vice versa.

Problem	Possible cause	Solution
High area counts and large RSD for blanks	The tubes of the TN-detector are defective.	The tube has some tarnish. Contact technical support.
NPOC results too high	In TC measurements, the high temperature can cause IC collection in the furnace. When NPOC is measured after TC, the collected IC will break free and have an effect on the area counts.	Make sure to measure two NPOC blanks to prevent defective NPOC measurements.
Double peaks for TN measurement	The ceramic wool layer in the top of the furnace tube is defective.	Replace the furnace tube.
High and low area counts within a sample	There is a problem in the O-ring of the furnace exit.	Replace the O-ring and analyze 15 blanks.
	There is a problem in the furnace tube.	A very small tear in the furnace tube can cause a leak. Replace the furnace tube and analyze 15 blanks.
	The halogen scrubber is not assembled correctly.	Remove the halogen scrubber and install it again.
Halogen scrubber should be replaced within days	The tubing of the condenser is not assembled correctly. There is liquid in the detector chambers.	Examine the tubing. Identify if the particle filter is moist. If the filter is wet, identify if the CO ₂ detector is moist. Clean the CO ₂ detector as necessary.
System shows "Leakage below IC/Waste port"	There is a problem in the IC/waste port.	<ul style="list-style-type: none"> • Replace the sensor. • Remove the sensor below the leakage container. Rinse the sensor with ethanol and then with ultra pure water. Remove all possible liquids from the sensor. Install the sensor again. <p><i>Note: Make sure to install the sensor below the threshold before assembling the analyzer.</i></p>
Memory effect	The furnace tube is dirty.	Remove the furnace tube. Wait until the furnace tube temperature decreases to ambient temperature. Rinse the quartz and the catalyst separately with ultra pure water. Assemble and install the furnace tube in the analyzer. Startup the furnace heat.
IC vessel filled with liquid	The IC pump is defective.	Contact technical support.

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)