
DOC022.97.90250

sensION™ + MM340

06/2020, Edition 4

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Table of Contents

- | | |
|---|---|
| 1 Specifications on page 3 | 7 Advanced operation on page 11 |
| 2 General information on page 3 | 8 Data logger on page 13 |
| 3 Installation on page 5 | 9 Maintenance on page 14 |
| 4 User interface and navigation on page 8 | 10 Troubleshooting on page 16 |
| 5 Startup on page 8 | 11 Replacement parts and accessories on page 17 |
| 6 Standard operation on page 9 | A Standard solutions on page 18 |

Section 1 Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions	35 x 20 x 11 cm (13.78 x 7.87 x 4.33 in.)
Weight	1100 g (2.43 lb)
Meter enclosure	IP42, indoor use
Power requirements (external)	100–240 V, 0.4 A, 47-63 Hz
Meter protection class	Class III
Installation category	Category I
Altitude requirements	Standard 2000 m (6562 ft) ASL (Above Sea Level)
Storage temperature	–15 to +65 °C (5 to +149 °F)
Operating temperature	0 to 40 °C (41 to 104 °F)
Operating humidity	< 80% (non-condensing)
Resolution	pH: 0.1/0.01/0.001, ORP: 0.1/1 mV, ISE: programmable, temperature: 0.1 °C (0.18 °F)
Measuring error (± 1 digit)	pH: ≤ 0.005, ORP: ≤ 0.2 mV, temperature: ≤ 0.2 °C (≤ 0.36 °F)

Specification	Details
Reproducibility (± 1 digit)	pH: ± 0.001, ORP: ± 0.1 mV, temperature: ± 0.1 °C (± 0.18 °F)
Data storage	330 results and last 9 calibrations
Connections	2 Combined or indicator probe: BNC connector (Imp. >10 ¹² Ω); 2 Reference electrodes: banana connector; A.T.C. type Pt 1000: banana or telephonic; 2 magnetic stirrers: RCA connector
Temperature correction	Manual, Pt 1000 temperature probe (A.T.C.), isopotential pH programmable, standard value 7.00
Measurement display lock	Continuous measurement, by stability
Display	Liquid crystal, backlit, 128 x 64 dots
Keyboard	PET with protective treatment
Certification	CE

Section 2 General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

2.1 Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

2.1.1 Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

2.1.2 Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.

Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

2.2 Certification

▲ CAUTION

This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.

Canadian Radio Interference-Causing Equipment Regulation, ICES-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

2.3 Product overview


The sensION™+ meters are used with probes to measure various parameters in water.

The sensION™+ MM340 meter has two measurement channels to measure pH, ORP (mV) or ISE (concentration) with a selective probe. Measurement data can be stored and transferred to a printer or PC.

2.4 Product components

Refer to [Figure 1](#) to make sure that all components have been received. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 1 Meter components





1 Electrolyte for the probe	5 Probe holder
2 Buffer solutions (pH 4.01, pH 7.00 and pH 10.01)	6 Rod with o-ring
3 Calibration beakers (with magnetic bar inside)	7 Probe (included with kits only)
4 Power supply	8 Meter

Section 3 Installation

3.1 Assemble the probe holder

Follow the numbered steps to assemble the probe holder and to connect the magnetic stirrer.


3.2 Connect to AC power

The meter can be powered by AC power with the universal power adapter.


1. Select the correct adapter plug for the power outlet from the adapter kit.
2. Connect the universal power adapter to the meter ([Figure 2](#)).
3. Connect the universal power adapter to an AC receptacle ([Figure 3](#)).
4. Turn the meter on.

Figure 2 Connector panel

1	Magnetic stirrer 1 connector, Channel 1	7	PC keyboard, mini DIN connector
2	Reference electrode (separated electrodes) connector, Channel 2	8	Temperature probe connector, Channel 1
3	Reference electrode (separated electrodes) connector, Channel 1	9	Combined pH electrode (or indicator) connector, Channel 1
4	Separated temperature probe connector, Channel 1	10	Combined pH electrode (or indicator) connector, Channel 2
5	RS-232 for printer or PC connector	11	Magnetic stirrer 2 connector, Channel 2
6	Power supply		

Figure 3 AC power connection


Note: Position the equipment so it is not difficult to disconnect power to the equipment.

Section 4 User interface and navigation

4.1 User interface

Keypad description


1 RETURN key: cancel or exit the current menu screen to the previous menu screen	5 DOWN key: scroll to other options, change a value
2 MEASUREMENT key: confirm the selected option	6 ON/OFF: turn on or turn off the meter
3 UP key: scroll to other options, change a value	7 LEFT key: change between Channel 1 and Channel 2, enter numbers and letters
4 RIGHT key: change between Channel 1 and Channel 2, enter numbers and letters	

4.2 Display description

The meter display shows the concentration, units, temperature, calibration status, operator ID, sample ID, date and time.

Figure 4 Single screen display

1 Sample ID	4 Measurement channel
2 Measurement unit and value (pH, ORP (mV) or ISE)	5 Sample temperature (°C or °F)
3 Measurement mode or time and date	6 Visual measurement timer

4.3 Navigation

Use the to return to the previous menu. Use the measure key to take a sample measurement or to confirm options. Use the arrow keys to scroll to other options or to change a value. To change the parameters use the arrow keys and . Refer to each task for specific instructions.

Section 5 Startup

5.1 Turn the meter on and off

NOTICE

Make sure that the probe is connected to the meter before the meter is turned on.

Push to turn on or turn off the meter. If the meter does not turn on, make sure that the AC power supply is properly connected to an electrical outlet.

5.2 Change the language

The display language is selected when the meter is powered on for the first time.

1. Use the ▲ or ▼ to select a language from a list.
2. Confirm with ✓. The measurement screen shows DATA OUTPUT.
3. Select Deactivated if no printer or PC is connected and confirm. Refer to [Select the data output](#) on page 11 for more information about Data Output.

Section 6 Standard operation

6.1 Calibration

▲ DANGER	
	<p>Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.</p>

6.1.1 Calibration settings

The calibration settings contain Calibration type, Calibration frequency and Display options.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select the following options:

Option	Description
Stability C.:	Criterion by stability—select Fast, Standard or Strict.
Calibration type	Calibration type—select Technical buffers, DIN19266 Buffers, User Buffers, To a X value, data introduction or Theoretical calibration. Refer to Calibration types on page 9 for more information.

Option	Description
Cal. frequency	Calibration reminder—can be set between 0–7 days (default daily). The display shows the remaining time to the new calibration. Refer to Set the calibration reminder on page 10 for more information.
Display mV	Display mV—select YES or NO to show mV.
Standard 220 mV	Standard 220 mV—calibrate with a 220 mV standard solution
To a X value	To a X value—use a specific standard solution. Adjust the value during the calibration.
Data introduction	Data introduction Compensation—enter the compensation value
Factory adjust	Factory adjust—enter the compensation value

6.1.1.1 Calibration types

Different calibration types can be selected.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select Calibration type.

Option	Description
Technical buffers	pH 2.00, 4.01, 7.00, 9.21 and 10.01 at 25 °C (77 °F)
DIN19266 Buffers	pH 1.679, 4.006, 6.865, 7.000, 9.180, 10.012 and 12.454
User Buffers	Selected when the technical or DIN19266 buffers are not used. Refer to Technical buffer solutions (DIN 19267) on page 18 for pH values of specific buffer sets at varying temperatures.
Calibration to a X value	To adjust manually any scale value of the measured pH.

Option	Description
Data introduction	Manual probe constant introduction.
Theoretical calibration	The probe calibration data is replaced at 25 °C (77 °F).

6.1.2 Calibration procedure

This procedure is for general use with liquid calibration solutions. Refer to the documents that are included with each probe for additional information.

Note: Solutions should be stirred during calibration. For more information about the stirring settings, refer to [Change the stirring settings](#) on page 12.

1. Pour the buffers or calibration solutions into the labeled calibration beakers.
2. From the main menu use the ▲ or ▼ and ◀ and ▶ to select the CALIBRATION parameter. Confirm.
3. If required select the Operator ID (1 to 10) and confirm.
4. Rinse the probe with deionized water and put the probe into the first calibration beaker. Be sure that there are no air bubbles in the membrane.
5. Push ✓ to Start calibration.
6. Push ✓ to measure the first calibration solution. The next calibration solution is shown.
7. Rinse the probe with deionized water and put the probe into the second calibration beaker. Be sure that there are no air bubbles in the membrane.
8. Push ✓ to measure the second calibration solution. The next calibration solution is shown.
9. Rinse the probe with deionized water and put the probe into the third calibration beaker. Be sure that there are no air bubbles in the membrane.
10. Push ✓ to measure the third calibration solution. When the calibration is good, the display briefly shows Calibration OK and then returns to the main menu.

Note: When a printer is connected the print menu opens and the result can be printed.

6.1.3 View the calibration data

Data from the most recent calibration can be shown.

1. From the main menu use the ▲ or ▼ to select DATA LOGGER. Confirm.
2. Select Display data.
3. Select Calibration data and confirm with ✓. The last calibration data is shown.
 - pH-the slope and offset values are shown alternating with the deviation (in %) and calibration temperature.
 - ORP-the measured mV value and calibration temperature are shown.
 - Conductivity-the cell constant and calibration temperature for each standard are shown.

6.1.4 Set the calibration reminder

The calibration reminder can be set between 0 to 23 hours or 1-7 days (default 1 day). The display shows the remaining time to the new calibration.

Note: When 0 days is selected, the calibration reminder is turned off.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select Cal. frequency and confirm.
4. Use the ◀ and ▶ to advance to the next step and use the ▲ or ▼ to change a value. Confirm. Push ✓ to start the calibration.

6.2 Sample measurements

Each probe has specific preparation steps and procedures for taking sample measurements.

1. From the main menu use the ▲ or ▼ and ◀ and ▶ to select MEASURE. Confirm.
2. Use the ▲ to change the following settings. Confirm every entry.

Option	Description
Resolution	Select the resolution: 1, 0.1, 0.01 (default) or 0.001
Measure	Stability—select By stability Criterion: Fast (variation < 0.02 pH in 6 s), Standard (variation < 0.01 pH in 6 s) or Strict (variation < 0.002 pH in 6 s). In continuous—enter the time interval for the In continuous Acquis. interval (data storage or printing data). By time—enter the time interval for data storage or printing data.
Display mV	Display mV—select YES or NO to show mV.
Limits	Limits—select YES or NO. YES: Enter the upper and lower limits. An acoustic warning appears when the measurement is out of limit. The report output shows an A next to the measured value when the measurement was out of limit.
Isopotential	Isopotential—change the Isopotential pH value in Data introduction. Select Calculate to calibrate the probe again.

3. Push ✓ to start the measurement.

Note: If the measurement is not stabilizing after 120 seconds, the meter turns automatically into the continuous measurement mode.

Section 7 Advanced operation

7.1 Change the measurement units

The measurement units can be changed individually for each channel.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Measurement units and confirm.
3. Select Channel 1 or Channel 2 and confirm.
4. Select ORP (mV), pH or ISE and confirm.

7.2 Use a sample ID

The sample ID tag is used to associate readings with a particular sample location. If assigned, stored data will include this ID.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Sample ID and confirm.
3. Use the ▲ or ▼ to select

Option	Description
Automatic	A consecutive number will be automatically assigned to every sample.
Manual	A keyboard or a barcode scanner is required to enter the sample ID name before taking a measurement (maximum 15 characters).

7.3 Select the data output

Data can be stored or transferred to a printer or to a PC.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Data Output and confirm.
3. Use the ▲ or ▼ to select

Option	Description
Deactivated	Select Deactivated if no printer or PC is connected.
For Printer	Select Dot matrix printer or Thermal printer.
For Computer	Select Terminal, LabCom or LabCom Easy. The LabCom Software controls several modules, pH and conductivity meters, automatic burettes, Samplers and so on from a computer. The LabCom Easy software gets measurements and calibration data from the meters

7.4 Change the date and time

The date and time can be changed from the Date / Time menu.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Date / Time and confirm.
3. Use the ◀ and ▶ to advance to the next step and use the ▲ or ▼ to change a value. Confirm.
The current date and time will be shown on the display.

7.5 Adjust the display contrast

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Display contrast and confirm.
3. Use the ◀ and ▶ to adjust the contrast of the display and confirm.

7.6 Adjust the temperature

The temperature measurement can be adjusted at 25 °C (77 °F) and/or 85 °C (185 °F) to increase accuracy.

1. Put the probe and a reference thermometer in a container of water at approximately 25 °C and allow the temperature to stabilize.
2. Compare the temperature read by the meter with that of the reference thermometer. The difference is the adjustment value for the meter.
Example: reference thermometer: 24.5 °C; meter: 24.3 °C.
Adjustment value: 0.2 °C.
3. Enter the adjustment value for the 25 °C reading:
 - a. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
 - b. Use the ▲ or ▼ to select Readjust temp. and confirm.
 - c. Select Channel 1 or Channel 2 and confirm. The temperature for channel 1 is measured by the conductivity cell and channel 2 is measured by the pH probe. If no Conductivity cell is connected, the temperature value must be entered manually or the temperature measured in channel 2 can be applied to the meter.
 - d. Use the ▲ or ▼ to select 25 °C and confirm.
 - e. Use the arrow keys to enter the adjustment value for 25 °C. Confirm.

4. Put the probe and a reference thermometer in a container of water at approximately 85 °C and allow the temperature to stabilize.
5. Compare the temperature from the meter with that of the reference thermometer. The difference is the adjustment value for the meter.
 - a. Use the ▲ or ▼ to select 85 °C and confirm.
 - b. Use the arrow keys to enter the adjustment value for 85 °C. Confirm.
 - c. Select Save changes and confirm.

7.7 Change the stirring settings

The stirring speed can be changed during calibration and during a measurement.

1. Use the ▲ or ▼ to change the stirring speed during calibration and during a measurement.

7.7.1 Turn the stirrer on or off

The stirrer 1 works with Channel 1 and 2 (Stirrer 1). A second stirrer can be connected to channel 2 (Stirrer 2). For the stirrer 2 activation refer to the following steps.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm with ✓.
2. Use the ▲ or ▼ to select Stirrer N.2 and confirm with ✓.
3. Use the ▲ or ▼ to select YES to turn the stirrer 2 on.
Note: Select NO to turn the stirrer 2 off.

7.8 Change the temperature units

The temperature units can be changed to Celsius or Fahrenheit.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Temperature units and confirm.
3. Use the ▲ or ▼ to select between Celsius or Fahrenheit and confirm.

Section 8 Data logger

8.1 Display data

The Display data log contains Measurement data, Electrode report and Calibration data. The stored data can be sent to a printer or to a PC. When the data log becomes full (400 data points), the oldest data point is deleted when a new data point is added.

1. Use the ◀ and ▶ to select Channel 1 or Channel 2 and confirm.
2. From the main menu use the ▲ or ▼ to select DATA LOGGER and confirm.
3. Use the ▲ or ▼ to select Display data and confirm.
4. Use the ▲ or ▼ to select

Option	Description
Measurement data	Measurement data—stores automatically each time a sample is measured
Electrode report	Electrode report—stores automatically the electrode history and measurement conditions
Calibration data	Calibration data—stores automatically the current calibration

8.2 Delete data

The entire measurement data or electrode report log can be deleted to remove data that has already been sent to a printer or PC.

1. From the main menu use the ▲ or ▼ to select DATA LOGGER and confirm.
2. Use the ▲ or ▼ to select Erase and confirm.
3. Use the ▲ or ▼ to select Measurement data or Electrode report and confirm. Confirm again to delete the data. The entire log is deleted at once.

8.3 Send data to a printer or to a computer

NOTICE

The data output (printer or PC) needs to be selected first, so that the Print menu is available (refer to [Select the data output](#) on page 11).

Note: Refer to [Report output](#) on page 13 to select the report output type.

1. From the main menu use the ▲ or ▼ to select DATA LOGGER and confirm.
2. Use the ▲ or ▼ to select Print and confirm. Select one of the following options and confirm with ✓ to print the data: Measurement data, Electrode data, Calibration data, Calibration report or Instrument condit.

8.4 Report output

NOTICE

The data output (printer or PC) needs to be selected first, so that the Type of report menu is available (refer to [Select the data output](#) on page 11).

Different report output types can be selected when a printer or a PC is connected.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Type of report and confirm.
3. When a printer or a computer and Terminal is connected, use the ▲ or ▼ to select

Option	Description
Reduced	Select Several or One sample as an output format

Option	Description
--------	-------------

Standard	Select Several or One sample as an output format. Select several: Users: The user name appears on the printed report (17 characters). Header: The company name can be added as a header (40 characters) and appears on the print report. Identify sensor: The sensor model and the sensor serial number can be added and appears on the print report.
-----------------	--

GLP	Select Several or One sample as an output format. Select several: Users: The user name appears on the printed report (17 characters). Header: The company name can be added as a header (40 characters) and appears on the print report. Identify sensor: The sensor model and the sensor serial number can be added and appears on the print report.
------------	--

4. When a computer is connected and LabCom Easy (refer to [Select the data output](#) on page 11 for more information) is selected, use the ▲ or ▼ to select

Option	Description
--------	-------------

Users	The user name appears on the printed report (17 characters).
--------------	--

Identify sensor	The sensor model and the sensor serial number can be added and appears on the print report.
------------------------	---

Section 9 Maintenance

▲ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

NOTICE

Do not disassemble the instrument for maintenance. If the internal components must be cleaned or repaired, contact the manufacturer.

9.1 Clean the instrument

NOTICE

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the display and accessories.



Clean the exterior of the instrument with a moist cloth and a mild soap solution.

9.2 Clean the probe

Clean the probe as needed. Refer to [Troubleshooting](#) on page 16 for more information about cleaning. Refer to the probe documentation for information about the probe maintenance.

9.3 Replace the magnetic stirrer

If the magnetic stirrer does not start, follow the numbered steps to replace the magnetic stirrer.

Section 10 Troubleshooting

Refer to the following table for common problem messages or symptoms, possible causes and corrective actions.

Table 1 Calibration warnings and errors

Error/Warning	Solution
Calibration out of range	Measured value out of range. Calibrate again. Connect a new probe.
Unknown buffer	Calibrate again.
Same buffer / standard. Re-calibrate	Examine the probe: Clean the probe (refer to Clean the probe on page 14 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter. Examine the buffer solution: Make sure that the buffer used matches the buffer specified in configuration; make sure of the temperature specification in configuration; use a new buffer solution.
Unstable measurement Time > 100 s (pH, EC and DO calibrations) Time > 240 s (ISE calibrations)	Calibrate again. Examine the probe: Clean the probe (refer to Clean the probe on page 14 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter. Make sure that the membrane and the diaphragm are properly immersed in the sample.
Temperature difference > 3.0 °C	Adjust the calibration solutions to the same temperature. Examine the temperature sensor.
Temperature out of range	Examine the temperature sensor. Connect a new probe.

Table 1 Calibration warnings and errors (continued)

Error/Warning	Solution
Outside allowable range	Offset or Slope outside the range. Examine the buffer solution: Make sure that the buffer used matches the buffer specified in configuration; make sure of the temperature specification in configuration; use a new buffer solution. Examine the probe: Clean the probe (refer to Clean the probe on page 14 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Signal too low / high (DO calibrations)	DO probe error. Examine the probe. Connect a different probe to make sure that if problem is with probe or meter. Use a new standard solution.
Cell constant over limits (EC calibrations)	Insert the probe into the appropriate standard and re-read. Examine the probe: Clean the probe (refer to Clean the probe on page 14 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Cell constant deviation error (EC calibrations)	Examine the probe: Clean the probe (refer to Clean the probe on page 14 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Not calibrated	There are no calibration data stored in the instrument. Perform calibration.

Table 2 Measurement warnings and errors

Error/Warning	Solution
-----	Measured value is out of range. Examine the probe: Clean the probe (refer to Clean the probe on page 14 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Unstable measurement (pH, EC and DO measurements) Time > 120 s	Make sure that the membrane and the diaphragm are properly immersed in the sample. Examine the temperature. Examine the probe: Clean the probe (refer to Clean the probe on page 14 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to make sure that if problem is with probe or meter.
10.389 mg/L >>>>> 00012 00:13	ISE measurement: The measured value is superior to the highest standard used in calibration. Measure again.
0.886 mg/L <<<<<< 00018 00:11	ISE measurement: The measured value is inferior to the lowest standard used in calibration. Measure again.

Section 11 Replacement parts and accessories

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
sensION+ PH3 Lab pH-meter with accessories, without probe	LPV2000.97.0002
sensION+ PH31 Lab pH-meter, GLP, with accessories, without probe	LPV2100.97.0002
sensION+ MM340 Lab pH & Ion-meter, GLP, 2 channels, with accessories, without probe	LPV2200.97.0002

Replacement parts (continued)

Description	Item no.
sensION+ EC7 Lab conductivity meter, with accessories, without probe	LPV3010.97.0002
sensION+ EC71 Lab conductivity meter, GLP, with accessories, without probe	LPV3110.97.0002
sensION+ MM374, 2 channel Lab meter, GLP, accessories, without probes	LPV4110.97.0002
sensION+ MM378, 2 channel Lab meter, GLP, accessories, without probes	LPV4130.97.0002

Consumables

Description	Item no.
pH buffer solution 4.01, 125 mL	LZW9460.99
pH buffer solution 7.00, 125 mL	LZW9461.97
pH buffer solution 10.01, 125 mL	LZW9470.99
pH buffer solution 4.01, 250 mL	LZW9463.99
pH buffer solution 7.00, 250 mL	LZW9464.97
pH buffer solution 10.01, 250 mL	LZW9471.99
pH buffer solution 4.01, 1000 mL	LZW9466.99
pH buffer solution 7.00, 1000 mL	LZW9467.97
pH buffer solution 10.01, 1000 mL	LZW9472.99
Electrolytic solution (KCl 3M), 125 mL	LZW9510.99
Electrolytic solution (KCl 3M), 250 mL	LZW9500.99
Electrolytic solution (KCl 3M), 50 mL	LZW9509.99
Electrolytic solution 0.1 M, 125 mL	LZW9901.99
Enzyme solution	2964349

Consumables (continued)

Description	Item no.
Pepsin Cleaning Solution	2964349
Electrode cleaning solution	2965249
0.1 N HCl solution	1481253
Ethanol, 95% (Grease, oils, fats)	2378900

Accessories

Description	Item no.
Thermal printer, USB, for sensION+ benchtop instruments	LZW8203.97
Thermal paper for printer LZW8203, bag with 4 rolls	LZM078
Power supply for sensION+ benchtop instruments, 230-115 VAC	LZW9008.99
LabCom Easy PC SW, for sensION+ GLP, CD, cable, USB adapter	LZW8997.99
LabCom PC SW, for sensION+ GLP, CD, cable, USB adapter	LZW8999.99
Magnetic stirrer with sensor holder, for sensION+ multimeters	LZW9319.99
3x50 mL printed beakers for pH calibration	LZW9110.97
Three-sensor holder, for sensION+ benchtop instruments	LZW9321.99
Radiometer probe holder for sensION+ benchtop instruments	LZW9325.99
Pyrex glass chamber, continuous flow measurements	LZW9118.99
PP protector, electrode storage	LZW9161.99
PTFE coated stirring bar, 20 x 6 mm	LZW9339.99

Section A Standard solutions

Technical buffer solutions (DIN 19267)

Refer to [Table 3](#) pH and ORP (mV) values of specific buffer sets at varying temperatures.

Table 3 pH, ORP (mV) and temperature values

Temperature		pH					mV
°C	°F						
0	32	2.01	4.01	7.12	9.52	10.30	—
10	50	2.01	4.00	7.06	9.38	10.17	245
20	68	2.00	4.00	7.02	9.26	10.06	228
25	77	2.00	4.01	7.00	9.21	10.01	220
30	86	2.00	4.01	6.99	9.16	9.96	212
40	104	2.00	4.03	6.97	9.06	9.88	195
50	122	2.00	4.06	6.97	8.99	9.82	178
60	140	2.00	4.10	6.98	8.93	9.76	160
70	158	2.01	4.16	7.00	8.88	—	—
80	176	2.01	4.22	7.04	8.83	—	—
90	194	2.01	4.30	7.09	8.79	—	—

Buffer solutions (DIN 19266)

Refer to [Table 4](#) for pH values of specific buffer sets at varying temperatures.

Table 4 pH and temperature values

Temperature		pH						
°C	°F							
5	32	1.668	4.004	6.951	7.087	9.395	10.245	13.207
10	50	1.670	4.000	6.923	7.059	9.332	10.179	13.003
20	68	1.675	4.001	6.881	7.016	9.225	10.062	12.627

Table 4 pH and temperature values (continued)

Temperature		pH						
°C	°F							
25	77	1.679	4.006	6.865	7.000	9.180	10.012	12.454
30	86	1.683	4.012	6.853	6.987	9.139	9.966	12.289
40	104	1.694	4.031	6.838	6.970	9.068	9.889	11.984
50	122	1.707	4.057	6.833	6.964	9.011	9.828	11.705
60	140	1.723	4.085	6.836	6.968	8.962	—	11.449
70	158	1.743	4.126	6.845	6.982	8.921	—	—
80	176	1.766	4.164	6.859	7.004	8.885	—	—
90	194	1.792	4.205	6.877	7.034	8.850	—	—

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)