
DOC026.97.80097

HIAC 9703+ Liquid Particle Counter

05/2013, Edition 4

User Manual
Manuel d'utilisation
Manual del usuario
Manual do Usuário
用户手册
取扱説明書
사용 설명서

English.....	3
Français.....	21
Español.....	39
Português.....	57
中文.....	75
日本語.....	93
한글.....	111

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	343 x 337 x 482 mm (13.5 x 13.3 x 19.0 in.)
Weight	10.7 kg (23.6 lb)
Power requirement	External Class III power adapter: 100–240 VAC, 50–60 Hz, 2.5 A input; 24 VDC, 3.75 A output
	Instrument: 24 VDC, 75 W maximum
Viscosity limits	< 50 cp
Coincidence loss limit	7,500
Operating temperature	5 °C to 40 °C (41 °F to 104 °F)
Operating humidity	0 to 80% relative humidity for temperatures up to 31 °C decreasing linearly to 50% RH at 40 °C
Maximum altitude	2,000 m (6,560 ft)
Sample bottle clearance	153 mm (6.02 in.)
Sample temperature	5 °C to 40 °C (41 °F to 104 °F)
Flow rate	10 to 100 mL/min (determined by the sensor calibration flow rate)
Flow rate accuracy	> 95%
Volume accuracy	> 95% of syringe size ¹
Tare volumes for probes	0.091 mL for 8.1 cm (3.2 in.)/1.2 mm (0.047 in.) ID
	0.172 mL for 15.4 cm (6.1 in.)/1.2 mm (0.047 in.) ID
	1.57 mL for 14 cm (5.5 in.)/large bore

Specification	Details
Sensor concentration limits (particle per ml)	HRLD150: 18,000
	HRLD150JA: 18,000
	HRLD400: 10,000
	HRLD600JS: 6,000
	MC05: 9,000
Wetted materials	PTFE Teflon [®] , TFE Teflon, PFA Teflon, Kel-F, 316 Stainless steel, Glass, Viton [®]
Number of channels	4 (2 for light scatter, 2 for light extinction)
Number of bins	65,536 (16,384 bins per channel)
Certification	CE

¹ Under ideal conditions, inaccuracy in flow rate is ≤1% of syringe size.

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements.

Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

⚠ DANGER
Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.
⚠ WARNING
Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.
⚠ CAUTION
Indicates a potentially hazardous situation that may result in minor or moderate injury.
NOTICE
Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol, if noted on the instrument, will be included with a danger or caution statement in the manual.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.

	This symbol indicates the presence of devices sensitive to Electrostatic Discharge (ESD) and indicates that care must be taken to prevent damage with the equipment.
	This symbol indicates a laser device is used in the equipment.
	This symbol indicates the need for protective eye wear.
	This symbol indicates a potential pinch hazard.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user.

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:


1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Product overview

The HIAC 9703+ liquid particle counter system measures particles in a variety of liquids with a syringe sampler and sensor. The instrument must connect to a computer for operation and data management. Refer to [Figure 1](#) and [Figure 2](#).


An external sensor is used with the HIAC 9703+ liquid particle counter system. Refer to [Figure 1](#), item 2. The external sensors are classified as Class 1 Laser Products. Refer to the sensor manual for FDA and IEC compliance.

Figure 1 HIAC 9703+ liquid particle counter

1 Syringe	4 Sample target
2 Sensor	5 Docking module
3 Sample probe	6 Front panel

Figure 2 Back panel


1 Drain fitting	4 USB connector
2 Service-only ports	5 Power connector
3 Sensor connector	

Product components


Make sure that all components have been received. Refer to [Figure 3](#), [Figure 4](#) and [Figure 5](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 3 HRLD components


1 HIAC 9703+ liquid particle counter	3 Cable, sensor
2 Sensor, HRLD	

Figure 4 MC05 components

1 HIAC 9703+ liquid particle counter	2 Sensor, MC05
--------------------------------------	----------------

Figure 5 Shipping kit

1 Nut	7 Probe, small bore
2 Hex Key, 1.5 mm	8 Cable, USB
3 Syringe, 10 ml	9 Power supply, 24V
4 Power cord, AC	10 O-rings
5 Stir Bar, Teflon	11 Tubing
6 Probe, large bore	12 Power cord, European

Installation

⚠ WARNING

Electrical shock hazard. Externally connected equipment must have an applicable country safety standard assessment.

⚠ CAUTION

Multiple hazards. Only qualified personnel should conduct the tasks described in this document.





Site location




Position the instrument on a dry, level surface. The location must be indoors and near an electrical outlet and a computer. Make sure the power switch is accessible.

Install the HRLD sensor

The sensor is held in place by the holding arm. A sampling probe and tubing must connect to the sensor for sample flow. A cable must connect to the sensor for signal transmission. To install the sensor, refer to the illustrated steps.

Note: Do not put Teflon tape on the fitting threads.








Install the MC05 sensor

The sensor is held in place by the holding arm. A sampling probe and tubing must connect to the sensor for sample flow. A cable must connect

to the sensor for signal transmission. To install the sensor, refer to the illustrated steps.

Note: Do not put Teflon tape on the fitting threads.


Connect the drain line

NOTICE


Dispose of waste in accordance with local, state and national environmental regulations.

The sample waste from the instrument must be collected in a compatible waste container. Monitor the liquid level in the waste container to prevent overflow.


Install or remove the docking module

Docking modules are available for single or multiple-container sampling.


Install the small-vial clamp

The optional small-vial clamp connects to the docking module with 2 screws.

Connect the power cord

A power supply unit is used to supply power to the instrument.

Remove the magnetic stirrer


The magnetic stirrer in the docking module can be removed if necessary.

1. Pull the tab on the docking module away from the instrument to remove the module.
2. Remove the four screws on the bottom plate of the module. Separate the bottom plate from the module enclosure.
3. Remove the screw that holds the magnet to the motor shaft. Hold the shaft with pliers so that the shaft does not move.
4. Remove the magnet and install the screw. Hold the shaft with pliers so that the shaft does not move.
5. Put the bottom plate on the module enclosure and install the four screws.

User interface

The buttons on the front panel allow the user to operate the syringe, sample probe and stirrer. Refer to [Figure 6](#), [Table 1](#) and [Table 2](#).

Figure 6 Front panel

1 Syringe controls	4 Stop button
2 Stirrer speed (RPM)	5 Sample probe controls
3 Power switch	6 LED status indicators

Table 1 Button descriptions

Button	Name	Description
	Load	Sends the syringe plunger to the top or bottom position.
	Clean	Flushes the system.
	Back-flush	Back-flushes the system when held for 3 seconds.

Table 1 Button descriptions (continued)

Button	Name	Description
	Stop	Stops the movement of the syringe and probe. Sampling is stopped. The stirrer continues to operate.
 	Set	Moves the sample probe up or down. Hold to increase speed.
	Lift/Lower	Sends the sample probe to the top or bottom position.

Table 2 LED status indicators

LED color	Description
Green	Instrument is powered on
Yellow	Instrument is sampling
Flashing yellow	Shows the boot up sequence

Operation

⚠ DANGER	
	Fire hazard. This product is not designed for use with flammable liquids.

⚠ WARNING	
	Explosion hazard. Make sure that the drain tube is free of all obstructions. If the drain tube has a blockage or is pinched or bent, high pressure can build up in the instrument.

⚠ CAUTION	
 	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Turn the instrument on or off


When power is applied to the instrument, the sampling probe automatically moves to the top position and the magnetic stirrer starts. The display on the front panel shows the stirrer speed in RPM.

1. Make sure that there are no obstructions to the movement of the probe.
2. Push the **ON/OFF** switch to turn the instrument on or off.
*Note: Push the **STOP** button to immediately stop the movement of the syringe and sample probe.*

Attach the syringe

The syringe drive mechanism must be at the bottom position to attach or remove the syringe.

Figure 7 Attach the syringe

1 Plunger bolt	3 Syringe drive
2 Syringe	4 Valve

1. Push the **ON/OFF** switch to turn the instrument on.
2. Push the **LOAD** button under the syringe compartment to set the plunger bolt at the lowest point. Refer to [Figure 7](#).
3. Loosen and remove the plunger bolt from the syringe drive.
4. Hold the syringe and push the bolt through the plunger hole.
5. Push the bolt in the syringe drive hole and turn the bolt by hand to tighten.
6. Pull the syringe barrel up and install the syringe threads in the valve.

7. Turn the syringe by hand to tighten.

Note: Do not overtighten.

Set the sample probe position

The lift mechanism for the sample probe can be set to automatically go to the correct position for sampling.

1. Put a sample container on the module. If a small vial is used, refer to [Put a vial in the small-vial clamp](#) on page 18.
2. Push the **SET** button to lower the sample probe to the correct position.
3. Push the **LIFT/LOWER** button to raise the sample probe to the top position.
4. Push the **LIFT/LOWER** button to automatically move the sample probe down and stop at the correct bottom position.
Note: This position is recognized by the external software.
The bottom position is saved until a new bottom position is set.


Flush the syringe

Flush the syringe when a new syringe is installed or when a new sample is used. To optimize the performance of the syringe:

- Do not operate the syringe for more than a few cycles without liquid.
 - Thoroughly flush the syringe with distilled or deionized water after use.
1. Push the **LIFT/LOWER** button to raise the sample probe.
 2. Put the container with the flushing liquid under the sample probe.
 3. Push the **LIFT/LOWER** button to lower the sample probe in the flushing liquid.
 4. Push the **CLEAN** button. The syringe fills with the flushing liquid and pushes the liquid out to the drain.
Note: The system can also be flushed with the external software.

Put a vial in the small-vial clamp

The small-vial clamp holds small vials under the center of the sample probe. To install the clamp assembly, refer to [Install the small-vial clamp](#) on page 14.

Stop all operations

The **STOP** button can be used to simultaneously stop movement of the syringe and sample probe. Use the buttons or external software to resume operations.

Connect the computer cable

A USB cable connects the instrument to a computer.

1. Install a compatible software application on the computer.
2. Turn on the instrument power.
3. Connect the USB cable to the USB port on the instrument and to the computer.

Operate the instrument by computer

The instrument must connect to a computer for configuration, calibration, data collection and data management. Refer to the software user manual for specific instructions.

Maintenance

▲ CAUTION

Personal injury hazard. Only qualified personnel should conduct the tasks described in this section of the manual.

NOTICE

Do not disassemble the instrument for maintenance. If the internal components must be cleaned or repaired, contact the manufacturer.

Clean the instrument

Clean the exterior of the instrument with a moist cloth and a mild soap solution and then wipe the instrument dry.

Clean the syringe

Frequent use and dirty sample material make it necessary to clean the syringe more often. Clean the syringe daily or weekly to make sure that it works properly. Use a mild detergent for standard applications. To find whether a different cleaning solution is compatible with the instrument, refer to [Specifications](#) on page 3.

1. Put a container with the cleaning solution on the docking module.
2. Push the **LIFT/LOWER** button to lower the sample probe.
3. Push the **LOAD** button to fill the syringe with the cleaning solution.
Note: If the plunger does not move easily in the syringe barrel, remove the plunger, wipe the plunger with alcohol and wet the plunger seal with deionized water.
4. Allow the solution to stay in the syringe for a minimum of 30 minutes.
5. Push the **LOAD** button to send the cleaning solution to the drain.
6. Push the **LIFT/LOWER** button to raise the sample probe.
7. Put a container with distilled or deionized water on the docking module.
8. Push the **LIFT/LOWER** button to lower the sample probe.
9. Push the **CLEAN** button to flush the syringe with the rinse water.
10. Repeat the flush with rinse water for a minimum of 10 cycles.

Troubleshooting

Back-flush the system

A large particle can reduce or stop the liquid flow. If this occurs, a flush in the reverse direction can push the particle out through the sample probe.

1. Put a small waste container under the sample probe.
2. Push the **LIFT/LOWER** button to lower the sample probe.
3. Push and hold the **BACK-FLUSH** button until the back-flush cycle starts.
The syringe fills with liquid from the drain tubing. When the syringe closes, the liquid is sent through the sensor and sample probe.
4. Repeat the back-flush cycle as necessary until a flush in the forward direction with a clean sample flows at a normal rate.

How to clean spills

▲ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

1. Obey all facility safety protocols for spill control.
2. Discard the waste according to applicable regulations.

Replacement parts and accessories

▲ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)