EZ Arsenic

For test kit 2822800 (EZ Arsenic)

0 to 500, 0 to 4000 ppb As

DOC326.98.00006


Additional copies available on www.hach.com

Test preparation


- For samples with sulfide greater than 15 ppb, follow the [Optional procedure for removing sulfide](#) before performing the test.
- Do not expose reacted strips to direct sunlight. Reaction products are photosensitive and may turn dark.
- Do not allow test strips to touch the reaction vessel solution. Test strips react with gases, not solution.
- Orient the test strip pad **paper side down** and **centered** over the hole in the black cap so the generated gases can make good contact with the pad.
- Two samples may be analyzed simultaneously with this kit.

DANGER: Hydrogen and arsine gasses are generated during the test. Work in a well-ventilated area away from open flames and other sources of ignition. Review Material Safety Data Sheets for safe handling, storage and disposal information.

0–500 ppb (0, 10, 25, 50, 100, 250, 500 ppb) test procedure

0–4000 ppb (0, 35, 75, 175, 500, 1500, 4000 ppb) test procedure

Required items

Description	Unit	Catalog no.
EZ Arsenic Reagent Set	each	2823200
Reaction Bottle	each	2800200
Cap, reaction bottle	each	4934800

Interferences

Refer to [Table 1](#) and [Table 2](#) for interference information. Interferences other than those listed are not likely.

Table 1 Interfering substances

Ion or Substance	Concentration
Acidity	< pH 5. Do not acid-preserve samples. If samples are below pH 5, adjust pH to between 5 and 6 before beginning test.
Antimony	> 250 ppb
Nitric acid	Interferes with the reduction step. Do not use samples preserved with nitric acid because low results will be observed. If samples must be preserved, use HCl or sulfamic acid to adjust sample to pH 2. Adjust to pH 5–7 before running the test.
Selenium	> 1 ppm
Sulfide	> 15 ppb (see Optional procedure for removing sulfide)
Tellurium	Likely to interfere, but not tested.

Table 2 Non-interfering substances

Ion or Substance	Highest concentration tested
Alkalinity	1000 ppm as CaCO ₃
Hardness	1000 ppm as CaCO ₃
Iron	100 ppm
Temperature	10 to 40 °C (50 to 104 °F)

Optional procedure for removing sulfide

If a rotten egg smell is detected after adding reagent #1, sulfide is present at interfering levels. Complete the following steps to remove the sulfide before beginning the test procedure:

- Tear off a small piece of cotton and form a ball the size of a pea.
- Saturate the cotton with a few drops of lead acetate. Squeeze the excess liquid out of the cotton, leaving it damp.
- Press the saturated cotton ball into the small opening of the reaction bottle cap from the bottom. Be sure that the cotton is firmly in place and that a gap remains between the cotton and the top surface of the cap.
- Insert the test strip as detailed in step 1 of the 0–500 or 0–4000 ppb test procedure and continue with the test.

Note: The lead acetate must not contact the test strip!

Always wear gloves or wash hands thoroughly after handling lead acetate.

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