
MC3051Pt-9

Combined Redox Electrode

Electrode de Platine Combinée

MC3051Pt-9 Combined Redox Electrode

Introduction

The MC3051Pt-9 Combined Platinum Electrode is ideal for a wide range of redox measurements. The special electrode design with the platinum sensing ring and the Ag/AgCl reference element built into the same unit ensures convenient handling during use.

Preparation for Measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole.
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.
If necessary, refill with KCl•Ag 3M KCl Solution, saturated with AgCl.
4. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the sensing element down.

Maintenance

1. Electrode contamination is a major cause of faulty measurements.
2. The electrode should be rinsed with distilled water after measurements.
3. Check frequently the level of filling solution.
4. In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc...
 - KS400 Pepsin in HCl Solution or RENOVO•X Xtra Strong Cleaning Solution: protein deposits (milk, cheese, serums...). Duration of treatment 1 to 2 hrs.
 - KS410 Thiourea solution: for porous plugs contaminated with sulphides or blocked by an AgCl precipitate. Duration of treatment, a few hours until the porous pin turns white.

- KS400 Pepsin in HCl Solution or RENOVO•X Xtra Strong Cleaning Solution: protein deposits (milk, cheese, serums...). Duration of treatment 1 to 2 hrs.
- RENOVO•N Normal Cleaning Solution: greasy or oily deposits...
- Chemical cleaning: the electrode can be cleaned with a sulphochromic solution, or with a boiling solution of hydrochloric acid.

Storage

Between measurements and up to 1 week: seal the filling hole with paraffin film or with the electrode stopper and immerse the electrode in KS110 3M KCl Solution.

Long-term storage: the electrode can be stored dry. Clean the electrode tip with distilled water and gently wipe it dry with soft tissue. Fit the electrode protection cap filled with KS110 3M KCl Solution, place the electrode in its box and store at room temperature.

Accessories

KCl•Ag 3M KCl Solution saturated with AgCl, 100 ml	S21M004
KS110 3M KCl Solution, 500 ml	C20C320
○ RENOVO•N Normal Cleaning Solution, 250 ml	S16M001
□ RENOVO•X Xtra Strong Cleaning Solution, 250 ml	S16M002
KS400 Pepsin in HCl Solution, 250 ml	C20C370
KS410 Thiourea Solution, 250 ml	C20C380
CL114 Electrode Cable, (1 m) with BNC plug	A94L114
CL116 Electrode Cable, (1 m) with type 7 plug	A94L116

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RADIOMETER ANALYTICAL S.A.
D31M090 • Radiometer Analytical S.A. • France • 9809A