

---

# Potassium

---

For water and wastewater

Tetraphenylborate Method

---

## Introduction

Potassium, one of the most abundant elements, is found in many minerals. Soils contain approximately 1 to 4% potassium. Concentrations of potassium in most drinking water is usually less than 20 mg/L; occasionally brines may contain more than 100 mg/L. The greatest areas of interest in measurement of potassium levels probably are medicine and agriculture, due to the importance of potassium as a mineral for plants and animals. Potassium salts, particularly potash, are common in fertilizers.

The Tetraphenylborate Method for determination of potassium in water is accurate, rapid, and inexpensive. In the reaction, a precipitate is formed and the resulting increase in turbidity is measured. All necessary reagents are packaged in three powder pillows to provide reagent stability, convenience and accuracy.

## Chemical reactions

Potassium combines with sodium tetraphenylborate to form potassium tetraphenylborate, a white precipitate. The precipitate remains in suspension in samples with low concentrations of potassium, causing an increase in turbidity.



Figure 1 Chemical reaction between potassium and tetraphenylborate

The sodium tetraphenylborate is contained in Potassium 3 Reagent Powder Pillows. Ammonium salts, magnesium and calcium interfere with the precipitation. Potassium 1 Reagent Powder Pillows and Potassium 2 Reagent Powder Pillows prevent these interferences.