




Method 8345

Direct Reading Method

MR (0 to 50 mg/L)

Scope and Application: For drinking and wastewater.

The estimated detection limit for program number 1510 is 0.6 mg/L.



1. Press the soft key under **HACH PROGRAM**.

Select the stored program number for mid range chlorine dioxide (ClO₂) by pressing **1510** with the numeric keys.

Press: **ENTER**

Note: Analyze samples immediately. See *Sample Collection, Storage and Preservation* following these steps.


Note: The Flow Cell and Sipper Modules can be used with this procedure. Use minimum volumes of 20 and 10 mL, respectively.




2. The display will show:

**HACH PROGRAM:
1510 Chlor: Dioxide,
MR**


The wavelength (λ), **360 nm**, is automatically selected.



3. Fill a sample cell to the 10-mL mark with deionized water (the blank).



4. Place the blank into the cell holder. Close the light shield.



5. Press the soft key under **ZERO**.


The display will show:

0.0 mg/L ClO₂

Note: For alternate concentration units, press the soft key under **OPTIONS**. Then press the soft key under **UNITS** to scroll through the available options. Press **ENTER** to return to the read screen.



6. Fill another sample cell to the 10-mL mark with sample (the prepared sample).



7. Place the prepared sample into the cell holder. Close the light shield. The results in mg/L ClO₂ (or chosen units) will be displayed.

Sample Collection, Storage and Preservation

Collect samples in clean plastic or glass bottles. Samples must be analyzed immediately. Chlorine dioxide is very volatile and unstable.

Accuracy Check

Standard Solution Method

Preparing chlorine dioxide standards is difficult and dangerous. In addition, **these standards are both explosive and volatile!** Only a trained chemist should prepare the standards using appropriate safety equipment and precautions. Hach does not recommend preparation of chlorine dioxide standards. If independent standard preparation is required, please see the instructions in *Standard Methods for the Examination of Water and Wastewater*, 18th ed., under the headings “Stock chlorine dioxide solution” and “Standard chlorine dioxide solution” (pg. 4–54). Prepare a 25.0-mg/L chlorine dioxide standard.

Method Performance

Precision

Standard: 25.0 mg/L ClO₂

Program	95% Confidence Limits
1510	24.5–25.5 mg/L ClO ₂

For more information on determining precision data and method detection limits, refer to Section 1.5.

Estimated Detection Limit

Program	EDL
1510	0.6 mg/L ClO ₂

For more information on derivation and use of Hach’s estimated detection limit, see Section 1.5.2. To determine a method detection limit (MDL) as defined by the 40 CFR part 136, Appendix B, see Section 1.5.1.

Sensitivity

Program Number: 1510

Portion of Curve	ΔAbs	ΔConcentration
0.010 Abs	0.010	0.28 mg/L
27.5 mg/L	0.010	0.25 mg/L
49.5 mg/L	0.010	0.22 mg/L

See Section 1.5.3 *Sensitivity Explained* for more information.

Calibration Standard Preparation

Preparing chlorine dioxide standards is difficult and dangerous. In addition, **these standards are both explosive and volatile!** Only a trained chemist should prepare the standards using appropriate safety equipment and precautions. Hach does not recommend preparation of chlorine dioxide standards. If independent standard preparation is required, please see the instructions in *Standard Methods for the Examination of Water and Wastewater*, 18th ed., under the headings “Stock chlorine dioxide solution” and “Standard chlorine dioxide solution” (pg. 4–54). Using the standards prepared above and the analysis procedure, generate a calibration curve.

Summary of Method

Chlorine dioxide, a yellow gas, can be measured directly in a water solution. This method uses a wavelength of 360 nm to increase the sensitivity of the test.

Safety

Good safety habits and laboratory techniques should be used throughout the procedure. Consult the *Material Safety Data Sheet* for information specific to the reagents used. For additional information, refer to Section 1.

CHLORINE DIOXIDE, continued

Pollution Prevention and Waste Management

For information on pollution prevention and waste management, refer to Section I.

REQUIRED REAGENTS AND STANDARDS

Description	Quantity Required		Unit	Cat. No
	Per Test			
Water, deionized	10 mL	4 liters		272-56

REQUIRED EQUIPMENT AND SUPPLIES

DR/4000 1-Inch Cell Adapter	1	each.....	48190-00
Sample Cells, 1-inch, matched pair.....	2	2/pkg.....	26126-02

OPTIONAL EQUIPMENT AND SUPPLIES

DR/4000 Carousel Module Kit		each.....	48070-02
DR/4000 Flow Cell Module Kit, 1-inch.....		each.....	48070-04
DR/4000 Flow Cell Module Kit, 1-cm.....		each.....	48070-05
DR/4000 Sipper Module Kit, 1-inch.....		each.....	48090-03



FOR TECHNICAL ASSISTANCE, PRICE INFORMATION AND ORDERING:
In the U.S.A. – Call toll-free 800-227-4224
Outside the U.S.A. – Contact the HACH office or distributor serving you.
On the Worldwide Web – www.hach.com; E-mail – techhelp@hach.com

HACH COMPANY
WORLD HEADQUARTERS
Telephone: (970) 669-3050
FAX: (970) 669-2932