



Method 8196


Esterification Method*

(0 to 2800 mg/L)

Scope and Application: For digester sludges.

The estimated detection limit for program number 3800 is 12 mg/L as acetic acid (HOAC).

* Adapted from *The Analyst*, 87, 949 (1962)



1. Press the soft key under **HACH PROGRAM**.

Select the stored program number for volatile acids by pressing **3 8 0 0** with the numeric keys.


Press: **ENTER**

Note: If samples cannot be analyzed immediately, see *Sample Collection, Preservation and Storage* following these steps.




2. The display will show: **HACH PROGRAM: 3800 Volatile Acids**

The wavelength (λ), **495 nm**, is automatically selected.




3. Pipet 0.5 mL of deionized water into a dry 25-mL sample cell (the blank).




4. Filter or centrifuge 25 mL of sample using the labware listed under **REQUIRED EQUIPMENT AND SUPPLIES**.

Note: Centrifuging is faster than filtration.




5. Pipet 0.5 mL of the filtrate or supernatant into another dry 25-mL sample cell (the prepared sample).


Note: For proof of accuracy, use 0.5 mL of a 500-mg/L volatile acid standard solution in place of the sample (see *Accuracy Check*).



6. Pipet 1.5 mL of ethylene glycol into each sample cell. Swirl to mix.




7. Pipet 0.2 mL of 19.2 N Sulfuric Acid Standard Solution into each cell. Swirl to mix.



8. Place both cells into a boiling water bath.


Note: You may boil the sample cells in a 600-mL beaker.

VOLATILE ACIDS, continued




9. Press the soft key under **START TIMER**.


Note: A 3-minute reaction period will begin.




10. When the timer beeps, cool the solutions to 25 °C (until the cell feels cold) with running tap water.




11. Using a pipet filler, pipet 0.5 mL of Hydroxylamine Hydrochloride Solution into each cell. Swirl to mix.




12. Using a pipet filler, pipet 2.0 mL of 4.5 N Sodium Hydroxide Standard Solution into each cell. Swirl to mix.



13. Add 10 mL of Ferric Chloride Sulfuric Acid Solution to each cell. Swirl to mix.




14. Add 10 mL of deionized water into each cell. Swirl to mix.




15. Immediately press the soft key under **START TIMER**. Another 3-minute reaction period will begin.

Note: During this 3-minute reaction period, complete Steps 16-17.



16. Blot each sample cell dry. Immediately place the blank in the cell holder. Close the light shield.




17. Press the soft key under **ZERO**.

The display will show:

0 mg/L HOAC

Note: For alternate concentration units, press the soft key under **OPTIONS**. Then press the soft key under **UNITS** to scroll through the available options. Press **ENTER** to return to the read screen.



18. When the timer beeps, immediately place the prepared sample in the cell holder. Close the light shield. Results in mg/L volatile acids as acetic acid will be displayed.

Sample Collection, Preservation and Storage

Collect samples in clean plastic or glass bottles. Analyze as soon as possible after collection. Samples can be stored for up to 24 hours by cooling to 4 °C (39 °F) or below. Warm samples to room temperature before analysis.

Accuracy Check

Standard Additions Method

- a. Leave the unspiked sample in the sample compartment. Verify that the units displayed are in mg/L. Select standard additions mode by pressing the soft keys under **OPTIONS, (MORE)** and then **STD ADD**.
- b. Press **ENTER** to accept the default sample volume (mL), 25.
- c. Press **ENTER** to accept the default standard concentration (mg/L), 62,500.
- d. Press the soft key under **ENTRY DONE**.
- e. Snap the neck off a Volatile Acid Voluette® Ampule Standard, 62,500 mg/L as acetic acid.
- f. Use the TenSette® Pipet to add 0.1 mL, 0.2 mL, and 0.3 mL of standard, respectively, to three 25-mL samples and mix each thoroughly.
- g. Analyze each standard addition sample as described above. Accept the standard additions reading by pressing the soft key under **READ** each time. Each addition should reflect approximately 100% recovery.
- h. After completing the sequence, the display will show the extrapolated concentration value and the “best-fit” line through the standard additions data points, accounting for matrix interferences.
- i. See Section 1.4.1 *Standard Additions* for more information.

VOLATILE ACIDS, continued

Standard Solution Method

Prepare a 500-mg/L volatile acid standard solution by pipetting 4.00 mL of a 62,500-mg/L Volatile Acid Standard Solution into a 500-mL Class A volumetric flask. Dilute to volume with deionized water. Prepare this solution daily. Perform the esterification procedure as described above.

To adjust the calibration curve using the reading obtained with the 500-mg/L standard solution, press the soft keys under **OPTIONS, MORE** then **STD: OFF**. Press **ENTER** to accept the displayed concentration, the value of which depends on the selected units. If an alternate concentration is used, enter the actual concentration and press enter to return to the read screen. See Section 1.5.5 *Adjusting the Standard Curve* for more information.

Method Performance

Precision

Standard 500 mg/L as acetic acid (HOAC)

Program	95% Confidence Limits
3800	492–507 mg/L HOAC

For more information on determining precision data and method detection limits, refer to Section 1.5.

Estimated Detection Limit

Program	EDL
3800	12 mg/L as HOAC

For more information on derivation and use of Hach's estimated detection limit, see Section 1.5.2. To determine a method detection limit (MDL) as defined by 40 CFR part 136, appendix B, see Section 1.5.1.

Sensitivity

Program Number: 3800

Portion of Curve	Δ Abs	Δ Concentration
Entire Range	0.010	23.4 mg/L

See Section 1.5.3 *Sensitivity Explained* for more information.

Calibration Standard Preparation

To perform a volatile acids calibration using the esterification method, use the 62,500-mg/L Volatile Acids Standard Solution (Cat. No. 14270-10) listed under optional reagents and standards. Prepare standards containing 313, 625, 938, 1250, 1563, 2188, and 2500 mg/L volatile acids as follows:

- a. Pipet 1.00, 2.00, 3.00, 4.00, 5.00, 7.00, and 8.00 mL of the 62,500-mg/L Volatile Acids Standard Solution into seven different 200-mL Class A volumetric flasks, using Class A glassware.
- b. Dilute to the mark with deionized water and mix thoroughly.
- c. Using the esterification method and the calibration procedure described in the *User-Entered Programs* section of the *DR/4000 Spectrophotometer Instrument Manual*, generate a calibration curve from the standards prepared above.

Summary of Method

The volatile acid test is designed specifically for determining volatile acids in digester sludges. The method is based on esterification of the carboxylic acids present in the sample and subsequent determination of the esters by the ferric hydroxamate reaction. All volatile acids present are reported as their equivalent mg/L as acetic acid.

Safety

Good safety habits and laboratory techniques should be used throughout the procedure. Consult the *Material Safety Data Sheet* for information specific to the reagents used. For additional information, refer to *Section 1*.

Pollution Prevention and Waste Management

For information on pollution prevention and waste management, refer to *Section 1*.

REQUIRED REAGENTS AND STANDARDS

	Cat. No.
Volatile Acids Reagent Set (90 tests)	22447-00
Includes: (1) 2039-53, (2) 2042-53, (1) 818-42, (1) 2040-53, (1) 2038-32	

Description	Quantity Required		Cat. No.
	Per Test	Unit	
Ethylene Glycol.....	3 mL	1000 mL.....	2039-53
Ferric Chloride-Sulfuric Acid Solution	20 mL	1000 mL.....	2042-53
Hydroxylamine Hydrochloride Solution, 100-g/L	1 mL	100 mL.....	818-42
Sodium Hydroxide Standard Solution, 4.5 N.....	4 mL	1000 mL.....	2040-53
Sulfuric Acid Standard Solution, 19.2 N.....	0.4 mL.....	100 mL MDB.....	2038-32
Water, deionized	20.5 mL.....	4 L.....	272-56

REQUIRED EQUIPMENT AND SUPPLIES

Cylinder, graduated, 10-mL	1	each.....	508-38
DR/4000 1-Inch Cell Adapter	1	each.....	48190-00
Filter Paper, folded, 12.5-cm.....	1	100/pkg.....	1894-57
Finger Cots	2	2/pkg.....	14647-02
Flask, Erlenmeyer, 50-mL.....	1	each.....	505-41

VOLATILE ACIDS, continued

REQUIRED EQUIPMENT AND SUPPLIES(continued)

Funnel, poly, 65-mm	1	each.....	1083-67
Hot Plate, 4" micro, 120 VAC	1	each.....	12067-01
Hot Plate, 4" micro, 240 VAC	1	each.....	12067-02
Pipet, Filler, safety bulb.....	1	each.....	14651-00
Pipet, TenSette®, 1.0- to 10.0-mL	1	each.....	19700-10
Pipet Tips, for 19700-10 Pipet	2	50/pkg.....	21997-96
Pipet, serological, 2-mL	1	each.....	532-36
Pipet, volumetric, Class A, 0.50-mL.....	1	each.....	14515-34
Pipet, volumetric, Class A, 10.00-mL.....	1	each.....	14515-38
Water Bath and Rack.....	1	each.....	1955-55

OPTIONAL REAGENTS AND STANDARDS

Description	Unit	Cat. No.
Volatile Acids Standard Solution, 10-mL Voluette® Ampule, 62,500 mg/L as HOAC	16/pkg	14270-10
Volatile Acids Standard Solution, 1000-mg/L as HOAC	100 mL	14205-42

OPTIONAL EQUIPMENT AND SUPPLIES

Ampule Breaker Kit	each.....	21968-00
Aspirator, Nalgene* vacuum pump	each.....	2131-00
Bottle, wash, 500-mL	each.....	620-11
Beaker, 600-mL	each.....	500-52
Centrifuge, Spinette, 115 VAC, 60 Hz	each.....	22413-00
Centrifuge Tubes, 15-mL	10/pkg.....	22787-39
Centrifuge Tube Caps.....	20/pkg.....	25852-20
Cylinder, graduated, mixing, 25-mL	each.....	1896-40
Cylinder, graduated, 100-mL	each.....	508-42
Cylinder, graduated, 250-mL	each.....	508-46
DR/4000 Carousel Module Kit	each.....	48070-02
DR/4000 Flow Cell Module Kit, 1-inch.....	each.....	48070-04
DR/4000 Flow Cell Module Kit, 1-cm.....	each.....	48070-05
DR/4000 Sipper Module Kit, 1-inch.....	each.....	48090-03
Filter Paper, 9-cm dia.	100/pkg.....	506-55
Flask, filtering, 500-mL	each.....	546-49
Flask, volumetric, Class A, 200-mL	each.....	14574-45
Funnel, Buchner	each.....	550-87
Pipet, TenSette, 0.1- to 1.0-mL	each.....	19700-01
Pipet Tips, for 19700-01 Pipet	50/pkg.....	21856-96
Pipet, serological, 5-mL	each.....	532-37
Pipet, volumetric, Class A, 1.00-mL.....	each.....	14515-35
Pipet, volumetric, Class A, 2.00-mL.....	each.....	14515-36
Pipet, volumetric, Class A, 3.00-mL.....	each.....	14515-03
Pipet, volumetric, Class A, 4.00-mL.....	each.....	14515-04
Pipet, volumetric, Class A, 5.00-mL.....	each.....	14515-37
Pipet, volumetric, Class A, 7.00-mL.....	each.....	14515-07
Pipet, volumetric, Class A, 8.00-mL.....	each.....	14515-08
Test Tube Holder	each.....	634-00
Tubing, rubber	3.6 m (12 ft).....	560-18
Tweezers, balance weight, plastic	each.....	14282-00

* Nalgene is a registered trademark of Nalge Nunc International



FOR TECHNICAL ASSISTANCE, PRICE INFORMATION AND ORDERING:
In the U.S.A. – Call toll-free 800-227-4224
Outside the U.S.A. – Contact the HACH office or distributor serving you.
On the Worldwide Web – www.hach.com; E-mail – techhelp@hach.com

HACH COMPANY
WORLD HEADQUARTERS
Telephone: (970) 669-3050
FAX: (970) 669-2932