
DOC022.98.80516

# Sample changer AS1000

10/2024, Edition 6

**User Manual**  
**Benutzerhandbuch**  
**Manuale utente**  
**Manuel d'utilisation**  
**Manual del usuario**  
**Manual do utilizador**  
**Нáвод к пoужитí**  
**Brugervejledning**  
**Gebbruikershandleiding**  
**Instrukcja obsługi**  
**Bruksanvisning**  
**Käyttöohje**  
**Ръководство за потребителя**  
**Felhasználói kézikönyv**  
**Manualul utilizatorului**  
**Kullanıcı Kılavuzu**  
**Návod na použitie**  
**Navodila za uporabo**  
**Korisnički priručnik**  
**Εγχειρίδιο χρήστη**

## Table of Contents

---

English .....	3
Deutsch .....	18
Italiano .....	33
Français .....	48
Español .....	63
Português .....	78
Čeština .....	93
Dansk .....	108
Nederlands .....	123
Polski .....	138
Svenska .....	153
Suomi .....	168
български .....	183
Magyar .....	199
Română .....	214
Türkçe .....	229
Slovenský jazyk .....	244
Slovenski .....	259
Hrvatski .....	274
Ελληνικά .....	289

# Table of Contents

- 1 Specifications on page 3
- 2 General information on page 4
- 3 Installation on page 9
- 4 Operation on page 14
- 5 Maintenance on page 15
- 6 Troubleshooting on page 16
- 7 Replacement parts and accessories on page 16

## Section 1 Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x H x D)	395 x 300 x 535 mm (15.55 x 11.81 x 21.06 in.)
Weight approximately	Approximately 10 kg (22.05 lb)
Enclosure	<b>Lifting module:</b> Enamelled steel and PP; <b>Rotating module:</b> PP
Pollution degree	2
Installation category	II
Protection class	I
Operating temperature	15 to 40 °C (59 to 104 °F)
Storage temperature	-5 to 40 °C (23 to 104 °F)
Relative humidity	20 to 80%, non-condensing
Maximum sound power level	60 dBA
Power requirements	<b>Instrument input:</b> 24 VDC, 2.5 A; <b>External plug-in power supply:</b> Input 100–240 VAC, 50-60 Hz, 1.5 A, Class I; output, 24 VDC, 2.5 A
Maximum voltage fluctuation	±10% of nominal voltage
Altitude	2000 m (6562 ft) maximum
Environmental conditions	Indoor use
Magnetic stirrer <sup>1</sup>	RCA connector, maximum voltage 12 V, maximum current 150 mA
Rotating module connection <sup>2</sup>	RJ45 telephone connector
Communications	RJ11 connector, RS232 inputs/outputs
Tray rotational speed	4 rpm (nominal)
Vertical movement (sensor holder)	88 mm (3.5 in.)
Tray positions	12, 20 or 30
Sample beakers (H x D)	<b>50 mL:</b> 96 x 29 mm (3.77 x 1.14 in.); <b>50 mL:</b> 88 x 60 mm (3.46 x 2.36 in.); <b>90 mL:</b> 75 x 43 mm (2.95 x 1.69 in.); <b>150 mL:</b> 88 x 60 mm (3.46 x 2.36 in.)
Certifications	<b>Safety:</b> IEC/EN 61010-1, IEC/EN 61010-2-081, CAN CSA C22.2 No. 61010-1, CAN CSA C22.2 No. 61010-2-081, ANSI/UL 61010-1, ANSI/UL 61010-2-081 <b>EMC:</b> IEC/EN 61326-1

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## Section 2 General information

In no event will the manufacturer be liable for damages resulting from any improper use of product or failure to comply with the instructions in the manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

### 2.1 Safety information

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

If the equipment is used in a manner that is not specified by the manufacturer, the protection provided by the equipment may be impaired. Do not use or install this equipment in any manner other than that specified in this manual.

#### 2.1.1 Use of hazard information

<b>▲ DANGER</b>
Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.
<b>▲ WARNING</b>
Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.
<b>▲ CAUTION</b>
Indicates a potentially hazardous situation that may result in minor or moderate injury.
<b>NOTICE</b>
Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## 2.1.2 Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	This symbol, if noted on the product, indicates a crush hazard. Keep hands and fingers clear.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

## 2.1.3 Certification

### Canadian Radio Interference-Causing Equipment Regulation, ICES-003, Class B:

Supporting test records reside with the manufacturer.

This Class B digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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### FCC Part 15, Class "B" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Move the equipment away from the device receiving the interference.
2. Reposition the receiving antenna for the device receiving the interference.
3. Try combinations of the above.

## 2.1.4 Chemical Safety

### ▲ DANGER

Chemical hazards. If this instrument is used to monitor a treatment process and/or chemical feed system for which there are regulatory limits and monitoring requirements related to public health, public safety, food or beverage manufacture or processing, it is the responsibility of the user of this instrument to know and abide by any applicable regulation and to have sufficient and appropriate mechanisms in place for compliance with applicable regulations in the event of malfunction of the instrument.

### NOTICE

Never use this instrument to do tests on living beings.

Normal operation of this device may require the use of chemicals.

- Do not use the instrument in samples that could cause a biological hazard.
- Observe all cautionary information printed on the original solution containers and safety data sheets prior to their use.
- Dispose of all consumed solutions in accordance with the local and national regulations and laws.
- Select the type of protective equipment suitable to the concentration and quantity of the dangerous material being used.

## 2.2 Product overview


The TitraLab AS1000 series is an automatic sample changer used in analytical laboratories with the TitraLab AT1000 series analyzers. There are three sampler models:

- AS1000.XX.20090: Sampler for TitraLab AT1000 series, 20 beakers, 90 mL
- AS1000.XX.30050: Sampler for TitraLab AT1000 series, 30 beakers, 50 mL
- AS1000.XX.12150: Sampler for TitraLab AT1000 series, 12 beakers, 50/150 mL

The sampler is PC-controlled with the TM1000 software. The instrument has a tray that holds the sample beakers. Refer to [Table 1](#) for instrument configurations. The rotating module turns the tray to change the sample position while the lifting module moves sensors and tubes up and down into each sample. The sampler has an integrated magnetic stirrer. Use the sampler for a wide variety of aqueous sample applications. Refer to the analyzer and software documentation for more information. Refer to [Figure 1](#).

**Note:** The type of tray, probe holder and beakers is dependent on the model. Refer to [Table 1](#).

**Figure 1 Sampler**

1 Lifting module	6 Tray
2 Probe and tubes holder	7 Rotating module
3 Sample beakers	8 Magnetic stirrer
4 Blue LED (Flashes when the tray is in position 1. For the remaining positions, the LED is on.)	9 Power switch
5 Green LED (ON/OFF indicator LED)	10 Rear panel

**Table 1 Instrument configurations**

Model	Tray positions	Beakers	Minimum sample volume (mL) <sup>3</sup>		Tubes <sup>4</sup>	Probes <sup>5</sup>
			Combined	Pt-Pt/ISE		
AS1000.XX.20090	20	90 mL	51	30	3	1
AS1000.XX.30050	30	50 mL	21	11	3	1
AS1000.XX.12150	12	50 mL	46	19	4	2
		150 mL	70	39	4	2

<sup>3</sup> Some probes can require a larger minimum sample volume. Make sure that the sensing part of the probe is fully into the sample during the measurement.


<sup>4</sup> Maximum number of tubes (suction tube not included).

<sup>5</sup> Maximum number of accepted probes.


## 2.3 Product components

Make sure that all components have been received. Refer to [Figure 2](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.


**Figure 2** Sampler components


AS1000.XX.12150

AS1000.XX.20090

AS1000.XX.30050

1 Lifting module	9 RS-232 adapter cable with USB/Serial adapter
2 Tray screw cap	10 TM1000 software
3 Tube holder	11 Tray
4 Rotating module	12 Sample beakers
5 Power cord	13 Magnetic stir bars
6 Power supply	14 Probe holder
7 Base with magnetic stirrer	15 Ring adapters for 50 mL beaker
8 Ethernet cables	

**Note:** The type of tray, probe holder and beakers is dependent on the model.


## 2.4 Instrument connections

Connect the power supply, the rotating module, the magnetic stirrer and the PC to the rear panel of the instrument. Refer to [Figure 3](#) for instrument connections.

External equipment that is intended for connection to the instrument must obey the applicable product standard for such equipment, e.g., IEC 60950-1 or IEC 62368-1 for IT-equipment. The safety of a system that includes external equipment is the responsibility of the user who assembles the system.

**Figure 3 Instrument connections**

1 Magnetic stirrer connection	3 External power supply connection	5 Connections label
2 Rotating module connection	4 PC connection	

## Section 3 Installation

### ⚠ DANGER

Explosion hazard. The instrument is not approved for installation in hazardous locations.

### ⚠ DANGER

Fire hazard. This product is not designed for use with flammable liquids.

### ⚠ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

### ⚠ CAUTION



Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.




### 3.1 Installation guidelines

- This instrument is for indoor use only.
- The plug of the mains supply cord or the input connector of the external power supply must be easily accessible so the power can be disconnected quickly in an emergency.
- Protective earth ground (PE) connection is required.
- Keep the instrument away from temperature extremes, including heaters, direct sunlight and other heat sources.
- Put the instrument on a stable and level surface in a well ventilated place.
- Make sure that there is at least 15 cm (6 in.) of space on all sides of the instrument to prevent electrical parts from overheating.
- Do not operate or keep the instrument in dusty, damp or wet locations.
- Always keep the surface of the instrument and all accessories dry and clean.

### 3.2 Assemble the instrument

Refer to the illustrated steps that follow.

### 3.3 Connect to AC power

#### ⚠ CAUTION

Electrical shock and fire hazards. Make sure that the supplied cord and non-locking plug meet the applicable country code requirements.

#### ⚠ WARNING

Fire hazard. Use only the external power supply that is specified for this instrument.

1. Connect the power cord to the power supply.
2. Connect the power supply to the instrument (refer to [Instrument connections](#) on page 8).
3. Connect the power cord to an electrical outlet. The electrical outlet must have a connection for protective earth ground (PE).

### 3.4 Install the probe holder and tubes

#### ▲ WARNING

Pinch hazard. Parts that move can pinch and cause injury. Do not touch moving parts.

#### ▲ WARNING

Chemical hazard. Make sure that the reagent tubes and all the sample beakers are in the correct position in the tray to avoid leaks and the potential escape of reagents. Make sure that the tubes and beakers are installed before using the instrument. Do not use the instrument if the tubes or beakers are not installed.

#### NOTICE


It is necessary to install the probes correctly. Make sure that:




- The sensing part of the probe is fully into the sample.
- The probe will not hit the magnetic stir bar when the holder moves down.
- When the holder is in the upper position, the probe will not hit the beakers when the tray moves.
- The probe adapter is in the correct position. Different sensor requires different probe adapter position.


#### Pre-requisites:

- Make sure that the instrument is connected to a PC with the TM1000 software installed. Refer to [Operate the instrument by computer](#) on page 14.
- Set the instrument power to on. The rod moves to the upper position and the tray turns to position 1. The blue LED flashes to show that the tray is in position 1.

Refer to the illustrated steps that follow to complete the probe holder and tube installation:

**6**

## Section 4 Operation

**▲ DANGER**


Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

### 4.1 Operate the instrument by computer

The instrument must connect to a computer for configuration. Refer to the TM1000 software help for specific information.

Refer to [Figure 4](#) for a recommended installation configuration.

**Figure 4 Recommended configuration**

## Section 5 Maintenance

### ⚠ DANGER

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

### ⚠ DANGER

Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.

### ⚠ WARNING

Multiple hazards. The technician must make sure that the equipment operates safely and correctly after maintenance procedures.

### NOTICE

Do not disassemble the instrument for maintenance. If the internal components must be cleaned or repaired, contact the manufacturer.

## 5.1 Maintenance schedule

Table 2 shows the recommended schedule of maintenance tasks. Refer to Table 2 to identify the maintenance items. Facility requirements and operating conditions may increase the frequency of some tasks.

Table 2 Maintenance schedule

Task	as necessary	3 months	1 year
Clean the instrument on page 15	X	X	
Clean spills on page 16	X		
Clean the probe on page 16	X		
Service maintenance (contact technical support)			X

## 5.2 Clean the instrument

### NOTICE

Never use flammable or corrosive solvents to clean any part of the instrument. Use of these solvents may degrade the environmental protection of the unit and may void the warranty.

1. Make sure that the probe holder is in the upper position. Refer to the TM1000 software.
2. Remove the probes and tubes from the holder. Clean the holder with a moist cloth if necessary.
3. Set the instrument to off.
4. Remove all of the sample beakers.
5. Remove the tray. Use tap water to clean the tray. Remove all stains and spills.
6. Disconnect all cables.
7. Remove the rotating module and then the lifting module from the base.
8. Clean the exterior surface with a moist cloth or with a mixture of water and mild detergent. Dry with a soft cloth.
9. Assemble the instrument. Refer to [Assemble the instrument](#) on page 10.

### 5.3 Clean spills

#### ▲ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

1. Obey all facility safety protocols for spill control.
2. Discard the waste according to applicable regulations.

### 5.4 Clean the probe

Refer to the probe documentation for information about the probe maintenance.

## Section 6 Troubleshooting

Problem	Possible cause	Solution
The blue indicator LED is off.	The tray was moved manually.	Start the instrument or send an order from the software.
The tray does not stop at position 1 or at any specified position.	The sampler does not sense the beaker or position 1 of the tray.	Restart the instrument. If the rod does not move to the upper position and the tray stops at position 1, contact technical support.
The sensor holder does not complete the vertical movement.	There is an internal damage	Contact technical support.
The beakers are not aligned. The sensors are not in the correct position. There is too much noise when the sampler operates.	The tray or one of the modules was not correctly installed.	Make sure that all instrument components are correctly installed.
The sampler does not operate.	There is a communication failure between the instrument and the PC.	Examine the instrument and instrument cables for damage. Make sure that the cables are correctly connected. Restart the instrument. If damage is found or the problem continues, contact technical support.
The green indicator LED is off.	Power is not supplied to the instrument.	Make sure that power is supplied to the instrument and that the instrument is set to on. If there is power and the instrument does not start contact technical support.

## Section 7 Replacement parts and accessories

#### ▲ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer. Use only the approved parts.

**Note:** Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Description	Item no.
Set of beakers, 50 mL (10x)	LZE108
Set of beakers, 150 mL (10x)	LZE109

## Section 7 Replacement parts and accessories (continued)

Description	Item no.
Set of beakers, 50 mL (20x)	LZE193
Set of beakers, 90 mL (30x)	LZE194
Holder, one probe	LZE191
Holder, two probes	LZE192
Magnetic stir bar, PTFE coated, 6 x 20 mm (10x)	LZE136
Magnetic stir bar, PTFE coated, 12 x 4.5 mm (30x)	LZE217
TM1000 Software with Ethernet cables	LZE127
Tube holder, 4 positions	LZE141
RS232 cable with a USB adapter for TitraLab 1000 Series	LZE195
Ring adapter for 50-mL beaker	LZE213
Power supply (without power cable)	LZE144
Power cable (EU)	YAA080
Power cable (US)	XLH055
Power cable (UK)	XLH057

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)