
Lachat QuikChem In-Line Module

05/2017, Edition 3

User Manual
Bedienungsanleitung
Manuel d'utilisation
Manual del usuario
Manuale utente

English.....	3
Deutsch.....	19
Français.....	35
Español.....	51
Italiano.....	67

Table of contents

Specifications on page 3	Operation on page 11
General information on page 3	Maintenance on page 12
Installation on page 6	Troubleshooting on page 16
User interface and navigation on page 10	Replacement parts and accessories on page 17

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	23 x 54 x 14 cm (9 x 21.3 x 5.5 in.)
Weight	4.5 kg (10 lb)
Power requirement	115 VAC, 60 Hz, 1.2 A
	230 VAC, 50 Hz, 0.8 A
Fuse	1.6 A, 250 VAC, slow blow
Operating temperature	5 to 40 °C (41 to 104 °F)
Pollution degree	2
Protection class	1
Installation category	II
Storage/operating humidity	Maximum relative humidity is 80% for temperatures up to 31 °C (87.8 °F), decreases linearly to 50% relative humidity at 40 °C (104 °F)
Location	Indoor or laboratory

Specification	Details
UV lamp ¹ —(most common)	354 nm (dark blue), Ushio type F6T5BLB (50717). Use with the cyanide determination.
	254 nm (clear), G6T5 (50775). Use with the total nitrogen or the total phosphorous determination.
Certifications	CE mark; listed by ETL to UL and CSA safety standards (cETLus) for North America.
Warranty	1 year

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

¹ Refer to the specific QuikChem method for lamp specifications for an application.

Use of hazard information

⚠ DANGER
Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.
⚠ WARNING
Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.
⚠ CAUTION
Indicates a potentially hazardous situation that may result in minor or moderate injury.
NOTICE
Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.

	This symbol indicates the need for protective eye wear.
	This symbol, when noted on the product, identifies the location of a fuse or current limiting device.
	This symbol identifies a risk of chemical harm and indicates that only individuals qualified and trained to work with chemicals should handle chemicals or perform maintenance on chemical delivery systems associated with the equipment.
	This symbol identifies the presence of a strong corrosive or other hazardous substance and a risk of chemical harm. Only individuals qualified and trained to work with chemicals should handle chemicals or perform maintenance on chemical delivery systems associated with the equipment.
	This symbol indicates that the marked item can be hot and should not be touched without care.
	This symbol indicates the presence of devices sensitive to Electrostatic Discharge (ESD) and indicates that care must be taken to prevent damage with the equipment.

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Electrostatic discharge (ESD) considerations

NOTICE	
	Potential Instrument Damage. Delicate internal electronic components can be damaged by static electricity, resulting in degraded performance or eventual failure.


Refer to the steps in this procedure to prevent ESD damage to the instrument:

- Touch an earth-grounded metal surface such as the chassis of an instrument, a metal conduit or pipe to discharge static electricity from the body.
- Avoid excessive movement. Transport static-sensitive components in anti-static containers or packages.
- Wear a wrist strap connected by a wire to earth ground.
- Work in a static-safe area with anti-static floor pads and work bench pads.

Product overview

This instrument contains a UV lamp and a heater block that prepare a sample for an analyzer. This instrument is for indoor use only. Refer to [Figure 1](#).

Figure 1 Product overview


1 Display panel	5 Power switch	9 Thumbscrews
2 UV lamp	6 Heater block	10 Access cover
3 Thermocouple	7 Power connector	11 UV lamp indicator
4 Heat rod	8 Fuse box	12 Sample in/out ports

Product components

Make sure that all components have been received. Refer to [Figure 2](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 2 Instrument components

1 Instrument	3 Clips with unions
2 Power cord	4 Heater block

Installation





▲ WARNING

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Install the heater block

Refer to the illustrated steps that follow.

Instrument placement

⚠ CAUTION

Ozone inhalation hazard. Under certain conditions, this instrument produces ozone concentrations above safe exposure limits. Plumb waste gases to a fume hood or to the building exterior in accordance with local, regional and national requirements.

This instrument is rated for an altitude of 2000 m (6562 ft) maximum. Use of this instrument at an altitude higher than 2000 m can slightly increase the potential for the insulation to breakdown, which can result in an electric shock hazard. The manufacturer recommends that users with concerns contact technical support.


Put the instrument in a position that makes it easy to disconnect the instrument from the power source.

Put the instrument in a position so that the fan ventilation opening on the left side is not blocked. Do not allow the power cord to contact hot surfaces on the QC8500 measurement channels.

There are two placement options:

- Put the instrument on top of the QuikChem System Core ([Figure 3](#)). Make sure that the heater controller is on the left side of the system.
- Put the instrument below the peristaltic pump if the space on top of the QuikChem system is not available. The peristaltic pump is located on the left side of the QuikChem system.

Figure 3 Coupler installation

1 Analyzer	3 Sample outlet tubes
2 Sample tubes to analyzer	4 Sample inlet tubes

Plumbing

Make sure to use the specified tubing size. Refer to the QuikChem method for details.

Prepare the reagent and standards

1. Refer to the QuikChem method for details about the reagent and standards preparation.
2. Read all system notes.
3. Connect the manifold as shown in the manifold diagram of the QuikChem method.

Connect the power cord

Items to collect:

- Surge-protected power strip
1. Connect the power cord to the instrument. The power cord must be rated for the supply voltage and current.
 2. Connect the power cord to a power strip.
 3. Turn the power strip on.
 4. Turn the instrument power switch on.

Start the peristaltic pump

Refer to the pump operation section in the QuikChem user manual for details on pump tube installation and peristaltic pump operation.


Look for leaks

1. After the manifold is installed and all pump tubes are securely attached, put all reagent lines in deionized water.
2. Make sure that there are no leaks on the manifold, the injection valve or the instrument.

User interface and navigation

Figure 4 shows the instrument display and keypad. Table 1 gives the key functions.

Figure 4 Display and keypad

1 PV LED display	3 UP arrow	5 Show/change
2 SV LED display	4 DOWN arrow	6 Enter

Table 1 Key functions

Key	Function
↩	Saves the change.
↻	Shows or changes the temperature setting.
▼	Decreases the value.
▲	Increases the value.

Operation

▲ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

Set the temperature

Note: Do not make the temperature higher than 60 °C (140 °F) without water or reagents in the tubing.

1. Push  until the display shows SP (set point).
2. Push the **UP** or **DOWN** arrows to raise or lower the set point. The display continues to show the SP mode.
3. When the temperature is set, push  to save the information in memory.
4. Push  once. The display shows the current temperature in Celsius. A red light shows the heater is on.
5. Let the temperature stabilize.

Find the correct timing

To find the correct timing in the Omnion FIA:

Items to collect:

- Stop watch or other timer
- Dye (provided with the system) - red dye gives the best results

Note: Due to differences in the viscosity and boiling point of each reagent, do not use water in place of reagents during the dye test.

1. Put some dye in a test tube.
2. Create a method in Omnion. Enter all valve timing parameters (cycle period, sample reaches first valve, load period and inject period) and sampler timing parameters.
3. Run the dye as a sample. Start the timer when the probe goes into the test tube.

4. Follow the dye through the system. When the sample gets to Port 6 at the injection valve, stop the timer.
5. Record the time, which is the sample reaches first valve time. Make sure to change this time in the valve timing of the method.
6. Save the method and run the dye again to measure the next parameter.
7. Run the dye and follow it through the system. When the dye gets to Port 6, start the timer. This is the instant when the valve toggles for the first time.
8. When the sample slug gets to Port 6, stop the timer.
9. Record this time and subtract 5 seconds to get the load period that is specified in the valve timing of the Omnion method.
 - The cycle period is specified in the method (cycle period = load period + inject period).
 - The inject period is automatically set by the software once the cycle period and load period are defined.

Calibrate the instrument

1. Put all reagent lines in their respective reagent bottles.
2. Put the standard solutions into the standard vials.
3. Refer to the analyzer user manual to start a calibration run.

Brackish timing in Omnion FIA

Some of the in-line methods make Brackish (or manual) integration timing necessary. The QuikChem method shows when the brackish timing is needed.

To learn more about how to set brackish timing, refer to the software manual or search for the keywords graphical events programming or brackish events in the Omnion help system.

Tips

- Do not use the stand-by speed of the peristaltic pump.
- Put reagent lines in deionized water to avoid wasting reagents.

- Make sure that the deionized water flows through the wash reservoir.
- The sample line is usually not a green-green pump tube. Refer to the QuikChem method for the correct pump tube.
- Set the peristaltic pump to 35.
- Get to the required temperature before a run is started.
- Allow the reagents to flow at normal speed through the entire system for at least 15 minutes before a run is started.
- Change the bubble trap membrane once a month or when the precision is poor or air bubbles are a constant problem.
- Look at the method support data in the QuikChem method and compare the data to the results. Make sure that the peak height and area are similar for standards of the same concentration.
- If the instrument and other standard FIA channels run simultaneously, the sample line for the standard channels can be changed from green-green to orange-orange. Cut the pump tube like any other sample line pump tube. Use a dye solution to make sure that the timing is correct. This change uses less sample.

Maintenance

⚠ DANGER	
	Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.
⚠ DANGER	
	Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.
⚠ WARNING	
 	Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

Maintenance schedule

Table 2 shows the recommended schedule of maintenance tasks. Facility requirements and operating conditions may increase the frequency of some tasks.

Table 2 Maintenance schedule

Task	Daily	Annually	As necessary
Clean the instrument on page 12	X		
Replace a fuse on page 13			X
Prepare the tubing on page 13		X	X
Replace the lamp on page 14			X
Replace the heater block on page 15			X

For the maintenance tasks specific to an individual method, refer to the QuikChem method.

Clean the instrument

NOTICE
Do not use any cleaning or decontamination procedures other than those recommended by the manufacturer. Other procedures can damage the instrument. Contact technical support to make sure that the procedures do not damage the instrument.

1. Set the instrument to off. Remove power from the instrument.
2. Clean up any spills before the instrument is cleaned.
 - a. Use baking soda to neutralize acidic spills.
 - b. Use water to remove the baking soda solution.
 - c. Dry the area with a lint-free cloth.
3. Clean all surfaces of the instrument with water.
4. Dry all surfaces of the instrument with a lint-free cloth.

Replace a fuse



⚠ DANGER

Fire hazard. Use the same type and current rating to replace fuses.

To replace a fuse, refer to the illustrated steps in [Figure 5](#).


Figure 5 Fuse replacement

Prepare the tubing

To prepare the tubing, refer to the QuikChem method and [Figure 6](#).

Figure 6 Tubing preparation

1 Tubing	3 Thermal tape, heater	5 Start and stop marks
2 Thermal tape, high temperature	4 Tubing for other connections ²	6 Tubing winds around the lamp/heater block

Replace the lamp

⚠ DANGER	
	Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.

² For lengths and other information, refer to the QuikChem method.

⚠ CAUTION

Burn hazard. Turn off power and allow the instrument to cool before this procedure.


Never disconnect the lamp when the power is turned on.

Disconnect the AC input before any part of the instrument is opened.

If the lamp has been recently used, let the instrument temperature cool. The time required depends upon the temperature set point.

1. Remove all power to the instrument.
2. Remove the top right cover panel.
3. Disconnect the lamp tubing.
4. Remove the lamp from the brackets.
5. Remove the aluminum thermal tape and tubing from the lamp. If the thermal tape is still in good condition after it is removed, use it on the new part.
6. To prepare the tubing, refer to [Prepare the tubing](#) on page 13.
7. To add tubing around the new lamp, refer to [Figure 7](#).
8. Install the lamp.
9. Connect the tubing.
10. Install the cover panel and secure it.
11. Apply power to the instrument.

Figure 7 Tubing replacement on the lamp

Replace the heater block

⚠ DANGER

Electrocution hazard. Remove power from the instrument before doing maintenance or service activities.

⚠ CAUTION

Burn hazard. Turn off power and allow the instrument to cool before this procedure.


Never disconnect the thermocouple and heat stick when the power is turned on. The result is uncontrolled heating.

Disconnect the AC input before any part of the instrument is opened.

If the heater has been recently used, let the instrument temperature cool. The time required depends upon the temperature set point.

1. Remove all power to the instrument.
2. Remove the top right cover panel.
3. Disconnect the heater block tubing.
4. Refer to the QuikChem method to make sure that the type and amount of tubing required for an application is correct.
5. To add tubing around the heater block, refer to [Prepare the tubing](#) on page 13 and [Figure 8](#).
6. To install the heater block, refer to [Install the heater block](#) on page 7.

Figure 8 Tubing replacement on the heater block

Troubleshooting

⚠ CAUTION

Burn hazard. Obey safe handling protocols during contact with hot sample components.

Decrease the temperature on the instrument to 70 °C (158 °F) or less before the pump tubes are examined. The typical amount of time to decrease the temperature from 140 °C (284 °F) to 70 °C (158 °F) is approximately 10 minutes.

General troubleshooting considerations

The in-line manifolds are different than the standard FIA manifolds. The main differences are:

- Sample line—The sample line used is usually not a green-green pump tube. Smaller pump tubes are chosen so that the sample (mixed with a digestion solution) is put on the heater and on the UV lamp for a longer period of time.
- Pre-valve sample preparation—The sample line from the sampler does not go directly to Port 6 of the injection valve. The sample line mixes with a digestion solution, gets heated, gets UV irradiated, then goes through the diffusion block or bubble tap membrane, and travels to Port 6 of the injection valve.

Table 3 shows the most common problems encountered:

Table 3 Troubleshooting guide

Problem	Solution
Air spikes in the FIAgram	Make sure that the timings agree with the information in Find the correct timing on page 11.
	Do a degass procedure on the reagents.
	Replace the membrane in the dialysis block or replace the bubble trap membrane.

Table 3 Troubleshooting guide (continued)

Problem	Solution
Poor recovery of a QC sample that consists of complex-cyanide	Make sure that the reagents are made fresh.
	Make sure that sample flows through all pump lines.
	Make sure that the temperature agrees with the QuikChem method.
	Look at the small hole on the top of the instrument to see if the UV lamp is turned on. If not, replace the lamp.
Poor or no recovery of a QC sample	Examine the pump tubes for flow.
	Examine the temperature setting.
	Make sure that the reagents are made correctly.
	For the cyanide determination, make sure that the acid line is put in the digestion solution. Examine the pump tubes for the acid, the sample and the acceptor.
	Read the QuikChem method notes to learn more about the expected recoveries.

Table 3 Troubleshooting guide (continued)

Problem	Solution
Air flow though the diffusion block gets very fast, then slows almost to a stop	Make sure that there are no blockages or crimps in the reagent and samples tubing.
	Make sure that the waste tubing flows freely. If not, there is excessive back pressure.
	Make sure that the instrument waste goes into a waste container that is applicable for that solution. Make sure that the tubing is above the liquid level and that the container is vented to prevent pressure inside the container. To prevent cavitation and vapor lock, the waste container must be at the same level as the instrument.
Imprecision when running replicates	Examine the reagent and water bath tubing for plugs or crimps. Small dirt particles or algae from a dirty water bath causes flow problems.
	Examine the timings in Find the correct timing on page 11.
	Replace the bubble trap membrane.
	Make sure that all pump tubes are in good condition and that there is good liquid flow.

Replacement parts and accessories

▲ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
Tubing, heater block	50028
Tubing, lamp	50728
UV lamp, 6 W	— ³
UV lamp, 8 W	— ³
Heater block only	25008 ⁴
Thermal tape, heater block	31182
Thermal tape, high temperature, 4-inch	31184
Fuse, 1.6 A, 250 V, slow blow	20205

³ Refer to the QuikChem method, and then contact technical support.

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)