
# ***TitraMaster 85***

Getting started  
Prise en main

# Contents

<i>Introduction</i> .....	3
<i>Notes to PC Administrator</i> .....	5
<i>Minimum configuration recommended</i> .....	6
<i>Installing TitraMaster 85</i> .....	8
<i>Getting help</i> .....	11
<i>Starting up</i> .....	12
<i>Main window</i> .....	16
<i>TitraMaster 85 icons</i> .....	18
<i>Quick tour</i> .....	20
<i>Uninstalling</i> .....	20

# Sommaire

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# Introduction

The TitraMaster 85 is a Windows-based PC program dedicated for use with any of the new range of workstations from Radiometer Analytical.

TitraMaster 85 exists in 4 versions:

- TitraMaster 85 Data Collector,
- TitraMaster 85 Data Collector FDA21CFR/11,
- TitraMaster 85 Bi-directional,
- TitraMaster 85 Bi-directional FDA21CFR/11.

With the **Data collector function**, the analysis is created and edited from the workstation. The analysis can be started either from the workstation or from TitraMaster 85. Curves and results are collected by TitraMaster 85. The Data collector version of TitraMaster 85 can be used with TitraLab® or MeterLab® compatible workstations. Concerning the ION450 and ION570, all methods can be run and the curves collected by TitraMaster 85, except the EC/pH Measurement method (pH and conductivity measurements performed simultaneously).

With the **Bi-directional function**, the application is created, edited and started from the PC and TitraMaster 85. On starting an application, this application is exported from the PC to the workstation. Then, TitraMaster 85 works as a Data collector version. Curves and results are collected by the software. The Bi-directional version of TitraMaster 85 can be used with TitraLab® or MeterLab® compatible workstations.

With the ION450, ION570, TIM840 and TIM845 workstations, TitraMaster 85 runs as a Data collector (the analysis is created and edited from the workstation only but you can run it either from the workstation or from the PC and the "Run" button of TitraMaster 85).

TitraMaster 85 automatically records all data obtained on up to 7 units interfaced via a single standard serial or USB PC com port. TitraMaster 85 is a powerful data archiving tool, which enables you to access your stored results and curves in a flash via user-programmable requests. You can consult the dedicated analytical database and filter and sort the archived results by date, method type etc. providing the information you need at a glance.

The TitraMaster 85 **FDA21CFR/11 version** offers all the necessary features to make your laboratory fulfil its obligations regarding the FDA21CFR/11 standard requested by the Food and Drug Administration. It means:

- An encrypted login and password of at least 6 characters each saved in a secure database.
- 3 hierarchical levels of user profiles.
- Fully traceable of results records.
- Totally secure result database for both archiving and restoring operations.
- Password protection before validation of any changes.
- Electronic signature for record certification and application certifications.
- Export of signed records in read-only format.
- Event viewer reporting all operations performed by operators using the software.

This Getting started instructions are common to all TitraMaster 85 versions. If you are using a FDA21CFR/11 version, please read the important chapter starting on the next page, otherwise skip to chapter “Minimum configuration recommended”.

Take advantage of new features and product updates by visiting the Radiometer Analytical [www.titration.com](http://www.titration.com) Resource Centre.

# Notes to PC Administrator

*This chapter deals with the FDA21CFR/11 versions of TitraMaster 85.*

Under Windows 2000 Pro or XP Pro, all the TitraMaster 85 application files are secured automatically when installing the application and this for any user profiles apart from the PC administrator or a user "With power". Therefore, it is required to use one of these two recent operating systems.

The PC Administrator must follow some rules in order to totally secure the TitraMaster 85 installation directory. If these rules are fulfilled, the PC Administrator will be the only person able to modify, handle or delete the TitraMaster 85 files.

1. The PC Administrator installs TitraMaster 85 on a local disk with a NTFS file system.
2. The PC Administrator decides which users can use TitraMaster 85 then creates a "User" account for each user in Windows. The "User" level means a person that can only read the files within the Windows Operating System.
3. Before starting TitraMaster 85, logon in Windows to one of the users created in step 2.  
To XP Pro Windows users:  
Do not open more than one TitraMaster 85 application as 2 TitraMaster 85 applications cannot run simultaneously on the same computer.
4. The directories or discs used for saving or exporting TitraMaster 85 files must be configured as read-only directories or discs (and must have an NTFS file system).

## **Note**

TitraMaster 85 with network accesses (examples: database archiving and analyse exportation to a network disk,...) must be configured by the PC Administrator.  
Contact us at Radiometer Analytical SAS for more information.

# Minimum configuration recommended

## 1 workstation only

Operating system: W2000 Pro (SP2) or XP Pro (SP2).

Internet Explorer version 6.0 or higher must be installed on your PC.

- Pentium 1 GHz processor or equivalent, 256 MB RAM,
- SVGA 1024 x 768 screen resolution, 17" colour screen display,
- 1 available serial port or 1 available USB port with 1 USB to serial converter,
- 1 CD-ROM drive,
- 1 PS2 or USB mouse type,
- 1 Windows compatible printer (optional),
- Hard disk, 150 MB of free space for software installation, 1 GB required for proper use with your database,
- Excel 97 (or better) installed on the PC to use Export to Excel function and table & customised printing formats.

## 2 to 7 workstations

Operating system: W2000 Pro (SP2) or XP Pro (SP2).

Internet Explorer version 6.0 or higher must be installed on your PC.

- Pentium 2 GHz processor or equivalent, 512 MB RAM,
- SVGA 1024 x 768 screen resolution, 19" colour screen display,
- 1 available serial port or 1 available USB port with 1 USB to serial converter,
- 1 CD-ROM drive,
- 1 PS2 or USB mouse type,
- 1 Windows compatible printer,
- Hard disk, 150 MB of free space for software installation, 1 GB required for proper use with your database,
- Excel 97 (or better) installed on the PC to use Export to Excel function and table & customised printing formats.

## TitraMaster 85 FDA21CFR/11 versions, 1 to 7 workstations

Same configuration as for 2 to 7 workstations with:

- **NTFS file system** required.

# Installing TitraMaster 85

## Update version from V1.x.x

If you want to install an update version of TitraMaster 85, first backup your database, uninstall TitraMaster 85 and remove manually the corresponding installation directory. Install the TitraMaster 85 full version then restore the database.

Proceed as follows to backup your files:

1. In TitraMaster 85, select **File > General setup**.
2. Click the **Database** tab then click the **Archive** icon.
3. In the **Save As** window, select the backup file destination folder and type a name for the backup file (DBZ file).
4. Click **Save** to close the **Save As** window.
5. Click **OK** to close the **General setup** window.  
Quit TitraMaster 85.

*To restore a database, consult the TitraMaster 85 help, topics "Handling the result database/Restoring a database".*

## Update version from V2.x.x

Use the TitraMaster 85 Upgrade CD-ROM or a file downloaded from the TitraLab Resource Centre. It will automatically upgrade your installed version.

We recommend you to backup your database before upgrading TitraMaster 85.

## Before starting installation

Refer to system requirements and make sure that your PC complies with the minimum configuration recommended (see previous chapter).

### **Internet Explorer 6.0 or higher must be installed on your PC.**

If you have a French or an English Internet Explorer, Internet Explorer 6.0 can be installed from the TitraMaster 85 CD-ROM as follows:

1. In Windows, click *Start\Settings\Control panel\Add/Remove Programs*.
2. Click the **Add new Programs** then **CD or Floppy** buttons.
3. Browse the CD-ROM to search for this file:  
/ie6setup.exe (Internet Explorer in English)  
/ie6fr/ie6setup.exe (Internet Explorer in French)
4. Close all application windows.
5. Run the ie6setup.exe file.
6. At end of installation, switch the computer Off then On.

### **For TitraMaster 85 FDA21CFR/11 users:**

TitraMaster 85 is installed by the PC Administrator on a local disk with NTFS file system.

TitraMaster 85 must be installed into a directory of the local disc and not into the root directory. Select preferably the installation directory proposed by default.

The TitraMaster 85 installation program creates in a Windows user group, a user account called **TM85USER** having the “With power” access rules. Do not modify the access rules of this user. Do not delete this user: if you delete **TM85USER**, you will have to uninstall then re-install TitraMaster 85.

**TM85USER** is using the TitraMaster 85 application to access disc in Read/Write mode irrespective of the user’s rights of the operator logged on Windows. If you (PC administrator) need to give additional rights to **TM85USER** (example: access to a network printer), please contact Radiometer Analytical SAS.

# Installation procedure

1. Close all applications and disable your antivirus application.
2. Insert the installation CD into a CD drive. Installation starts automatically (AUTORUN). If it does not, use your Windows Explorer to locate Setup.exe and run it.
3. Select the language of the installation software.
4. Click **OK** (information message about languages) then **Next**. Installation starts with the **Welcome** dialogue box. Click **Next**.
5. Accept the License agreement and click the **Next** button.
6. In the **Customer information** dialogue box, enter your name, the name of your company and the TitraMaster 85 serial number indicated on a label in the TitraMaster 85 installation CD box (xxxVxxxNxxx). Enter this number correctly as you will be asked for it, together with the serial number of the workstation if you call us for servicing purposes. Click **Next**.
7. Enter the software access code (4 sets of 5 characters). This code is found on a label in the TitraMaster 85 installation CD box. Click **Next** then **OK**.

## **Important**

Keep carefully your serial number and access code in a safe place. If you loose them, you will be not able to re-install the software!

8. Review your installation settings. Accept the installation directory proposed by default clicking **Next** or change it clicking the **Change** button.
9. Click the **Next** button. The setup program creates the destination directory, transfers the files into the directory and to the C:\Windows\System directory and creates an icon for starting the TitraMaster 85 in the **Start** folder.
10. Click the **OK** button (information message about the embedded software of the workstation).
11. Click the **Finish** button at the end of installation. TitraMaster 85 will be installed in French if Windows is French and in English in all other cases.

12. Keep the installation CD into the CD drive and start TitraMaster 85 to end the installation.  
Installation is complete, you can remove the installation CD.

**Important:**

Make sure that you have downloaded the latest version of the workstation (and ABU52 or ABU65 burette if relevant) embedded software. Use FlashLoader Light to download this software into the workstation.

Flashloader Light and the workstation embedded software are available on your TitraMaster85 CD-ROM.

If you change your Windows operating system on your PC, you must uninstall then reinstall TitraMaster 85.

## Getting help

Click  to gain access to the TitraMaster 85 on line help.


You can also press F1 or select the command **? > Contents** in the TitraMaster 85 window to display the contents menu of the help file. From here, select the item of interest or search for a certain topic.

# Starting up

## Connecting TIM/ION workstations

Connect the **PC/PRINTER** socket on the workstation to a 9-pin male serial port COM1 on the PC using cable part no. A95X501. You can also connect the **PC/PRINTER** socket on the workstation to an USB port on the PC using cable part no. A95X501 and an USB to serial port converter.


PC serial or USB port

**Note:** *If you are using a USB to serial Adapter, Windows assigns by default the USB port to COM5 to COM8. You have to enter your configuration panel to manually configure your adapter COM port setup to have TitraMaster 85 connected to COM1 or COM4 (and not above COM4).*

If necessary, connect the printer to a parallel or USB port on the PC. Connect the workstation to the mains using the line cord part no. A95S001 (230 V) or A95S002 (115 V). Switch on the workstation (**O/I** switch).

**Note:** *You can connect up to 7 workstations to a PC serial or USB port. The workstations are interconnected as follows:*

PC serial or USB port

**Important:** RS232 cable length must not exceed 20 m (length from the PC to the last workstation connected).

# Starting TitraMaster 85 for the first time

1. Select **Start > TitraMaster 85** and select Supervisor.

2. Type "default" (use lower case letters).


3. Type "password" (use lower case letters).

4. Click the **Change password** button.

Type a new login (6 characters min.).

Type a new password (6 characters min.).

Re-enter this password and click **OK** (\*).

5. If you have connected on switched on your workstation, check that this option is selected (\*\*).


## (\*) **Concerning passwords**

In TitraMaster 85, any user, irrespective of his profile, must change password every at least 90 days. Passwords cannot be reused.

## (\*\*) **Note to portable PC users only.**

At first start, uncheck this option (no "✓" must be displayed) then select the **File > Scan USB -> RS232 ports** command from the TitraMaster 85 main window in the menu bar. Then refresh the workstation list.


6. TitraMaster 85 scans the PC's serial and USB ports and detects the workstations connected and on-line.

Workstation(s) detected.  
TitraMaster 85 assigns a name to each one. "V05.1" is the software embedded version.

In the "User" box, the name of the user connected is displayed.

7. Click  then  to create a new user account at the Supervisor level.

### Important


The Supervisor is the only person who can create and unlock user accounts as logins and passwords are unique and strictly confidential. It is important to create a second Supervisor user account to ensure a replacement when the first Supervisor is absent.

For further information refer to the Help file, topic "Handling users - Creating a user account".

7. Select **File > General setup** to configure your printouts and backup file.  
For further information, refer to the Help file, topics "Printouts - Preparing printouts" and "Handling database - Archiving a database".

# Main window

1 **Function toolbar:** see next page.

2 **Workstation list:** shows the workstation connected to the PC. The status of the workstation is shown by a sunny, cloudy, stormy or question mark icon.

The term "Free" indicates that the workstation shows that the workstation is ready for operation.

3 **Run button (available in the Bi-directional version only):** click to run the application of the workstation selected. The application is selected in the **Application** tab of the **Run** window.

4 **Log on/Log off button:** Click to log off and allow a different user to use TitraMaster 85.

5 **Run window of the workstation selected:** displays the current analysis in progress on this workstation. Using TitraMaster 85 Bi-directional, displays the applications created for this workstation.

6 **Users bar:** to log off, select None from the list. To log back on, select your id from the bar and enter your login and password.

## The Function toolbar

**Log on users:** access the list of users. Log on to a user with or without automatic detection of the connected workstations.

**Handling users:** access the list of users. Create, edit or delete a user account. *Supervisor access only.*


**Handling workstations:** access the list of workstations. View workstation data and enter notes.

**Creating and editing sample preparations:** create or edit a sample preparation to be run by a SAC850/SAC950 sample changer. *Supervisor access only.*

**Creating and editing applications (available in the Bi-directional version only):** create or edit an application. *Supervisor access only.*

**Applications (available in the Bi-directional version only):** consult the dedicated application database, filter and sort the archived applications.

**Results:** consult the dedicated result database, filter and sort the archived results.

**Check:** consult the control cards of a set of pre-selected analyses.

**Help:** get on-line help.

**Quit:** quit TitraMaster 85.

# TitraMaster 85 icons

## General icons

Create new user or request to sort archived data.

Delete.

Print.

Display or hide information.

Enter an equation.

Open on-line help.

Enter edit mode to modify, e.g. the workstation name.

Modify the request name.

Recalculate a result.

Reset a result searching request.

Open the Traceability form of an analysis.

Display table in list or icon format.

Icon of the **Run** window. Open the sample stack of the workstation current sequence or the workstation application selected (\*).

Waiting icon: buffer result for a calibration cycle which is not completed.

Empty all the analyses from the **Run** window. Analyses are not deleted from the result database.

Copy the selected application and create a new application or create a new version of this application (\*).

Icon of the **Run** window. Export the application to a selected workstation (\*).

Archive the selected application to an .xml file (\*).

Restore an archived application from an .xml file (\*).

Icon of the **Run** window of a TIM580/585. Start the “Stand by” mode of the selected application(\*).

Icon of the **Run** window. Start calibrating an electrode of the selected application(\*).

Icon of the **Run** window. Start calibrating a reagent of the selected application(\*).

Analysis is running.

(\* Available in the TitraMaster 85 Bi-directional version only.

Application file (set of method or sequence of methods parameters) (\*).

Unlock a selected workstation or a user account (\*).

This application has been signed (certified) (\*)(\*\*).

Signing icon. To sign an analysis or an application (\*\*).

## Result indicator icons

Result accepted by the instrument and the user. Analysis or calibration cycle result is calculated with this test result.

Recalculated result.

This recalculated result has been signed (certified) (\*\*).

This analysis has been signed (certified) (\*\*).

Analysis stopped. Result not saved.

Result lies outside the preset acceptance limits. Test result rejected by instrument but accepted by user (Supervisor level only). Analysis or calibration cycle result is calculated with this test result.

Result accepted and no errors have occurred (green).

Result accepted, with an acceptance error (red).

Result accepted by instrument but rejected by user. Analysis or calibration cycle result is calculated without this test result.

Result calculated using an entered equation. Result lies within acceptance limits (green) or outside acceptance limits (red).

Result rejected by instrument and user. Analysis/calibration cycle result calculated without test result.

Error has occurred.

Internal procom error. No result found. Restart TitraMaster 85.

(\*) Available in TitraMaster 85 Bi-directional version only.

(\*\*) Available in TitraMaster 85 FDA21CFR/11 version only.

## Quick tour

The Quick tour has been developed to guide you through the initial use of TitraMaster 85. This nature of this tour depends on you version of TitraMaster 85 installed.

It is highly recommended to run this Quick tour to familiarise yourself with the TitraMaster 85 features.

Select the command **? > Tutorial** in the TitraMaster 85 window then click a demonstration topic to start the tour:

[Restoring the demonstration database](#)

[Creating an application](#)

[Running an application \(Bi-directional mode\)](#)

[Running an analysis \(Data collector mode\)](#)

[Viewing results/curves in deferred time](#)

[Reprocessing curves](#)

[Recalculating results](#)

[Handling the database](#)

[TitraMaster 85 and the FDA21CFR/11 Standard](#)

## Uninstalling

Close all applications and disable your Anti Virus application then select the **Start > Settings > Control Panel > Add/Remove Program** of Windows. Next select “TitraMaster 85 “ and click “Add/Remove”. All the files that have been created and copied by the TitraMaster 85 setup program are removed. The TitraMaster 85 icon placed in the “Start” folder is removed. The files created after the TitraMaster 85 installation are not affected by the uninstallation program (for example, the compressed database (\*.dbz files), the result database, the application database, the user list).

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