
DOC022.97.90252

sensION™ + EC71

03/2013, Edition 3

User Manual
Manuel d'utilisation
Manual del usuario
Manual do Usuário

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Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions	35 x 20 x 11 cm (13.78 x 7.87 x 4.33 in.)
Weight	1100 g (2.43 lb)
Meter enclosure	IP42
Power requirements (external)	100–240 V, 0.4 A, 47-63 Hz
Meter protection class	Class II
Storage temperature	–15 to +65 °C (5 to +149 °F)
Operating temperature	0 to 40 °C (41 to 104 °F)
Operating humidity	< 80% (non-condensing)
Measuring error (± 1 digit)	Conductivity: ≤ 0.1 %, Salinity/TDS: ≤ 0.5 %, temperature: ≤ 0.2 °C (0.36 °F)
Reproducibility (± 1 digit)	Conductivity: ± 0.1 %, Salinity/TDS: ± 0.1 %, temperature: ± 0.1 °C (0.18 °F)
Data storage	330 results and last 9 calibrations
Connections	Conductivity probe with built-in Pt1000 sensor: telephonic connector; magnetic stirrer: RCA connector RS232C for printer or PC: telephonic connector; external PC keyboard: mini DIN connector
Temperature correction	Manual, Pt 1000 temperature probe (A.T.C.), NTC 10 kΩ probe
Measurement display lock	Continuous measurement, by stability and by time
Display	Liquid crystal, backlit, 128 x 64 dots
Keyboard	PET with protective treatment
Certification	CE

General information

Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol, if noted on the instrument, will be included with a danger or caution statement in the manual.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/98/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.

Product overview


The sensION™+ meters are used with probes to measure various parameters in water.

The sensION™+ EC71 meter measures conductivity, salinity, TDS and temperature. Measurement data can be stored and transferred to a printer or PC.

Product components

Refer to [Figure 1](#) to make sure that all components have been received. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 1 Meter components





1 Standard solutions (147 $\mu\text{S}/\text{cm}$, 1413 $\mu\text{S}/\text{cm}$ and 12.88 mS/cm)	5 Rod with o-ring
2 Calibration beakers (with magnetic bar inside)	6 Probe (included with kits only)
3 Probe holder	7 Meter
4 Power supply	

Installation

Assemble the probe holder

Follow the numbered steps to assemble the probe holder and to connect the magnetic stirrer.

Connect to AC power

⚠ DANGER

Electrocution hazard. If this equipment is used outdoors or in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used to connect the equipment to its main power source.

The meter can be powered by AC power with the universal power adapter.

1. Select the correct adapter plug for the power outlet from the adapter kit.
2. Connect the universal power adapter to the meter (Figure 2).
3. Connect the universal power adapter to an AC receptacle (Figure 3).
4. Turn the meter on.

Figure 2 Connector panel




Figure 3 AC power connection

User interface and navigation

User interface

Keypad description


1 RETURN key: cancel or exit the current menu screen to the previous menu screen	5 DOWN key: scroll to other options, change a value
2 MEASUREMENT key: confirm the selected option	6 ON/OFF: turn on or turn off the meter
3 UP key: scroll to other options, change a value	7 LEFT key: change the measuring unit, enter numbers and letters
4 RIGHT key: change the measuring unit, enter numbers and letters	

Display description

The meter display shows the concentration, units, temperature, calibration status, operator ID, sample ID, date and time.

Figure 4 Single screen display

1 Sample ID	4 Sample temperature (°C or °F)
2 Measurement unit and value (conductivity, TDS or salinity)	5 Visual measurement timer
3 Measurement mode or time and date	

Navigation

Use the ↶ to return to the previous menu. Use the measure key ✓ to take a sample measurement or to confirm options. Use the arrow keys ▲▼ to scroll to other options or to change a value. To change the parameters use the arrow keys ◀ and ▶. Refer to each task for specific instructions.

Start-up

Turn the meter on and off

NOTICE
Make sure that the probe is connected to the meter before the meter is turned on.

Push ⏻ to turn on or turn off the meter. If the meter does not turn on, make sure that the AC power supply is properly connected to an electrical outlet.

Change the language

The display language is selected when the meter is powered on for the first time.

1. Use the ▲ or ▼ to select a language from a list.
2. Confirm with ✓. The measurement screen shows DATA OUTPUT.
3. Select Deactivated if no printer or PC is connected and confirm. Refer to [Select the data output](#) on page 9 for more information about Data Output.

Standard operation

Calibration

⚠ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

Calibration settings

The calibration settings contain Calibration type, Calibration frequency and Display options.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select the following options:

Option	Description
Calibration type	Calibration type—select Molar Standards, Demal Standards, NaCl Standards, Calibration to a X value, data introduction or Theoretical calibration. Refer to Calibration types for more information.

Option	Description
Cal. frequency	Calibration reminder—can be set between 0–99 days (default 15 days). The display shows the remaining time to the new calibration. Refer to Set the calibration reminder on page 9 for more information.
Display Ω	Display Resistivity—select YES to show simultaneously conductivity and resistivity.

Calibration types

Different calibration types can be selected.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select Calibration type.

Option	Description
Molar Standards	147 µS/cm, 1413 µS/cm, 12.88 mS/cm and 111.8 mS/cm at 25 °C (77 °F)
Demal Standards	1049 µS/cm, 12.85 mS/cm and 111.31 mS/cm at 25 °C (77 °F)
NaCl Standards	1014.9 µS/cm 25 °C (77 °F)
Calibration to a X value	To adjust manually any scale value of the measured conductivity.
Data introduction	Manual probe constant introduction.
Theoretical calibration	The probe calibration data is replaced with $C=1.000\text{ cm}^{-1}$.

Calibration procedure

This procedure is for general use with liquid calibration solutions. Refer to the documents that are included with each probe for additional information.

Note: Solutions should be stirred during calibration. For more information about the stirring settings, refer to [Change the stirring settings](#) on page 10.

1. Pour the buffers or calibration solutions into the labeled calibration beakers.
2. From the main menu use the ▲ or ▼ and ◀ and ▶ to select the CALIBRATION parameter. Confirm.
3. If required select the Operator ID (1 to 10) and confirm.
4. Rinse the probe with deionized water and put the probe into the first calibration beaker. Be sure that there are no air bubbles in the probe measuring chamber.
5. Push ✓ to Start calibration.
6. Push ✓ to measure the first calibration solution.
The next calibration solution is shown.
7. Rinse the probe with deionized water and put the probe into the second calibration beaker. Be sure that there are no air bubbles in the probe measuring chamber.
8. Push ✓ to measure the second calibration solution.
The next calibration solution is shown.
9. Rinse the probe with deionized water and put the probe into the third calibration beaker. Be sure that there are no air bubbles in the probe measuring chamber.
10. Push ✓ to measure the third calibration solution.
When the calibration is good, the display briefly shows Calibration OK and then returns to the main menu.
Note: When a printer is connected the print menu opens and the result can be printed.

View the calibration data

Data from the most recent calibration can be shown.

1. From the main menu push .
2. Use the ▼ to view the last calibration data. The meter returns automatically to the measurement screen after 3 seconds.

Set the calibration reminder

The calibration reminder can be set between 0 to 99 days (default 15 days). The display shows the remaining time to the new calibration.

Note: When 0 days is selected, the calibration reminder is turned off.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select Cal. frequency and confirm.
4. Use the ◀ and ▶ to advance to the next step and use the ▲ or ▼ to change a value. Confirm.
Push ✓ to start the calibration.

Advanced operation

Use a sample ID

The sample ID tag is used to associate readings with a particular sample location. If assigned, stored data will include this ID.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Sample ID and confirm.
3. Use the ▲ or ▼ to select

Option	Description
Automatic	A consecutive number will be automatically assigned to every sample.
Manual	A keyboard or a barcode scanner is required to enter the sample ID name before taking a measurement (maximum 15 characters).

Select the data output

Data can be stored or transferred to a printer or to a PC.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Data Output and confirm.

- Use the ▲ or ▼ to select

Option	Description
Deactivated	Select Deactivated if no printer or PC is connected.
For Printer	Select Dot matrix printer or Thermal printer.
For Computer	Select Terminal, LabCom or LabCom Easy. The LabCom Software controls several modules, pH and conductivity meters, automatic burettes, Samplers and so on from a computer. The LabCom Easy software gets pH and conductivity data from a PC.

Change the date and time

The date and time can be changed from the Date / Time menu.

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select Date / Time and confirm.
- Use the ◀ and ▶ to advance to the next step and use the ▲ or ▼ to change a value. Confirm.
The current date and time will be shown on the display.

Adjust the display contrast

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select Display contrast and confirm.
- Use the ◀ and ▶ to adjust the contrast of the display and confirm.

Adjust the temperature

The temperature measurement can be adjusted at 25 °C (77 °F) and/or 85 °C (185 °F) to increase accuracy.

- Put the probe and a reference thermometer in a container of water at approximately 25 °C and allow the temperature to stabilize.
- Compare the temperature read by the meter with that of the reference thermometer. The difference is the adjustment value for the meter.

Example: reference thermometer: 24.5 °C; meter: 24.3 °C.
Adjustment value: 0.2 °C.

- Enter the adjustment value for the 25 °C reading:
 - From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
 - Use the ▲ or ▼ to select Readjust temp. and confirm.
 - Use the ▲ or ▼ to select 25 °C and confirm.
 - Use the arrow keys to enter the adjustment value for 25 °C. Confirm.
- Put the probe and a reference thermometer in a container of water at approximately 85 °C and allow the temperature to stabilize.
- Compare the temperature from the meter with that of the reference thermometer. The difference is the adjustment value for the meter.
 - Use the ▲ or ▼ to select 85 °C and confirm.
 - Use the arrow keys to enter the adjustment value for 85 °C. Confirm.
 - Select Save changes and confirm.

Change the stirring settings

The magnetic stirrer can be turned on and the stirring speed can be changed in the Stirring menu.

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select Stirring and confirm.
- To turn on/off the stirrer push ✓.
- When the stirrer is turned on, use the ▲ or ▼ to change the stirring speed in %.

Note: Use the ▲ or ▼ to change the stirring speed during calibration and during a measurement.

Change the temperature units

The temperature units can be changed to Celsius or Fahrenheit.

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select Temperature units and confirm.

- Use the ▲ or ▼ to select between Celsius or Fahrenheit and confirm.

Data logger

Display data

The Display data log contains Measurement data, Electrode report and Calibration data. The stored data can be send to a printer or to a PC. When the data log becomes full (400 data points), the oldest data point is deleted when a new data point is added.

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select DATA LOGGER and confirm.
- Use the ▲ or ▼ to select Display data and confirm.
- Use the ▲ or ▼ to select

Option	Description
Measurement data	Measurement data—stores automatically each time a sample is measured
Electrode report	Electrode report—stores automatically the electrode history and measurement conditions
Calibration data	Calibration data—stores automatically the current calibration

Delete data

The entire measurement data or electrode report log can be deleted to remove data that has already been sent to a printer or PC.

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select DATA LOGGER and confirm.
- Use the ▲ or ▼ to select Erase and confirm.
- Use the ▲ or ▼ to select Measurement data or Electrode report and confirm. Confirm again to delete the data.
The entire log is deleted at once.

Send data to a printer or to a computer

NOTICE

The data output (printer or PC) needs to be selected first, so that the Print menu is available (refer to [Select the data output](#) on page 9).

Note: Refer to [Report output](#) on page 11 to select the report output type.

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select DATA LOGGER and confirm.
- Use the ▲ or ▼ to select Print and confirm. Select one of the following options and confirm with ✓ to print the data: Measurement data, Electrode data, Calibration data, Calibration report or Instrument condit.

Report output

NOTICE

The data output (printer or PC) needs to be selected first, so that the Type of report menu is available (refer to [Select the data output](#) on page 9).

Different report output types can be selected when a printer or a PC is connected.

- From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
- Use the ▲ or ▼ to select Type of report and confirm.
- When a printer or a computer and Terminal is connected, use the ▲ or ▼ to select

Option Description

Reduced Select Several or One sample as an output format

Option	Description
Standard	Select Several or One sample as an output format. Select several: Users: The user name appears on the printed report (17 characters). Header: The company name can be added as a header (40 characters) and appears on the print report. Identify sensor: The sensor model and the sensor serial number can be added and appears on the print report.
GLP	Select Several or One sample as an output format. Select several: Users: The user name appears on the printed report (17 characters). Header: The company name can be added as a header (40 characters) and appears on the print report. Identify sensor: The sensor model and the sensor serial number can be added and appears on the print report.

4. When a computer is connected and LabCom Easy (refer to [Select the data output](#) on page 9 for more information) is selected, use the ▲ or ▼ to select

Option	Description
Users	The user name appears on the printed report (17 characters).
Identify sensor	The sensor model and the sensor serial number can be added and appears on the print report.

Maintenance

⚠ CAUTION

Multiple hazards. Do not disassemble the instrument for maintenance or service. If the internal components must be cleaned or repaired, contact the manufacturer.

⚠ CAUTION

Personal injury hazard. Only qualified personnel should conduct the tasks described in this section of the manual.

Clean the instrument

NOTICE

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Clean the exterior of the instrument with a moist cloth and a mild soap solution.

Clean the probe

Clean the probe as needed. Refer to [Troubleshooting](#) on page 14 for more information about cleaning. Refer to the probe documentation for information about the probe maintenance.



Use the cleaning agents listed in [Table 1](#) for contaminations on the conductivity probe.

Table 1 Cleaning agents for the conductivity probe

Contamination	Cleaning agent
Grease, oils, fats	Electrode cleaning solution
Limescale	0.1 N HCl solution

Replace the magnetic stirrer

If the magnetic stirrer does not start, follow the numbered steps to replace the magnetic stirrer.

Troubleshooting

Refer to the following table for common problem messages or symptoms, possible causes and corrective actions.

Table 2 Calibration warnings and errors

Error/Warning	Solution
UNSTABLE READING Time t> 100 s	Calibrate again. Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol.
Check temperature. Check the cell.	Make sure the standard temperature is constant during the calibration. Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol.
STANDARD TEMPERATURE < 15 °C or > 35 °C	Calibrate again. The standard temperature must be between 15 and 35 °C.
Difference C > 30 %	Calibrate again.
C < 0.05 or > C > 50.00 °C	Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol; connect a different probe to verify if problem is with probe or meter.
SAME BUFFERS	Examine the standard solution: Use new standard solution.
Buffer or cell in poor conditions	
BUFFER NOT RECOGNIZED	
Difference C > 30 %	
SAME STANDARDS	

Table 3 Measurement warnings and errors

Error/Warning	Solution
23.2 °C 728 µS/cm @ 25 °C	Calibrate again.
EC out of range	Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol; connect a different probe to verify if problem is with probe or meter.
Out of range °C	Examine the temperature sensor. Connect a different probe to verify if problem is with probe or meter.
TC = 0 Impossible to measure salinity	Modify TC
Time > 150 s	Examine the temperature. Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol; connect a different probe to verify if problem is with probe or meter.

Replacement parts and accessories

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
sensION+ PH3 Lab pH-meter with accessories, without probe	LPV2000.97.0002
sensION+ PH31 Lab pH-meter, GLP, with accessories, without probe	LPV2100.97.0002

Replacement parts and accessories (continued)

Description	Item no.
sensION+ MM340 Lab pH & Ion-meter, GLP, 2 channels, with accessories, without probe	LPV2200.97.0002
sensION+ EC7 Lab conductivity meter, with accessories, without probe	LPV3010.97.0002
sensION+ EC71 Lab conductivity meter, GLP, with accessories, without probe	LPV3110.97.0002
sensION+ MM374, 2 channel Lab meter, GLP, accessories, without probes	LPV4110.97.0002

Consumables

Description	Item no.
Conductivity standard solution 147 µS/cm, 125 mL	LZW9701.99
Conductivity standard solution 1413 µS/cm, 125 mL	LZW9711.99
Conductivity standard 12.88 mS/cm, 125 mL	LZW9721.99
Conductivity standard 147 µS/cm, 250 mL	LZW9700.99
Conductivity standard 1413 µS/cm, 250 mL	LZW9710.99
Conductivity standard 12.88 mS/cm, 250 mL	LZW9720.99
Enzyme solution	2964349
Pepsin Cleaning Solution	2964349
Electrode cleaning solution	2965249
0.1 N HCl solution	1481253
Ethanol, 95% (Grease, oils, fats)	2378900

Accessories

Description	Item no.
Thermal printer, RS232, for sensION+ benchtop	LZW8201.99
Thermal paper for printer LZW8201, bag with 5 rolls	LZW9117.99
Dot-impact printer, RS232, for sensION+ instruments	LZW8200.99
Standard paper for printer LZW8200, bag with 10 rolls	LZW9000.99
Ribbon for printer LZW8200, 3 units	LZW9001.99
LabCom Easy PC SW, for sensION+ GLP, CD, cable, USB adapter	LZW8997.99
LabCom PC SW, for sensION+ GLP, CD, cable, USB adapter	LZW8999.99
RS232 cable, for sensION+ benchtop instruments	LZW9135.99
RS232 cable for sensION+ benchtop, with USB adapter	LZW9135USB.99
Magnetic stirrer with sensor holder, for sensION+ MM benchtop	LZW9319.99
3x50 mL printed beaker for benchtop conductivity calibration	LZW9111.99
Three-sensor holder, for sensION+ benchtop instruments	LZW9321.99
Holder and clamp for three sensors	LZW9155.99
Pyrex glass chamber, continuous flow measurements	LZW9118.99
PP protector, electrode storage	LZW9161.99

Standard solutions

Conductivity standard solutions

Refer to [Table 4](#) for conductivity values of standard solutions at varying temperatures.

Table 4 Conductivity and temperature values

Temperature		Conductivity (EC)			
°C	°F	µS/cm	µS/cm	mS/cm	mS/cm
15.0	59	119	1147	10.48	92.5
16.0	60.8	122	1173	10.72	94.4
17.0	62.6	125	1199	10.95	96.3
18.0	64.4	127	1225	11.19	98.2
19.0	66.2	130	1251	11.43	100.1
20.0	68	133	1278	11.67	102.1
21.0	69.8	136	1305	11.91	104.0
22.0	71.6	139	1332	12.15	105.4
23.0	73.4	142	1359	12.39	107.9
24.0	75.2	145	1386	12.64	109.8
25.0	77	147	1413	12.88	111.8
26.0	78.8	150	1440	13.13	113.8
27.0	80.6	153	1467	13.37	115.7
28.0	82.4	156	1494	13.62	—
29.0	84.2	159	1522	13.87	—
30.0	86	162	1549	14.12	—
31.0	87.8	165	1581	14.37	—
32.0	89.6	168	1609	14.62	—
33.0	91.4	171	1638	14.88	—
34.0	93.2	174	1667	15.13	—
35.0	95	177	1696	15.39	—

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)