M231Pt2 and M241Pt2-8

Double Platinum Electrodes


Electrodes Double Platine

Operating Instructions

Mode d'Emploi


M231Pt2

2 banana plugs
2 fiches banane

M241Pt2-8

BNC plug
Fiche BNC

View from the bottom
Vue de dessous

The design principles of the M231Pt2 and M241Pt2-8 Electrodes

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Double Platinum Electrodes

Operating Instructions

The M231Pt2 and M241Pt-8 are Double Platinum Electrodes ideal for Karl Fischer titrations and other imposed current potentiometric titrations.

Preparing for measurement

1. Remove the protection cap from the electrode.
2. Before starting a measurement, rinse the electrode with a suitable solvent.
3. Check that the electrode body is not damaged.

Cleaning and maintenance

Platinum wires

The sensing part of the electrode can be polished with finishing abrasive bands, part no. B90X550. While polishing, take great care not to bend the platinum wires.

In order to have reproducible and accurate results with the electrode, treatments are required. These treatments will depend on the last measurement performed:

Mild Cleaning

- Place the electrode overnight in RENOVO•N Normal Cleaning Solution.

Degreasing

- Remove mineral oils by rinsing the electrode with trichlorethylene and dabbing dry with a soft tissue.

Chemical cleaning

- The electrode can be cleaned with a sulphochromic solution, or with a boiling solution of hydrochloric acid.

Note:

If you are using an ethanol-based Karl Fischer solvent (example: E-Solvent), once a week, immerse the Pt-Pt electrode for 2 min. in a 10% v/v TritonX-100 solution. Then rinse with dry methanol and gently wipe. This treatment allows the electrode to recover full efficiency after a few minutes operation.

10% v/v TritonX-100 is available from reagent manufacturers or can be prepared by diluting 10 ml of TritonX-100 in 100 of deionised water.

Storage

- Replace the protection cap, and store the electrode in it's box at room temperature.

Accessories

Description	Part No.
BA03 Finishing Abrasive Bands	B90X550
○RENOVO•N Normal Cleaning Solution, 250 ml	S16M001
Plug Adapter 2 banana - BNC	A94P801
Electrode Plug Adapter, banana - Type 7	809-147

Specifications

Temperature range:	-10 to 100°C
Sensing element:	2 platinum wires
Diameter of Pt wire:	1 mm for M231Pt 0.5 mm for M241Pt-8
Distance between the Pt wires:	5 mm for M231Pt 4 mm for M241Pt-8
Connection:	2 x banana plugs for M231Pt 1 BNC for M241Pt-8

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