
DOC022.98.80554

# PTM450

03/2016, Edition 1

**User Manual**  
**Bedienungsanleitung**  
**Manuale utente**  
**Manuel d'utilisation**  
**Manual del usuario**  
**Manual do utilizador**  
**Návod k použití**  
**Brugsanvisning**  
**Gebbruikershandleiding**  
**Instrukcja obsługi**  
**Bruksanvisning**  
**Käyttäjän käsikirja**  
**Ръководство на потребителя**  
**Használati útmutató**  
**Manual de utilizare**  
**Kullanım Kılavuzu**  
**Návod na obsluhu**  
**Navodila za uporabo**  
**Korisnički priručnik**  
**Εγχειρίδιο χρήσης**

English .....	3
Deutsch .....	11
Italiano .....	19
Français .....	27
Español .....	35
Português .....	43
Čeština .....	51
Dansk .....	59
Nederlands .....	67
Polski .....	75
Svenska .....	83
Suomi .....	91
български .....	99
Magyar .....	107
Română .....	115
Türkçe .....	123
Slovenský jazyk .....	131
Slovenski .....	139
Hrvatski .....	147
Ελληνικά .....	155

## Table of contents

[Specifications](#) on page 3

[General information](#) on page 3

[Installation](#) on page 5

[Operation](#) on page 8

[Maintenance](#) on page 9

[Troubleshooting](#) on page 10

[Replacement parts and accessories](#) on page 10

## Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x H x D)	110 x 45 x 60 mm (4.33 x 1.77 x 2.36 in.)
Weight approximately	250 g (0.55 lb)
Operating temperature	15 to 35 °C (59 to 95 °F)
Storage temperature	-5 to 40 °C (23 to 104 °F)
Relative humidity	20 to 80%, non-condensing
Power requirements	<b>Instrument input:</b> 24 VDC 0.2 mA; <b>External plug-in power supply:</b> Input, 100–240 VAC, 50-60 Hz, 1.5 A, Class I; output, 24 VDC, 2.5 A
Output <sup>1</sup>	Minimum 20 mV; maximum 1250 mV
Impedance	10 kOhms
Wavelengths	3-position switch: 525 nm ±10 nm (green); 610 nm ±10 nm (red); OFF
Response time	< 1 second
Connections	Local BNC cable, 1.20 m (3.94 ft)
Certifications	<b>Safety:</b> IEC/EN 61010-1; <b>EMC:</b> IEC/EN 61326-1
Warranty	1 year (EU: 2 years)

## General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

## Safety information

### NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

<sup>1</sup> The signal output is not sensitive to daylight.

## Use of hazard information

<b>▲ DANGER</b>
Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.
<b>▲ WARNING</b>
Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.
<b>▲ CAUTION</b>
Indicates a potentially hazardous situation that may result in minor or moderate injury.
<b>NOTICE</b>
Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This is the safety alert symbol. Obey all safety messages that follow this symbol to avoid potential injury. If on the instrument, refer to the instruction manual for operation or safety information.
	This symbol indicates that a risk of electrical shock and/or electrocution exists.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.


## Product overview

The PTM450 is a photo-colorimetric titration module used with the OPT300 photo-colorimetric probe. Use the module and the probe to do photo-colorimetric titrations with the Titalab AT1000 series titrator.

Connect the module to the Titalab AT1000 series titrator with the legacy adapter and follow the application working procedure to complete a titration. Refer to the AT1000 series titrator and OPT3000 probe documentation for additional information.

Refer to [Figure 1](#).

**Figure 1 Product overview**


1 LED in connector (from probe)	5 BNC connector (to titration system)
2 LED out connector—525 nm (from probe)	6 External power supply connection
3 Titration module	7 LED out connector—610 nm (from probe)
4 3-position switch: Off, 525 nm, 610 nm	

## Product components

Make sure that all components have been received. Refer to [Figure 2](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

**Note:** The power adapter is not supplied with the instrument. The power adapter is sold separately. Refer to [Replacement parts and accessories](#) on page 10.

**Figure 2 Product components**

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## Installation

### ▲ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

## Installation guidelines

- This instrument is for indoor use only.
- The power supply connector on the rear panel must be easily accessible so the power can be disconnected quickly in case of emergency.
- Keep the instrument away from temperature extremes, including heaters, direct sunlight and other heat sources.

## Install the instrument

### ⚠ DANGER

Electrocution hazard. If this equipment is used outdoors or in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used for connecting the equipment to its main power source.

### ⚠ WARNING

Fire hazard. Use only the external power supply that is specified for this instrument.


### ⚠ CAUTION


Electrical shock and fire hazards. Make sure that the supplied cord and non-locking plug meet the applicable country code requirements.




Prepare the probe and the titration system. Refer to the probe and titration system documentation. Refer to the illustrated steps that follow.

**1**

**2**

# Operation

**⚠ DANGER**

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

## Do a titration

Do the steps that follow to complete a titration.

1. Prepare the titration system. Refer to the titration system and the application working procedure documentation.
2. Make sure that the titration module and the probe are correctly installed and connected.
3. Install the probe on the electrode holder.
4. Connect the titration module to the legacy adapter of the titration system. Follow the instructions that show on the TitraLab AT1000 to set up the legacy adapter.
5. Select the applicable wavelength on the module 3-position switch: 525 nm or 610 nm based on your application working procedure.

Usually, select the wavelength that is non-absorbing for the coloured solution at the titration start then absorbing for the color after the end point. Refer to [Table 1](#) and [Table 2](#).

6. Use the titration system and the instructions on the application working procedure to complete the titration.

**Table 1 Wavelength selection based on indicator**

Indicator	Non complexed form	Complexed form	Wavelength
Calmagite	Red	Blue	610 nm
Eriochrome black T	Red	Blue	610 nm

**Table 2 Wavelength selection based on pH indicator**


pH indicator	Colour = f (pH)	Wavelength
Helianthin	red = <3.1 orange = 3.1 to 4.4 yellow = > 4.4	525 nm
Phenolphthalein	no colour = <8.3 pink = 8.3 to 10 purple = >10.0	610 nm

## Example of a titration curve

In the example, complexometric titration with EDTA at pH = 10 is used to determine the total hardness of a mineral water. The equivalent point is found with colorimetric detection. The change in color (wine red to blue) causes a potential jump of 600 to 800 mV and an accurate and reproducible definition of the equivalent point (standard deviation around 1% on the 3 tests). Refer to [Figure 3](#).

**Figure 3 Titration curve**

1 Potential in mV	2 Volume in mL
-------------------	----------------

## Maintenance

### Clean the instrument

Clean the exterior of the instrument with a moist cloth, and then wipe the instrument dry.

### Clean the probe

The OPT3000 probe has a low maintenance level. Clean the probe with deionised water and dry with a soft tissue. Do not use organic solvents to clean the probe. Do not put the probe in organic solvent.

Replace the probe tip only if the mirror is damaged (e.g. scratches). Refer to the probe documentation for additional information.

## Troubleshooting

Problem	Possible cause	Solution
No output signal	The probe is defective. The connections are not correct.	<ol style="list-style-type: none"> <li>1. Remove the probe from the solution and make sure that the probe transmits light (red for the 610 nm position and green for the 525 nm position).</li> <li>2. Examine the titration system connections and power connections. Make sure that the 3-position switch is set to 525 nm or 610 nm. On the titration system, examine the probe connection. Disconnect and connect the probe again. Do a direct measurement to examine the output signal.</li> </ol>
There is noise in the output signal or the output signal is not stable	The probe is not correctly installed on the electrode holder. There is dirt on the probe mirror tip.	<ol style="list-style-type: none"> <li>1. Make sure that the probe is firmly attached in the sensor holder. A small vibration can cause not stable output signals.</li> <li>2. Make sure that the probe tip is correctly installed in the probe body.</li> <li>3. Make sure that there are no unwanted materials in the optical path of the probe.</li> </ol>
The output signal decreases and becomes low	There is dirt on the probe tip. The optical fibre is damaged.	<ol style="list-style-type: none"> <li>1. Make sure that there are no unwanted materials in the optical path of the probe. Clean the probe.</li> <li>2. Interchange the two LED Out fibre connections (525 nm and 610 nm).</li> <li>3. Examine the mirror surface and the optical fibre. Clean or polish the mirror surface and the optical fibre. Clean the probe.</li> <li>4. Replace the probe tip.</li> </ol>

## Replacement parts and accessories

**Note:** Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Description	Item no.
Photocolorimetric Titration Module PTM450	X61T005
Photocolorimetric probe OPT300	E71T001
Mirror tip for OPT300 probe	E91T001
Power supply (without power cable)	LZE144
Power cable (EU)	YAA080
Power cable (US)	XLH055
Power cable (UK)	XLH057

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)