
*M241Pt*  
*Platinum Electrode*

Operating Instructions

Mode d'Emploi

**M241Pt**

*The design principles of the M241Pt Electrode*  
*Schéma de principe de l'électrode M241Pt*

# Platinum Electrode

## Operating Instructions

---

The M241Pt Platinum Electrode is ideal for a wide range of applications. It can be used together with a standard Ag/AgCl reference electrode or, depending on the application, with a reference electrode which has an alternative liquid junction or reference system.

### Preparing for measurement

1. Remove the protection cap from the electrode.
2. Before starting a measurement, rinse the electrode with distilled water.
3. Check that the electrode body and sensing element are not damaged.

## **Cleaning and maintenance**

### **Platinum plate**

The sensing part of the electrode can be polished with finishing abrasive bands, part no. B90X550. While polishing, take great care not to bend the platinum plate.

In order to have reproducible and accurate results with the electrode, a few physical and electrochemical treatments are required. These treatments will depend on the last measurement performed:

### **Mild Cleaning**

- Place the electrode overnight in RENOVO•N Normal Cleaning Solution.

### **Degreasing**

- Remove mineral oils by rinsing the electrode with trichlorethylene and dabbing dry with a soft tissue.
- Remove organic oils by cathodic electrolysis in a hot alkaline solution, e.g. 15 g/l NaOH + 50 g/l Na<sub>2</sub>CO<sub>3</sub> at 50°C.

### **Degassing**

- Hydrogen produced during the electrolysis (see above) must be removed. Dip the electrode in water at 60°C and flush with pure nitrogen.

### **Chemical cleaning**

- The electrode can be cleaned with a sulphochromic solution, or with a boiling solution of hydrochloric acid.

### **Storage**

- Replace the protection cap, and store the electrode in its box at room temperature.

## Accessories

<b>Description</b>	<b>Part No.</b>
BA03 Finishing Abrasive Bands	B90X550
○RENOVO•N Normal Cleaning Solution, 250 ml	S16M001

## Specifications

Temperature range:	0 to 100°C
Sensing element:	M241Pt : platinum plate 5 x 5 mm
Connection:	Banana plug

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