
MC6091Ag-9
Combined Silver Electrode
Electrode d'Argent Combinée

MC6091Ag-9
Combined Silver and
Hg/Hg₂SO₄ Electrode

Introduction

The MC6091Ag-9 is a combined silver and Hg/Hg₂SO₄ Electrode. It can be used for the determination of halides either by direct measurements or by potentiometric titration.

Preparation for Measurement

1. Remove the protection cap from the electrode and any seals covering the filling hole.
2. Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.
3. Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.
If necessary, refill with KS160 Saturated K₂SO₄ Solution.
4. In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the sensing element down.

Maintenance

Silver rod

The sensing part of the electrode can be polished with finishing abrasive bands, part no. B90X550 or with silver polish.

Cleaning

The electrode must be carefully cleaned before use, but in many cases a thorough rinsing with distilled water will suffice. Owing to the high chemical resistance of the electrode, most cleaning chemicals and solvents may be used (the electrode tip and stem consist of pure silver and glass).

- To remove most electrode deposits which cannot be removed with pure water: flush with dilute nitric acid and/or ammonium hydroxide followed by distilled water.
- In case of severe contamination, clean the electrode with silver polish and then flush with acetone and distilled water.

Hg/Hg₂SO₄ Reference Electrode

1. The electrode should be rinsed with distilled water after measurements.
2. Check frequently the level of filling solution.
3. In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc...
 - KS400 Pepsin in HCl Solution or RENOVO•X Xtra Strong Cleaning Solution: protein deposits (milk, cheese, serums...). Duration of treatment 1 to 2 hrs.
 - KS410 Thiourea Solution: for porous pins contaminated with sulphides or blocked by an AgCl precipitate. Duration of treatment, a few hours until the porous plug turns white.
 - RENOVO•N Normal Cleaning Solution : greasy or oily deposits...

Storage

Between measurements: leave the electrode in KS160 Saturated K₂SO₄ Solution.

Long-term: seal the filling hole with paraffin film or with the electrode stopper. Cover the silver rod and porous pin with the protection cap filled with KS160 Saturated K₂SO₄ Solution.

Accessories

Finishing Abrasive Bands	B90X550
KS160 Saturated K ₂ SO ₄ Solution, 500 ml	C20C500
○ RENOVO•N Normal Cleaning Solution, 250 ml	S16M001
□ RENOVO•X Xtra Strong Cleaning Solution, 250 ml	S16M002
KS400 Pepsin in HCl Solution, 250 ml	C20C370
CL114 Electrode Cable, (1 m) with BNC plug	A94L114
CL116 Electrode Cable, (1 m) with type 7 plug	A94L116

Specifications

Temperature range:	0 to 60°C
Sensing element:	Silver rod
Reference element:	Hg/Hg ₂ SO ₄

In case of contact with mercury, wash exposed areas immediately with plenty of water.

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RADIOMETER ANALYTICAL SAS
D31M089 • Radiometer Analytical SAS • France • 2005-12C