
DOC022.97.80305

Lachat Micro Dist[©]

05/2014, Edition 3

User Manual
Manuel d'utilisation
Manual del usuario

English	3
Français	28
Español.....	54

Table of Contents

Specifications on page 3	Maintenance on page 20
General information on page 3	Troubleshooting on page 22
Installation on page 7	Replacement parts and accessories on page 23
User interface and navigation on page 8	Appendix on page 24
Operation on page 9	

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	72 x 17 x 15 cm (28.3 x 6.7 x 5.9 in.)
Enclosure	Aluminum
Weight	20 kg (45 lb)
Safety	Laboratory 61010-1, 61010-2-010 (UL, CSA, EN safety standards)
Pollution degree/installation category	2; II
Protection class	I
Power requirements	100–115 VAC, 50/60 Hz, 9–11 A 230 VAC, 50/60 Hz, 5 A
Fuse	100–115 VAC: 15 A, 250 VAC, slow-blow 230 VAC: 6.3 A, 250 VAC, slow-blow
Operating environment	Temperature: 5 to 40 °C (41 to 104 °F) Humidity: 92.5%, non-condensing maximum
Sample size	6.0 mL capacity per Micro Dist tube
Sample capacity	21 Micro Dist tubes
Tube material	Analytical-grade pure polypropylene
Temperature range	30 to 135 °C (86 to 275 °F)
Temperature accuracy	±0.5 °C (0.9 °F)
Temperature variation (hole-to-hole)	Less than ±0.5 °C (0.9 °F)
Controller set point error	Less than 1% of the temperature set point
Certifications	100–115 VAC model: UL and CSA safety standards (cETLus mark) 230 VAC model: CE

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.
	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	This symbol, when noted on the product, identifies the location of a fuse or current limiting device.
	This symbol indicates that the marked item can be hot and should not be touched without care.

Certification

Canadian Radio Interference-Causing Equipment Regulation, ICES-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.


Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Product overview


This instrument is used with Micro Dist® tubes to remove volatile analytes from liquid or solid samples with distillation. Refer to [Figure 1](#) and [Figure 2](#). This instrument is for indoor use only.

Figure 1 Instrument overview

1 Heater block	3 Power switch	5 Drain plate (under the heater block)
2 Display and keypad	4 Power connector	

Figure 2 Micro Dist tube overview


1 Tube cap	4 6.0 mL mark	7 Stop ring
2 Cap membrane	5 Break point	8 Sample tube
3 Collector tube	6 Collector tube membrane	

Product components

Make sure that all components have been received. Refer to [Figure 3](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 3 Instrument components

1 Collector tube rack	8 Pipette
2 Micro Dist ¹	9 Timer
3 Assembly press	10 Heat-resistant gloves, size large
4 Instructional CD	11 Tube cap (10x)
5 Sample tube rack	12 User-fill collector tube ² (10x)
6 Heat-resistant gloves, size medium	13 Cap membranes (20x)
7 Power cord	14 Sample tube (10x)

¹ A17102 includes only the Micro Dist block and a power cord.

² Not for use with the Ammonia-1 or -2 Method. Refer to [Accessories](#) on page 24.

Installation

Assemble the tube racks

Refer to the assembly instructions supplied with the tube racks.

Connect the power cord

DANGER

Electrocution hazard. Protective Earth Ground (PE) connection is required.

▲ DANGER

Electrocution and fire hazard. Attach the power cord so that it cannot touch hot surfaces.

▲ DANGER

Electrical shock and fire hazards. Replacement power cords must:

- Have the correct plug style for the outlet connection
- Have a rating sufficient for the supply voltage and current. Refer to the requirements in the Specifications section.
- Meet or exceed local electrical code requirements

▲ WARNING

Personal injury hazard. Instruments or components are heavy. Use assistance to install or move.

NOTICE

Make sure that there is at least 15 cm (6 in.) of clearance around the instrument to prevent heat damage to objects near the instrument and so there is sufficient air flow around the instrument.


This instrument is rated for an altitude of 2000 m (6562 ft) maximum. Use of this instrument at an altitude higher than 2000 m can slightly increase the potential for the insulation to breakdown, which can result in an electric shock hazard. The manufacturer recommends that users with concerns contact technical support.

1. Put the instrument on a stable, level, heat-resistant surface. Do not put the instrument in direct sunlight.
2. Connect the power cord to the power connector. Refer to [Figure 1](#) on page 6.
3. Connect the power cord to an electrical outlet with protective earth ground.

User interface and navigation

[Figure 4](#) shows the instrument display and keypad. [Table 1](#) gives the key functions.

Figure 4 Display and keypad

1 Temperature of the heater block (°C)

2 Heating light

Table 1 Key functions

Key	Function
	Shows or changes the temperature setting.
	Increases the value.
	Decreases the value.
ENTER	Saves the change.

Operation

⚠ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

⚠ WARNING

Gas inhalation hazard. Operate the instrument in a fume hood to prevent exposure to hazardous gas.

⚠ WARNING

Explosion hazard. To prevent the increase of pressure in the tubes during distillation, do not cover the ventilation opening in the tube caps.

⚠ WARNING

Explosion hazard. Hot sample can exit the sample tube explosively if foam or debris touches the membrane during distillation.

⚠ CAUTION

Fire hazard. This product is not designed for use with flammable liquids.

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Refer to [Troubleshooting](#) on page 22 if foam or debris touches the membrane during distillation.

Hazardous gas information

The Micro Dist methods do not typically release hazardous gases, but a hazardous situation can occur in a fault condition. Refer to [Table 2](#) for the gases that can be released from the Micro Dist.

Table 2 Hazardous gas information for specific methods

Method	Gas	Permissible exposure limit concentration (per mg/m ³)
Cyanide-1	Cyanide (HCN)	5 mg
Cyanide-2	Cyanide (HCN)	5 mg
Cyanide-3 (WAD)	Cyanide (HCN)	5 mg
Cyanide-5	Cyanide (HCN)	5 mg
Phenolics-1	Phenol (C ₆ H ₅ OH)	19 mg
Sulfide-1	Hydrogen sulfide (H ₂ S)	15 mg
Sulfide-2	Hydrogen sulfide (H ₂ S)	15 mg
Ammonia-1	Ammonia (NH ₃)	35 mg
Ammonia-2	Ammonia (NH ₃)	35 mg
Tritium-1	Tritium (³ H)	—


Preparation notes

- The instrument takes approximately 40 minutes to increase to the set temperature.
- For very low-concentration methods, such as ammonia, flush the collector tubes with deionized water before use.
- Use only ammonia collector tubes with the ammonia Micro Dist methods. Typically, the ammonia collector tubes contribute an ammonia blank of 0.015 mg or 15 ppb N/L in the 6.0 mL of distillate. Make sure to distill both calibration standards and all samples with the ammonia tubes.
- Use user-fill tubes the same day they are filled with trapping solution.
- Keep pre-fill tubes sealed in the package until ready to use. Use pre-fill tubes within 2 weeks after the sealed package is opened. Record the date the package was opened on the package.
- Do not put collector tubes in a refrigerator as cold causes the collector tube membrane to decrease in size. Trapping solution can then get around the collector tube membrane.
- For the cyanide methods, always distill the standards with the samples.

- The manufacturer recommends that the user distills a set of standards at least once a day or once every batch of samples. Then, use the Micro Dist standards as the calibration standards for the method used to determine the samples.

General procedure

Set the temperature

1. Turn the power switch on.
2. Push the loop key.
The display flashes between "SP" and the temperature setting.


3. Use the UP and DOWN arrows to select the temperature. Refer to [Method information](#) on page 16 for the temperature setting.

4. Push Enter.

5. Push the loop key.
The display shows the temperature of the heater block.

Prepare the collector tubes

If pre-filled collector tubes are used, do only step 1.


1. Put one collector tube in the rack with the M end up for each sample to distill.

2. Prepare the trapping solution if applicable. Refer to [Trapping solutions](#) on page 17. Put trapping solution in each collector tube.


3. Put a cap membrane over a collector tube. Put a tube cap on the collector tube. Refer to [Figure 2](#) on page 6.

4. Do step 3 for each collector tube.


Assemble the tubes

1. Put the necessary number of sample tubes in the sample tube rack (up to 21).


2. Put sample in each sample tube. Refer to [Method information](#) on page 16 for the quantity.


3. Add the releasing agent, neutralizing agent and/or buffer solution to one sample tube if applicable.


Refer to [Releasing agents](#) on page 18, [Neutralizing agents](#) on page 19 and [Buffer solutions](#) on page 19.

4. Assemble and compress the tube. Refer to [Compress the tube](#) on page 12.

5. Put the tube in the collector tube rack with the sample tube down.


6. Do steps 3–5 for each sample tube.


Compress the tube

Refer to [Figure 5](#) to compress the tube. Hold the middle of the collector tube securely with one hand while the tube is compressed.

Figure 5 Compress the tube

Distill the samples


1. When the instrument temperature is at the temperature setting, put on heat-resistant gloves and safety glasses.


2. Push each tube assembly into a hole in the heater block. Push until the stop ring on the sample tube hits the heater block. Refer to [Figure 2](#) on page 6.

3. Take off the gloves. Set and start the timer. Refer to [Method information](#) on page 16 for the distillation time.

4. When the time has expired, put on the heat-resistant gloves.


5. Remove one tube assembly from the heater block.


6. Within 4 seconds, pull down and twist the sample tube to remove it from the collector tube.


7. Discard the hot sample and the sample tube.

8. Invert the collector tube so that the D end is up.


9. Put the collector tube in the collector tube rack.


10. Do steps 5-9 until all the tubes are removed from the heater block. Then, take off the gloves.


Collect the distilled samples


1. Let the tubes cool for at least 10 minutes.


2. Remove a tube from the rack.


3. Hold the tube horizontally. Slowly roll the collector tube to collect drops on the sides of the tube.


4. Hold the tube vertically with the D end up. Tap the sides of the collector tube as necessary to move drops down to the M end.

5. Break off the D end of the tube with hands. Discard the D end of the tube.


6. Put a cap on the tube.


7. Put Parafilm® over the ends of the tube if the sample will not be used immediately.

Dilute the distilled samples


For Tritium-1, add distilled sample and dilution water to the scintillation vials as necessary. Add scintillation cocktail and count. Refer to the instructions from the manufacturer.
For the other methods:


1. Remove the cap from a tube.


2. Add deionized water to the 6.0 mL mark on the tube.

3. Put the cap on the tube.

4. Gently shake the tube side to side to mix. Do not invert the tube.

5. Put the tube in the rack with the M end down.

6. Do steps 1–5 for each tube.

Measure the sample

Measure the sample with the QuikChem Method or equivalent. Refer to [Methods of determination](#) on page 20.

There is enough distillate in each tube for two 25-second sample periods.

Method information

Method	Sample type	Range	Quantity	Temperature	Time
Cyanide-1	Free and Combined Cyanide/L	0.002–100 mg	Liquid sample or standard: 6.0 mL or 6 g Soil or sludge ¹ : 0.5 to 1.0 g mixed with 5 mL of deionized water <i>Note: The amount of cyanide in the soil or sludge in 6 mL should be between about 0.012 and 600 µg.</i>	120 °C	30 min
Cyanide-2	Free and Combined Cyanide/L in caustic extract (1.25 M NaOH)	0.5–50 mg			
Cyanide-3 (WAD)	Free and Weak Acid Dissociable (WAD) Cyanide/L	0.002–0.50 mg			
Cyanide-5	Free and Combined Cyanide/L in the presence of sulfide	0.05–50 mg			
Phenolics-1	Phenolics/L	0.002–50.0 mg	Liquid sample or standard: 6.0 mL or 6 g Soil or sludge ¹ : 0.5 to 1.0 g mixed with 5 mL of deionized water Adjust to approximately pH 4 with 1 M NaOH or 10% H ₂ SO ₄	130 °C	90 min
Sulfide-1	Acid-soluble sulfide	0.2–50 mg/kg	Standard: 5.5 mL or 5.5 g Liquid sample: 1.0 g mixed with 5 mL of deionized water Soil or sludge ¹ : 0.5 to 1.0 g mixed with 5 mL of deionized water	120 °C	30 min
Sulfide-2	Acid-soluble sulfide/L	0.01–10.0 mg	Liquid sample or standard: 6.0 mL or 6 g Soil or sludge ¹ : 0.5 to 1.0 g mixed with 5 mL of deionized water		

Method	Sample type	Range	Quantity	Temperature	Time
Ammonia-1: Phenate Method/ISE	N/L	0.010–2.0 mg	Liquid sample or standard: 6.0 mL or 6 g Soil or sludge ¹ : 0.5 to 1.0 g mixed with 5 mL of deionized water	120 °C	30 min
Ammonia-2: Nesslerization Method	N/L	0.010–2.0 mg	Liquid sample or standard: 6.0 mL or 6 g Soil or sludge ¹ : 0.5 to 1.0 g mixed with 5 mL of deionized water		
Tritium-1	Volatile tritiated compounds <i>Note: If the sample is primarily oil, do not distill, but count directly with the scintillation counter.</i>	—	Liquid sample or standard: 6.0 mL or 6 g	120 °C	90 min

¹ Refer to [Distill a solid sample](#) on page 23 to identify the best weight of sample to add to each tube.

Trapping solutions

Methods not shown do not use a trapping solution.

Method	Trapping solution ¹	Preparation
Cyanide-1	1.5 mL or 1.55 g of 1.00 M standardized NaOH solution	—
Cyanide-2		—
Cyanide-3 (WAD)		
Cyanide-5	1.5 mL or 1.55 g of trapping solution	Dissolve 0.80 g of PbCO ₃ in 1 L of 1.00 M standardized NaOH solution.
Sulfide-1	2.0 mL or 2.0 g of 0.043 M zinc acetate solution	Combine: <ul style="list-style-type: none">• 8.78 g zinc acetate dihydrate• 0.10 g concentrated HCl• 880.0 g deionized water• 43.2 g of a 37% formaldehyde solution Mix well until dissolved.
Sulfide-2	1.5 mL or 1.55 g of 1.00 M standardized NaOH solution	—
Ammonia-1: Phenate Method/ISE	1.0 mL or 1.0 g of 0.016 M standardized sulfuric acid solution	—
Ammonia-2: Nesslerization Method	2.0 mL or 2.0 g of 0.13 M boric acid solution	Dissolve 8.0 g of boric acid in 995.0 g of deionized water.

¹ Other concentrations of standardized solution can be used. Adjust the amount added to each tube accordingly.

Releasing agents

Methods not shown do not use a releasing agent.

Method	Releasing agent	Preparation
Cyanide-1	0.75 mL of 7.11 M sulfuric acid / 0.79 M magnesium chloride solution	Refer to Prepare the releasing agent—Cyanide-1, -2 or -5 on page 18.
Cyanide-2		
Cyanide-5		
Cyanide-3 (WAD)	0.75 mL of 0.50 M zinc acetate/0.52 M sodium acetate/0.87 M acetic acid	Refer to Prepare the releasing agent—Cyanide-3 (WAD) on page 18.
Sulfide-1	0.45 mL of 9.0 M sulfuric acid solution <i>Note:</i> When an unknown matrix is initially distilled, identify the amount of concentrated sulfuric acid necessary to decrease the sample pH to 1 or less. The volume necessary may be greater than 0.45 mL.	Refer to Prepare the releasing agent—Sulfide-1 or -2 on page 18.
Sulfide-2		

Prepare the releasing agent—Cyanide-1, -2 or -5

Makes 200 mL of releasing agent.

1. In a fume hood, put a 500 mL beaker on a top-loading balance.
2. Set the balance to zero.
3. Add 110.8 g of deionized water.
4. Add and fully dissolve 32.2 g of magnesium chloride hexahydrate ($\text{MgCl}_2 \cdot 6 \text{ H}_2\text{O}$).
5. Slowly add 139 g of concentrated sulfuric acid in increments of 40 g. Swirl and let cool between additions.
6. Let the solution cool to ambient temperature.
7. Put the solution in a closed container.

Prepare the releasing agent—Cyanide-3 (WAD)

Makes 1 L of releasing agent.

1. In a fume hood, put a 1 L container on a top-loading balance.
2. Set the balance to zero.
3. Add 950 g of deionized water.
4. Add and fully dissolve 70.97 g of sodium acetate trihydrate ($\text{NaC}_2\text{H}_3\text{O}_2 \cdot 3 \text{ H}_2\text{O}$).
5. Add and fully dissolve 100 g of zinc acetate monohydrate ($\text{ZnC}_2\text{H}_3\text{O}_2 \cdot \text{H}_2\text{O}$).
6. Add 52.5 g of glacial acetic acid ($\text{C}_2\text{H}_3\text{O}_2$).
7. Put the solution in a closed container.

Prepare the releasing agent—Sulfide-1 or -2

Makes 300 mL of releasing agent.

1. In a fume hood, put a 500 mL beaker on a top-loading balance.
2. Set the balance to zero.
3. Add 150 g of deionized water.
4. Slowly add 276 g of concentrated sulfuric acid in increments of 40 g. Swirl and let cool between additions.
5. Let the solution cool to ambient temperature.
6. Put the solution in a closed container.

Neutralizing agents

Methods not shown do not use a neutralizing agent.

Method	Neutralizing agent
Cyanide-2	200 µL of concentrated H ₂ SO ₄

Buffer solutions

Methods not shown do not use a buffer solution.

Method	Buffer solution	Preparation
Ammonia-1: Phenate Method/ISE	Samples preserved with 2 mL concentrated sulfuric acid per liter: 1.0 mL of 0.025 M borate buffer (Na ₂ B ₄ O ₇). Buffer at pH 9.5 to decrease the hydrolysis of cyanates and organic nitrogen compounds.	Refer to Prepare the buffer solution—Ammonia-1 or -2 on page 19.
	Samples not preserved with acid: 0.75 mL of 0.011 M borate buffer (Na ₂ B ₄ O ₇) Buffer at pH 9.5 to decrease the hydrolysis of cyanates and organic nitrogen compounds.	
Ammonia-2: Nesslerization Method	0.75 mL of 0.011 M borate buffer (Na ₂ B ₄ O ₇) Buffer at pH 9.5 to decrease the hydrolysis of cyanates and organic nitrogen compounds.	
Tritium-1	1 or 2 drops of 10 M NaOH	—

Prepare the buffer solution—Ammonia-1 or -2

1. Put a 1 L container on a top-loading balance.
2. Set the balance to zero.
3. Add 964 g of deionized water.
4. Add and fully dissolve sodium borate (Na₂B₄O₇). Refer to [Table 3](#) for quantity.
5. Add sodium hydroxide (NaOH). Refer to [Table 3](#) for quantity.
6. Mix well.
7. Put the solution in a closed container.

Table 3 Borate buffer information

Concentration	Sodium borate	Sodium hydroxide
0.025 M borate buffer	5 g (anhydrous)	22 g (solid)
	9.5 g (decahydrate, 10 H ₂ O)	22 g (solid)
0.011 M borate buffer	2.17 g (anhydrous)	36 mL of 0.1 M solution
	4.2 g (decahydrate, 10 H ₂ O)	36 mL of 0.1 M solution

Methods of determination

Method	Matrix and chemistry	Solution	QuikChem Method
Cyanide-1	Waters, solids, strong acid dissociable (SAD), total cyanide	Original sample volume, but now in 0.25 M NaOH	10-204-00-1-X or 10-204-00-1-X2
Cyanide-2	Caustic extracts, SAD, total cyanide		10-204-00-1-W
Cyanide-3 (WAD)	Waters, solids, weak acid dissociable (WAD)		10-204-00-1-WX
Cyanide-5	Waters, solids, SAD (total cyanide) in the presence of sulfide	Original sample volume, but now in 0.25 M NaOH and 0.8 mM Pb	10-204-00-1-X
Phenolics-1	Waters, solids, 4-AATP	—	10-210-00-1-X
Sulfide-1	Waters, iodometric determination	Original sample volume, but now in sulfide (ZnS)	No QuikChem Method Use an EPA Method: <ul style="list-style-type: none">• SW-846 Method 9030, Revision 1 (Dec. 1987)• Method 376.1
Sulfide-2	Waters/MTB colorimetric determination	Original sample volume, but now in 0.25 M NaOH	10-116-29-1-X
Ammonia-1: Phenate Method/ISE	Waters, phenate colorimetric or ISE determination	Original sample volume, but now in 0.003 M sulfuric acid	10-107-06-1-X
Ammonia-2: Nesslerization Method	Waters, solids, Nesslerization	Original sample volume, but now in 0.043 M boric acid solution	No QuikChem Method. Use standard methods.
Tritium-1	Contaminated waters, solids	—	—

Maintenance

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

NOTICE

Do not disassemble the instrument for maintenance. If the internal components must be cleaned or repaired, contact the manufacturer.

Clean the instrument

⚠ CAUTION

Burn hazard. Turn off power and allow the instrument to cool before this procedure.

⚠ CAUTION

Fire hazard. Do not use flammable agents to clean the instrument.

⚠ CAUTION

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

NOTICE

To prevent instrument damage, do not use cleaning agents such as turpentine, acetone or similar products to clean the instrument or the display.

Keep the instrument clean to get continuous and accurate operation.

1. Turn the power switch off.
2. Disconnect the power cord from the electrical outlet.
3. When the instrument is cool, clean the surface of the instrument with a soft, moist cloth and a weak soap solution.
4. Let the instrument fully dry before use.

Clean a spill

If a large amount of liquid gets in the heater block:

1. Disconnect the power cord from the electrical outlet.
2. Let the instrument cool.
3. Rinse the holes with water. Water drains out of the bottom of the instrument through the drain plate.
4. Let the instrument fully dry before use.

Replace the fuse

⚠ DANGER


Fire hazard. Use the same type and current rating to replace fuses.

Fuse failure is typically caused by an instrument problem. If the fuse continues to fail, contact Customer Service.

Note: 100–115 VAC model: Although there are positions for two fuses, only the upper fuse in the fuse drawer is connected in the circuit. The lower fuse is used as a replacement and is not necessary.

Refer to [Figure 6](#) to replace the fuse.

Figure 6 Replace the fuse

Troubleshooting

Problem	Possible cause	Solution
Low or not consistent analyte recoveries	Tubes are not pushed fully into the heater block	Push tubes into the heater block until the stop ring on the sample tube hits the heater block.
	Too much NaOH was added to the samples	Lower the pH with sulfuric acid before distillation. <i>Note: If too much NaOH is added to samples, the amount of HCN in the samples is decreased.</i>
	Sample tubes were not removed immediately after distillation	Immediately remove the sample tubes from the tube assembly after distillation.
	Too much liquid was put in the sample tube which caused distilled sample to be pulled back into the sample tube	Use a calibrated pipette to dispense sample, releasing agent, neutralizing agent and/or buffer solution into the sample tube.
	A method expects a higher NaOH concentration so the samples give high recoveries	For method Cyanide-1, the method of determination which follows the distillation must anticipate both 6.0 mL samples and 0.25 M NaOH. For example, when using USEPA Method 335.2 and scaling down the reagent quantities, but the solution is diluted to 100 mL, color development will not occur. See the scaled-down version of this method for an example.
	A method expects a lower NaOH concentration so the samples give low recoveries	Some "auto-analyzer" methods expect a different concentration of NaOH. If the undistilled standards used in these methods are not made up in 0.25 M NaOH, then there will be a matrix effect when the Micro Distilled samples, which are always in 0.25 M NaOH, are determined.
Large blanks in cyanide	Contaminated MgCl ₂	Use only 99.9999% pure or other ultra-pure grade. Use high quality grade (e.g., HPLC) of sulfuric acid when possible.
Sample boils too intensely in the sample tube	—	Add one Hengar granule or other small (ca. 1 mm diameter) boiling stone to the sample tube.

Problem	Possible cause	Solution
Results from the Micro Dist methods do not agree with other methods	Analytical method is not calibrated with the Micro Dist standards	Calibrate the analytical method with Micro Dist standards. If the recoveries get better for reference samples, use Micro Dist standards. <i>Note: Analytical methods used after the distillation must be scaled down in proportion to the smaller Micro Dist sample.</i>
Temperature of the block takes hours to get to the temperature set point	Instrument failure	When it takes more than 30–40 minutes for the instrument to increase to the temperature setting, contact Technical Support.
Collector tube cannot be broken	—	Use pruning shears or poultry scissors to cut the collector tube at the break point.
Foam and/or debris on the membrane	Soil or sludge sample is high in organics (i.e., grease or oils)	Use a smaller sample size or less deionized water. Refer to Distill a solid sample on page 23. <i>Note: Foam wets the hydrophobic membrane which causes the membrane to let distilled sample back into the sample tube.</i>

Distill a solid sample

Use different weights to identify the best weight of sample to add to each tube for each matrix type. If the sample is high in organic content, use only 0.5 g or less of sample. If the sample is low in organic content, use up to 1 g of sample.

Soil and sludge samples are diluted with deionized water. Calculate the amount of sample in mg/kg after analysis:

$$A = C \times ((1 \text{ L} / 1000 \text{ mL}) \times V)$$

$$B = A \text{ divided by } (W \times 1000 \text{ g/kg})$$

Where:

A = mg of analyte in the sample

B = mg of analyte per kg of sample (mg/kg)

C = determined concentration (mg/L)

V = dilution volume in tube (mL)

W = weight of the original sample (g)

If foam or debris continues to get on the collector tube membrane after the sample is decreased to 0.5 g:

1. Add 4–5 mL of deionized water instead of 6 mL to the sample tube.
2. Add activated charcoal to the sample tube so that it covers the surface of the solid.
3. Fill the empty space in the sample tube with glass wool.

Note: For soil or organic samples that contain cyanide, Biobeads™, manufactured by BIO-RAD, part number SM-2, are effective in laboratories.

Do a test on a known standard and a spiked sample that puts foam or debris on the membrane to identify if these procedures correct the problem.

Replacement parts and accessories

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Quantity	Item no.
Assembly press	1	17023L
Fuse, 15 A, 250 VAC, slow-blow	1	17085
Fuse, 6.3 A, 250 VAC, slow-blow	1	20210
Micro Dist, digital, 110 VAC	1	A17102
Pipette, 0.2–1.0 mL	1	BBP078
Pipette tips for BBP078	100	BBP079
Power cord, 115 VAC, 15 A	1	17040L
Power cord, 230 VAC, 10 A	1	31196
Tube, Micro Dist user-fill	10	A17017
Tube rack, 24 holes	1	2497903
Tube rack, 60 holes	1	2497901

Accessories

Description	Quantity	Item no.
Tube, Micro Dist user-fill	50	A17517
Tube, Micro Dist user-fill	100	A17117
Tube, Micro Dist user-fill for Ammonia-1 or -2 ¹	50	A17517A
Tube, Micro Dist user-fill for Ammonia-1 or -2 ¹	100	A17117A
Tube, Micro Dist pre-filled, for Cyanide-1, -2, -3 and Sulfide-2	21	A17001
Tube, Micro Dist pre-filled, for Phenolics-1	21	A17002
Tube, Micro Dist pre-filled, for Tritium-1	10	A17014

¹ Specially prepared for determination of NH₃.

Appendix

Cyanide-1 Method

USEPA Method 335.2 scaled-down

For 6.0 mL of micro distilled sample, the scale-down factor for this batch colorimetric method is 6/50 or 0.120. The changes suggested in apparatus or technique are [added in brackets]. The scale-down gives the same concentrations of reagents, sample and analyte from Section 8.7, p 335.2–4.

8.7 Withdraw 6.0 mL of the solution from the flask [collector tube] and transfer to a 20 mL erlenmeyer flask. Add 1.80 mL of sodium phosphate solution (7.6) and mix.

8.7.1 Pyridine-barbituric acid method: Add 0.24 mL of chloramine T (7.12) and mix. After one to two minutes, add 0.60 mL of pyridine-barbituric acid solution (7.13.1) and mix. [Add 3.36 mL of deionized water] and mix again. Allow eight minutes for color development, then read absorbance at 578 nm in a 1 cm cell within 15 minutes.

8.7.2 Pyridine-pyrazolone method: Add 0.060 mL of chloramine T (7.12) and mix. After one to two minutes add 0.60 mL of pyridine-pyrazolone solution (7.13.1) and mix. [Add 3.54 mL deionized water] and mix again. After 40 minutes read absorbance at 620 nm in a 1 cm cell.

Note: Some distillates can contain compounds that have a chlorine demand. One minute after the addition of chloramine T, measure for residual chlorine with KI-starch paper. If the test is negative, add an additional 0.060 mL of chloramine T. After one minute, measure the sample again.

Note: More than 0.060 mL of chloramine T will prevent the color from developing with pyridine-pyrazolone.

Equivalency of releasing agents

Method	Concentrations of releasing agents
Micro Dist Method Cyanide-1	0.75 mL of 7.11 M H ₂ SO ₄ is added to 6.0 mL of sample. This gives a sample concentration in H ₂ SO ₄ : $7.11 \text{ M H}_2\text{SO}_4 \times 0.75 \text{ mL} / (6.0 \text{ mL} + 0.75 \text{ mL})$ $= 0.79 \text{ M H}_2\text{SO}_4$ This same solution is also 0.79 M in MgCl ₂ · 6 H ₂ O. This gives a concentration of: $0.79 \text{ M MgCl}_2 \cdot 6 \text{ H}_2\text{O} \times 0.75 \text{ mL} / (6.75 \text{ mL total})$ $= (0.088 \text{ M MgCl}_2 \cdot 6 \text{ H}_2\text{O})$
USEPA Method 335.2	50 mL of 9 M (18 N) H ₂ SO ₄ is added to 500 mL of sample. This gives a sample concentration in H ₂ SO ₄ : $0.9 \text{ M H}_2\text{SO}_4 \times 50 \text{ mL} / (500 \text{ mL} + 50 \text{ mL} + 20 \text{ mL})$ $= 0.79 \text{ M H}_2\text{SO}_4$ 20 mL of a 510 g MgCl ₂ · 6 H ₂ O/L (2.5 M) solution is also added to 500 mL of sample. This gives a concentration of: $2.5 \text{ M MgCl}_2 \cdot 6 \text{ H}_2\text{O} \times 20 \text{ mL} / (570 \text{ mL total})$ $= (0.088 \text{ M MgCl}_2 \cdot 6 \text{ H}_2\text{O})$
USEPA Method 335.4	5 mL of 9 M (18 N) H ₂ SO ₄ is added to 50 mL of sample. This gives a sample concentration in H ₂ SO ₄ : $9 \text{ M H}_2\text{SO}_4 \times 5 \text{ mL} / (50 \text{ mL} + 5 \text{ mL} + 2 \text{ mL})$ $= 0.79 \text{ M H}_2\text{SO}_4$ 2 mL of a 510 g MgCl ₂ · 6 H ₂ O/L (2.5 M) solution is also added to 50 mL of sample. This gives a concentration of: $2.5 \text{ M MgCl}_2 \cdot 6 \text{ H}_2\text{O} \times 2 \text{ mL} / (57 \text{ mL total})$ $= (0.088 \text{ M MgCl}_2 \cdot 6 \text{ H}_2\text{O})$

Equivalency of trapping solutions

Method	Concentrations of trapping solutions
Micro Dist Method Cyanide-1	In the Micro Dist Method, 1.5 mL of 1.0 M NaOH becomes 6.0 mL solution after distillation and dilution to the mark: $1.0 \text{ M NaOH} \times (1.5 \text{ mL} / 6.0 \text{ mL}) = 0.25 \text{ M NaOH}$
USEPA Method 335.2	In the macro distillation, 50 mL of 1.25 M NaOH becomes 250 mL solution after distillation and dilution to the mark: $1.25 \text{ M NaOH} \times (50 \text{ mL} / 250 \text{ mL}) = 0.25 \text{ M NaOH}$
USEPA Method 335.4	In the macro distillation, 50 mL of 0.25 M NaOH stays undiluted.

Sulfide-2 Method

Comparison of the Micro Dist Sulfide-2 Method and the SW-846 Method 9030B, Revision 2

Both methods use a trapping solution which is 0.043 M in zinc acetate and contains 1.60% formaldehyde. The stock 37% formaldehyde solution contains 10–15% methanol. 2.0 mL of this trapping solution is in each Sulfide-1 collector tube.

Sources

Method	Sources
Cyanide-1	<ul style="list-style-type: none"> U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 335.2. EPA-600, 4-79-020, Revised March 1983. APHA-AWWA-WEF Standard Methods for the Examination of Water and Wastewater, 20th ed. (1998), Method 4500-CN-N., pp 4-0 to 4-34. U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Waste, Physical / Chemical Methods (SW-846), Method 9010B, Revision 2, December 1996, "Total and Amenable Cyanide: Distillation". U.S. Environmental Protection Agency, Methods for the Determination of Inorganic Substances in Environmental Samples, Method 335.4. EPA/600/R-93, August 1993. U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Waste, Physical / Chemical Methods (SW-846), Method 9213, Revision 0, December 1996, "Potentiometric determination of cyanide in aqueous samples and distillates with ion-selective electrodes"
Cyanide-2	—
Cyanide-3 (WAD)	APHA-AWWA-WPCF Standard Methods for the Examination of Water and Wastewater, 20th ed., Method 4500-CN-N., pp 4-0 to 4-34.
Cyanide-5	—
Phenolics-1	<ul style="list-style-type: none"> U.S. Environmental Protection Agency, Methods for the Determination of Inorganic Substances in Environmental Samples, Method 420.4. EPA/600/R-93/100, Revised August 1993. U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Waste, Physical / Chemical Methods (SW-846), Method 9065, Revision 0, September 1986, "Phenolics."
Sulfide-1	<ul style="list-style-type: none"> U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 376.1 (Titrimetric, Iodine), EPA-600, 4-79-020, Revised March 1983. U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 376.2 (Colorimetric, Methylene Blue), EPA-600, 4-79-020, Revised March 1983. U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9030B, Revision 2, December 1996, "Acid-soluble and acid-insoluble sulfides: distillation". U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9215, Revision 0, December 1996, "Potentiometric determination of sulfide in aqueous samples and distillates with ion-selective electrode"
Sulfide-2	<ul style="list-style-type: none"> U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 376.1 (Titrimetric, Iodine), EPA-600, 4-79-020, Revised March 1983. U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9030B, Revision 2, December 1996, "Acid-soluble and acid-insoluble sulfides: distillation". U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9215, Revision 0, December 1996, "Potentiometric determination of sulfide in aqueous samples and distillates with ion-selective electrode"
Ammonia-1: Phenate Method/ISE	APHA-AWWA-WEF Standard Methods for the Examination of Water and Wastewater, 20th ed. (1998), Method 4500-NH3 D.to H., pp. 4-106 to 4-112.

Method	Sources
Ammonia-2: Nesslerization Method	—
Tritium-1	<ul style="list-style-type: none"> <li data-bbox="296 198 1035 285">• DOE Method RP580 "Water distillation from soil and aqueous matrices using a microdistillation system for tritium determination", from "DOE Methods for evaluating environmental and waste management samples", DOE/EM-0089T, Rev 2, April 1995, U.S. Dept. of Energy. <li data-bbox="296 285 1023 357">• S. Stieg, "A Miniature Membrane Tube for Rapid Parallel Distillation of Cyanide, Phenolics, Ammonia, Sulfide, Methylmercury, and Tritium from Waters and Solids", American Environmental Laboratory, v9, no. 10, Nov-Dec 1997, pp10-11.

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Fuentes

Método	Fuentes
Cianuro-1	<ul style="list-style-type: none">U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 335.2. EPA-600, 4-79-020, Revised March 1983.APHA-AWWA-WEF Standard Methods for the Examination of Water and Wastewater, 20th ed. (1998), Method 4500-CN-N., pp 4-0 to 4-34.U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Waste, Physical / Chemical Methods (SW-846), Method 9010B, Revision 2, December 1996, "Total and Amenable Cyanide: Distillation".U.S. Environmental Protection Agency, Methods for the Determination of Inorganic Substances in Environmental Samples, Method 335.4. EPA/600/R-93, August 1993.U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Waste, Physical / Chemical Methods (SW-846), Method 9213, Revision 0, December 1996, "Potentiometric determination of cyanide in aqueous samples and distillates with ion-selective electrodes"
Cianuro -2	—
Cianuro -3 (WAD)	APHA-AWWA-WPCF Standard Methods for the Examination of Water and Wastewater, 20th ed., Method 4500-CN-N., pp 4-0 to 4-34.
Cianuro -5	—
Fenólicos-1	<ul style="list-style-type: none">U.S. Environmental Protection Agency, Methods for the Determination of Inorganic Substances in Environmental Samples, Method 420.4. EPA/600/R-93/100, Revised August 1993.U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Waste, Physical / Chemical Methods (SW-846), Method 9065, Revision 0, September 1986, "Phenolics."
Sulfuro -1	<ul style="list-style-type: none">U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 376.1 (Titrimetric, Iodine), EPA-600, 4-79-020, Revised March 1983.U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 376.2 (Colorimetric, Methylene Blue), EPA-600, 4-79-020, Revised March 1983.U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9030B, Revision 2, December 1996, "Acid-soluble and acid-insoluble sulfides: distillation".U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9215, Revision 0, December 1996, "Potentiometric determination of sulfide in aqueous samples and distillates with ion-selective electrode"
Sulfuro -2	<ul style="list-style-type: none">U.S. Environmental Protection Agency, Methods for Chemical Analysis of Waters and Wastes, Method 376.1 (Titrimetric, Iodine), EPA-600, 4-79-020, Revised March 1983.U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9030B, Revision 2, December 1996, "Acid-soluble and acid-insoluble sulfides: distillation".U.S. Environmental Protection Agency, Test Methods for Evaluating Solid Wastes (SW-846) Method 9215, Revision 0, December 1996, "Potentiometric determination of sulfide in aqueous samples and distillates with ion-selective electrode"

Método	Fuentes
Amoníaco--1: Método del fenato/ISE	APHA-AWWA-WEF Standard Methods for the Examination of Water and Wastewater, 20th ed. (1998), Method 4500-NH ₃ D.to H., pp. 4-106 to 4-112.
Amoníaco—2: Método de nesslerización	—
Tritio-1	<ul style="list-style-type: none"> DOE Method RP580 "Water distillation from soil and aqueous matrices using a microdistillation system for tritium determination", from "DOE Methods for evaluating environmental and waste management samples", DOE/EM-0089T, Rev 2, April 1995, U.S. Dept. of Energy. S. Stieg, "A Miniature Membrane Tube for Rapid Parallel Distillation of Cyanide, Phenolics, Ammonia, Sulfide, Methylmercury, and Tritium from Waters and Solids", American Environmental Laboratory, v9, no. 10, Nov-Dec 1997, pp10-11.

**HACH COMPANY World Headquarters**

P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
Tel. (970) 669-3050
(800) 227-4224 (U.S.A. only)
Fax (970) 669-2932
orders@hach.com
www.hach.com

HACH LANGE GMBH

Willstätterstraße 11
D-40549 Düsseldorf, Germany
Tel +49 (0) 2 11 52 88-320
Fax +49 (0) 2 11 52 88-210
info@hach-lange.de
www.hach-lange.de

HACH LANGE Sàrl

6, route de Compois
1222 Vésenaz
SWITZERLAND
Tel. +41 22 594 6400
Fax +41 22 594 6499