
M231Pt-9/M251Pt-9
Platinum Electrodes
Electrodes de Platine

M231Pt-9 / M251Pt-9 Platinum Electrodes

Introduction

The M231Pt-9 and M251Pt-9 Platinum Electrodes can be used with a standard Ag/AgCl reference electrode or, depending on the application, with a reference electrode which has an alternative liquid junction or reference system.

Preparation for Measurement

1. Remove the protection cap from the electrode.
2. Before starting a measurement, rinse the electrode with distilled water. Take care not to wet the screw cap.
3. Check that the electrode body is not damaged.

Maintenance

Platinum rod

The sensing part of the M231Pt-9 electrode can be polished with finishing abrasive bands, part no. B90X550. While polishing, take great care not to fold the platinum rod.

In order to have reproducible and accurate results with the electrode, a few physical and electrochemical treatments are required. These treatments will depend on the last measurement performed:

Mild cleaning

Place the electrode overnight in RENOVO•N Normal Cleaning Solution or place the electrode 15 minutes in a Triton X-100 solution.

Degreasing

- Remove mineral oils by rinsing the electrode with trichlorethylene and dabbing dry with a soft tissue.
- Remove organic oils by cathodic electrolysis in a hot alkaline solution, e.g. 15 g/l NaOH + 50 g/l Na₂CO₃ at 25°C.

Degassing

- Hydrogen produced during the electrolysis of the alkaline solution (see above) must be removed. Dip the electrode in water at 60°C and flush with pure nitrogen.

Chemical cleaning

- The electrode can be cleaned with a sulphochromic solution, or with a boiling solution of hydrochloric acid.

Storage

Replace the protection cap, and store the electrode in its box at room temperature.

Accessories

BAO3 Finishing Abrasive Bands	B90X550
○RENOVO•N	
Normal Cleaning Solution, 250 ml	S16M001
CL114 Electrode Cable, (1 m) with BNC plug	A94L114
CL116 Electrode Cable, (1 m) with type 7 plug	A94L116
CL111 Electrode Cable, (1 m) with banana plug	A94L111

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