
DOC022.97.80156

# senSIOn™ + MM378

07/2013, Edition 1

**User Manual**  
**Manuel d'utilisation**  
**Manual del usuario**  
**Manual do Usuário**

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## Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions	35 x 20 x 11 cm (13.78 x 7.87 x 4.33 in.)
Weight	1100 g (2.43 lb)
Meter enclosure	IP42
Power requirements (external)	100–240 V, 0.4 A, 47-63 Hz
Meter protection class	Class II
Pollution degree	2
Installation category	Category II
Altitude requirements	Standard 2000 m (6562 ft) ASL (Above Sea Level)
Storage temperature	–15 to +65 °C (5 to +149 °F)
Operating temperature	0 to 40 °C (41 to 104 °F)
Operating humidity	< 80% (non-condensing)
Resolution	pH: 0.1/0.01/0.001, ORP: 0.1/1 mV, ISE: programmable, temperature: 0.1 °C (0.18 °F), EC: variable, Resistivity: variable, NaCl: variable, TDS: variable, DO: 0.01 mg/L (from 0.00 to 19.99 mg/L) 0.1 mg/L (from 20.0 to 60.0 mg/L) 0.1% (from 0.0 to 19.9%) 1% (from 20 to 600%)
Measuring error (± 1 digit)	pH: ≤ 0.005, ORP: ≤ 0.2 mV, temperature: ≤ 0.2 °C (≤ 0.36 °F), EC: ≤ 0.5%, resistivity: ≤ 0.5%, NaCl: ≤ 0.5%, TDS: ≤ 0.5%, DO: ≤ 0.5%
Reproducibility (± 1 digit)	pH: ± 0.001, ORP: ± 0.1 mV, temperature: ± 0.1 °C (± 0.18 °F), EC: ± 0.1%, resistivity: ± 0.1%, NaCl ± 0.1%, TDS ± 0.1%, DO ± 0.2%
Data storage	330 results and last 9 calibrations

Specification	Details
Connections	1 Combined or indicator probe: BNC connector (Imp. >10 <sup>12</sup> Ω); 1 Reference electrode: banana connector; 2 A.T.C. type Pt 1000 (or NTC 10 kΩ probe): 1 banana and 1 telephonic connector; 2 magnetic stirrers: RCA connector  Conductivity probe with built-in Pt1000 sensor (or NTC 10 kΩ probe): telephonic connector RS232C for printer or PC: telephonic connector; external PC keyboard: mini DIN connector Dissolved Oxygen probe: BNC connector, with built-in NTC 30 kΩ probe: banana connector
Temperature correction	<b>Channel 1 pH:</b> Pt 1000 temperature probe (A.T.C.), NTC 10 kΩ probe, manual, isopotential pH programmable, standard value 7.00, <b>EC:</b> Pt 1000 (or NTC 10KΩ) temperature probe (A.T.C.), linear function, TC=0.00 to 9.99%/Temperature. TRef: 20°C (68 °F), 25°C (77 °F) or other temperature (between 0 and 35 °C), non-linear function for natural waters (UNE EN 2788) <b>Channel 2 DO:</b> NTC 30 KΩ probe, manual
Measurement display lock	Continuous measurement, by stability, by time
Display	Liquid crystal, backlit, 128 x 64 dots
Keyboard	PET with protective treatment
Certification	CE

## General information

Revised editions are found on the manufacturer's website.

## Safety information

### NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

## Use of hazard information

### ▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

### ▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

### ▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. <i>Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.</i>

## Product overview

The sensION™+ meters are used with probes to measure various parameters in water.

The sensION™+ MM378 meter has two measurement channels to measure pH, ORP (mV), conductivity, ISE (concentration) or dissolved oxygen. Channel 1 measures one or two parameters individually or simultaneously. Up to two probes can be connected to Channel 1. Channel 2 measures dissolved oxygen. Measurement data can be stored and transferred to a printer or PC.

## Certification

### Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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## FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.


Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

## Product components

Refer to [Meter components](#) on page 5 to make sure that all components have been received. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

**Figure 1 Meter components**





1 pH probe electrolyte solution (KCl)	7 Power supply
2 DO probe filling solution (included with kits only)	8 Probe holder
3 pH buffer standard solutions (pH 4.01 and pH 7.00)	9 Rod
4 Calibration beakers (with magnetic bar inside)	10 Probes (included with kits only)
5 Conductivity standard solution (1413 $\mu\text{S}/\text{cm}$ )	11 Meter
6 Calibration beaker (with magnetic bar inside)	

# Installation

## Assemble the probe holder

Follow the numbered steps to assemble the probe holder and to connect the magnetic stirrer.

## Connect to AC power

### ⚠ DANGER


Electrocution hazard. If this equipment is used outdoors or in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used to connect the equipment to its main power source.

The meter can be powered by AC power with the universal power adapter.


1. Select the correct adapter plug for the power outlet from the adapter kit.
2. Connect the universal power adapter to the meter ([Figure 2](#)).
3. Connect the universal power adapter to an AC receptacle ([Figure 3](#)).
4. Turn the meter on.

**Figure 2 Connector panel**

1 Magnetic stirrer 1 connector, Channel 1	7 PC keyboard, mini DIN connector
2 Temperature probe connector, Channel 2	8 Conductivity probe connector, Channel 1
3 Reference electrode (separated electrodes) connector, Channel 1	9 Temperature probe connector, Channel 1
4 Separated temperature probe connector, Channel 1	10 BNC combined (or indicator) pH probe connector, Channel 1
5 RS-232 for printer or PC connector	11 BNC DO probe connector, Channel 2
6 Power supply	12 Magnetic stirrer 2 connector, Channel 2

**Figure 3 AC power connection**


**Note:** Position the equipment so it is not difficult to disconnect power to the equipment.

# User interface and navigation

## User interface

### Keypad description


1 RETURN key: cancel or exit the current menu screen to the previous menu screen	5 DOWN key: scroll to other options, change a value
2 MEASUREMENT key: confirm the selected option	6 ON/OFF: turn on or turn off the meter
3 UP key: scroll to other options, change a value	7 LEFT key: change between Channel 1 and Channel 2, enter numbers and letters
4 RIGHT key: change between Channel 1 and Channel 2, enter numbers and letters	

### Display description

The meter display shows the concentration, units, temperature, calibration status, operator ID, sample ID, date and time.

Figure 4 Single screen display

1 Sample ID	5 Measurement channel
2 Measurement unit and value (conductivity)	6 Sample temperature (°C or °F)
3 Measurement unit and value (pH, ORP/mV or ISE)	7 Reference temperature
4 Measurement mode or time and date	8 Visual measurement timer

### Navigation

Use the to return to the previous menu. Use the measure key to take a sample measurement or to confirm options. Use the arrow keys to scroll to other options or to change a value. To change the parameters use the arrow keys and . Refer to each task for specific instructions.

### Startup

#### Turn the meter on and off

**NOTICE**

Make sure that the probe is connected to the meter before the meter is turned on.

Push to turn on or turn off the meter. If the meter does not turn on, make sure that the AC power supply is properly connected to an electrical outlet.

## Change the language

The display language is selected when the meter is powered on for the first time.

1. Use the ▲ or ▼ to select a language from a list.
2. Confirm with ✓. The measurement screen shows DATA OUTPUT.
3. Select Deactivated if no printer or PC is connected and confirm. Refer to [Select the data output](#) on page 13 for more information about Data Output.

## Standard operation

### Setup for dissolved oxygen

Before a calibration or a measurement for dissolved oxygen, the probe must be polarized and the settings for atmospheric pressure and salinity must be entered.

#### Polarize the electrode

If the probe or batteries are disconnected, connect the probe or install the batteries and wait for polarization:

Disconnect time	Polarization time
< 5 minutes	10 minutes
5 to 15 minutes	45 minutes
> 15 minutes	6 hours

## Calibration

### ⚠ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

#### Calibration settings

The calibration settings contain Calibration type, Calibration frequency and Display options.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select the following options:

Option	Description
<b>Stability C.:</b>	Criterion by stability—select Fast, Standard or Strict.
<b>Calibration type</b>	<b>pH:</b> Calibration type—select Technical buffers, DIN19266 Buffers, User Buffers, To a X value, data introduction or Theoretical calibration. <b>EC:</b> Calibration type—select Molar Standards, Demal Standards, NaCl Standards, Calibration to a X value, data introduction or Theoretical calibration. Refer to <a href="#">Calibration types</a> for more information. <b>DO:</b> Calibration type—select Water saturated air, Water saturated air & Zero, To a X value or Theoretical calibration.
<b>Cal. frequency</b>	Calibration reminder—can be set between 0-23 hours (pH) or 0-7 days (EC). The display shows the remaining time to the new calibration. Refer to <a href="#">Set the calibration reminder</a> on page 12 reminder for more information.
<b>Display mV</b>	Display mV—select YES or NO to show mV.
<b>Display nA</b>	Display nA—select YES or NO to show DO current in nA (nanoampere).

#### Calibration types

Different calibration types can be selected.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select Calibration type.

Option	Description
<b>Technical buffers</b>	pH 2.00, 4.01, 7.00, 9.21 and 10.01 at 25 °C (77 °F)

Option	Description
<b>DIN19266 Buffers</b>	pH 1.679, 4.006, 6.865, 7.000, 9.180, 10.012 and 12.454
<b>User Buffers</b>	Selected when the technical or DIN19266 buffers are not used. Refer to <a href="#">Technical buffer solutions (DIN 19267)</a> on page 20 for pH values of specific buffer sets at varying temperatures.
<b>Calibration to a X value</b>	To adjust manually any scale value of the measured pH or conductivity.
<b>Data introduction</b>	Manual probe constant introduction.
<b>Theoretical calibration</b>	<b>pH:</b> The probe calibration data is replaced at 25 °C (77 °F). <b>EC:</b> The probe calibration data is replaced with $C=1.000 \text{ cm}^{-1}$ <b>DO:</b> calibration is based on theoretical values (no physical calibration)
<b>Molar Standards</b>	147 µS/cm, 1413 µS/cm, 12.88 mS/cm and 111.8 mS/cm at 25 °C (77 °F)
<b>Demal Standards</b>	1049 µS/cm, 12.85 mS/cm and 111.31 mS/cm at 25 °C (77 °F)
<b>NaCl Standards</b>	1014.9 µS/cm 25 °C (77 °F)
<b>DO Standards</b>	1 point calibration: water saturated air 2 point calibration: water saturated air & zero

## Calibration procedure

This procedure is for general use with liquid calibration solutions. Refer to the documents that are included with each probe for additional information.

**Note:** Solutions should be stirred during calibration. For more information about the stirring settings, refer to [Change the stirring settings](#) on page 14.

1. Pour the buffers or calibration solutions into the labeled calibration beakers.
2. From the main menu use the ▲ or ▼ and ◀ and ▶ to select the CALIBRATION parameter. Confirm.
3. If required select the Operator ID (1 to 10) and confirm.

4. Push ✓ to Start calibration.
5. Push ✓ to measure the first calibration solution. The next calibration solution is shown.
6. Rinse the probe with deionized water and put the probe into the second calibration beaker. Be sure that there are no air bubbles in the membrane.
7. Push ✓ to measure the second calibration solution. The next calibration solution is shown.
8. Rinse the probe with deionized water and put the probe into the third calibration beaker. Be sure that there are no air bubbles in the membrane.
9. Push ✓ to measure the third calibration solution. When the calibration is good, the display briefly shows Calibration OK and then returns to the main menu.  
**Note:** When a printer is connected the print menu opens and the result can be printed.

## View the calibration data

Data from the most recent calibration can be shown.

1. From the main menu use the ▲ or ▼ to select DATA LOGGER. Confirm.
2. Select Display data.
3. Select Calibration data and confirm with ✓. The last calibration data is shown.
  - pH-the slope and offset values are shown alternating with the deviation (in %) and calibration temperature.
  - ORP-the measured mV value and calibration temperature are shown.
  - Conductivity-the cell constant and calibration temperature for each standard are shown.
  - DO-the current (in nA) and temperature are shown.

## Set the calibration reminder

**pH:** The calibration reminder can be set between 0 to 23 hours or 1-7 days (default 1 day). **EC:** The calibration reminder can be set between 0 to 99 days (default 15 days). The display shows the remaining time to the new calibration.

**Note:** When 0 days is selected, the calibration reminder is turned off.

1. From the main menu use the ▲ or ▼ to select CALIBRATION. Confirm.
2. Use the ▲ to enter the calibration menu.
3. Use the ▲ or ▼ to select Cal. frequency and confirm.
4. Use the ◀ and ▶ to advance to the next step and use the ▲ or ▼ to change a value. Confirm.  
Push ✓ to start the calibration.

## Sample measurements

Each probe has specific preparation steps and procedures for taking sample measurements.

1. From the main menu use the ▲ or ▼ and ◀ and ▶ to select MEASURE. Confirm.
2. Use the ▲ to change the following settings. Confirm every entry.

Option	Description
<b>Resolution</b>	Select the resolution: 1, 0.1, 0.01 (default) or 0.001
<b>Measure</b>	Stability—select By stability Criterion: Fast (variation < 0.02 pH in 6 s), Standard (variation < 0.01 pH in 6 s) or Strict (variation < 0.002 pH in 6 s). In continuous—enter the time interval for the In continuous Acquis. interval (data storage or printing data). By time—enter the time interval for data storage or printing data.
<b>Display mV</b>	Display mV—select YES or NO to show mV.
<b>Limits</b>	Limits—select YES or NO. YES: Enter the upper and lower limits. An acoustic warning appears when the measurement is out of limit. The report output shows an A next to the measured value when the measurement was out of limit.

Option	Description
<b>Isopotential</b>	Isopotential—change the Isopotential pH value in Data introduction. Select Calculate to calibrate the probe again.
<b>TC</b>	TC—select Linear or Natural waters. Linear: enter a value in %/°C (default 2.00%/Temperature). Natural waters: Non-linear for natural waters, according to EN27888
<b>Tref</b>	Reference temperature—select between 20 or 25 °C or Other temperature.
<b>Pressure</b>	Barometric pressure—select a pressure between 600 and 1133 mBar (450 mmHg and 850 mmHg). Preset to 1013 mBar
<b>Sal.Corr.</b>	Salinity correction—enter a correction factor between 00.0 and 45.0 g/L (preset to 1.0 g/L)
<b>Display nA</b>	Display nA—select YES or NO to show DO current in nA (nanoampere)

**Note:** The salinity value can be calculated from the conductivity value. Refer to [Table 7](#) on page 23

3. Push ✓ to start the measurement.

**Note:** If the measurement is not stabilizing after 120 seconds, the meter turns automatically into the continuous measurement mode.

## Advanced operation

### Change the measurement units

The measurement units can be changed individually for each channel.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Measurement units and confirm.
3. Select Channel 1 and confirm.

- Select Parameter 1 and then Parameter 2 and select one of the following options:

Option	Description
Parameter 1	mV, pH, ISE or Disabled
Parameter 2	EC, NaCl, TDS, $\Omega$ or Disabled

- Select Channel 2 and confirm with . Select % or mg/L and confirm.

## Use a sample ID

The sample ID tag is used to associate readings with a particular sample location. If assigned, stored data will include this ID.

- From the main menu use the or to select SYSTEM. Confirm.
- Use the or to select Sample ID and confirm.
- Use the or to select

Option	Description
<b>Automatic</b>	A consecutive number will be automatically assigned to every sample.
<b>Manual</b>	A keyboard or a barcode scanner is required to enter the sample ID name before taking a measurement (maximum 15 characters).

## Select the data output

Data can be stored or transferred to a printer or to a PC.

- From the main menu use the or to select SYSTEM. Confirm.
- Use the or to select Data Output and confirm.
- Use the or to select

Option	Description
<b>Deactivated</b>	Select Deactivated if no printer or PC is connected.

Option	Description
<b>For Printer</b>	Select Dot matrix printer or Thermal printer.
<b>For Computer</b>	Select Terminal, LabCom or LabCom Easy. The LabCom Software controls several modules, pH and conductivity meters, automatic burettes, Samplers and so on from a computer. The LabCom Easy software gets measurements and calibration data from the meters

*Note: Make sure that LabCom or LabCom Easy software has version 2.1 to be compatible with the instrument.*

## Change the date and time

The date and time can be changed from the Date / Time menu.

- From the main menu use the or to select SYSTEM. Confirm.
- Use the or to select Date / Time and confirm.
- Use the and to advance to the next step and use the or to change a value. Confirm.  
The current date and time will be shown on the display.

## Adjust the display contrast

- From the main menu use the or to select SYSTEM. Confirm.
- Use the or to select Display contrast and confirm.
- Use the and to adjust the contrast of the display and confirm.

## Adjust the temperature

The temperature measurement can be adjusted at 25 °C (77 °F) and/or 85 °C (185 °F) to increase accuracy.

- Put the probe and a reference thermometer in a container of water at approximately 25 °C and allow the temperature to stabilize.
- Compare the temperature read by the meter with that of the reference thermometer. The difference is the adjustment value for the meter.

Example: reference thermometer: 24.5 °C; meter: 24.3 °C.  
Adjustment value: 0.2 °C.

3. Enter the adjustment value for the 25 °C reading:
  - a. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
  - b. Use the ▲ or ▼ to select Readjust temp. and confirm.
  - c. Select Channel 1 or Channel 2 and confirm.
  - d. Use the ▲ or ▼ to select 25 °C and confirm.
  - e. Use the arrow keys to enter the adjustment value for 25 °C. Confirm.

Option	Description
--------	-------------

**Channel 1** Select Parameter 1 or Parameter 2. **Parameter 1** the temperature is measured by the pH probe. **Parameter 2** the temperature is measured by the conductivity cell. If one of the parameters is not used, the meter automatically assigns the temperature of the other parameter. If no A.T.C. is detected, the temperature value must be entered manually

**Channel 2** The temperature for Channel 2 is measured by the DO probe or manually entered.

4. Put the probe and a reference thermometer in a container of water at approximately 85 °C and allow the temperature to stabilize.
5. Compare the temperature from the meter with that of the reference thermometer. The difference is the adjustment value for the meter.
  - a. Use the ▲ or ▼ to select 85 °C and confirm.
  - b. Use the arrow keys to enter the adjustment value for 85 °C. Confirm.
  - c. Select Save changes and confirm.

## Change the stirring settings

The stirring speed can be changed during calibration and during a measurement.

1. Use the ▲ or ▼ to change the stirring speed during calibration and during a measurement.

## Turn the stirrer on or off

The stirrer 1 works with Channel 1 and 2 (Stirrer 1). A second stirrer can be connected to channel 2 (Stirrer 2). For the stirrer 2 activation refer to the following steps.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm with ✓.
2. Use the ▲ or ▼ to select Stirrer N.2 and confirm with ✓.
3. Use the ▲ or ▼ to select YES to turn the stirrer 2 on.  
*Note: Select NO to turn the stirrer 2 off.*

## Change the temperature units

The temperature units can be changed to Celsius or Fahrenheit.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Temperature units and confirm.
3. Use the ▲ or ▼ to select between Celsius or Fahrenheit and confirm.

## Data logger

### Display data

The Display data log contains Measurement data, Electrode report and Calibration data. The stored data can be send to a printer or to a PC. When the data log becomes full (400 data points), the oldest data point is deleted when a new data point is added.

1. Use the ◀ and ▶ to select Channel 1 or Channel 2 and confirm.
2. From the main menu use the ▲ or ▼ to select DATA LOGGER and confirm.
3. Use the ▲ or ▼ to select Display data and confirm.
4. Use the ▲ or ▼ to select

Option	Description
Measurement data	Measurement data—stores automatically each time a sample is measured

Option	Description
<b>Electrode report</b>	Electrode report—stores automatically the electrode history and measurement conditions
<b>Calibration data</b>	Calibration data—stores automatically the current calibration

## Delete data

The entire measurement data or electrode report log can be deleted to remove data that has already been sent to a printer or PC.

1. From the main menu use the ▲ or ▼ to select DATA LOGGER and confirm.
2. Use the ▲ or ▼ to select Erase and confirm.
3. Use the ▲ or ▼ to select Measurement data or Electrode report and confirm. Confirm again to delete the data.  
The entire log is deleted at once.

## Send data to a printer or to a computer

### NOTICE

The data output (printer or PC) needs to be selected first, so that the Print menu is available (refer to [Select the data output](#) on page 13).

**Note:** Refer to [Report output](#) on page 15 to select the report output type.

1. From the main menu use the ▲ or ▼ to select DATA LOGGER and confirm.
2. Use the ▲ or ▼ to select Print and confirm. Select one of the following options and confirm with ✓ to print the data: Measurement data, Electrode data, Calibration data, Calibration report or Instrument condit.

## Report output

### NOTICE

The data output (printer or PC) needs to be selected first, so that the Type of report menu is available (refer to [Select the data output](#) on page 13).

Different report output types can be selected when a printer or a PC is connected.

1. From the main menu use the ▲ or ▼ to select SYSTEM. Confirm.
2. Use the ▲ or ▼ to select Type of report and confirm.
3. When a printer or a computer and Terminal is connected, use the ▲ or ▼ to select

### Option Description

**Reduced** Select Several or One sample as an output format

**Standard** Select Several or One sample as an output format. Select several: **Users:** The user name appears on the printed report (17 characters). **Header:**The company name can be added as a header (40 characters) and appears on the print report. **Identify sensor:**The sensor model and the sensor serial number can be added and appears on the print report.

**GLP** Select Several or One sample as an output format. Select several: **Users:** The user name appears on the printed report (17 characters). **Header:**The company name can be added as a header (40 characters) and appears on the print report. **Identify sensor:**The sensor model and the sensor serial number can be added and appears on the print report.

4. When a computer is connected and LabCom Easy (refer to [Select the data output](#) on page 13 for more information) is selected, use the ▲ or ▼ to select

### Option Description

**Users** The user name appears on the printed report (17 characters).

**Identify sensor** The sensor model and the sensor serial number can be added and appears on the print report.

## Maintenance

### **▲ WARNING**

Multiple hazards. Do not disassemble the instrument for maintenance or service. If the internal components must be cleaned or repaired, contact the manufacturer.

### **▲ CAUTION**

Personal injury hazard. Only qualified personnel should conduct the tasks described in this section of the manual.

## Clean the instrument

### **NOTICE**

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Clean the exterior of the instrument with a moist cloth and a mild soap solution.



## Clean the probe

Clean the probe as needed. Refer to [Troubleshooting](#) on page 18 for more information about cleaning. Refer to the probe documentation for information about the probe maintenance.

## Replace the magnetic stirrer

If the magnetic stirrer does not start, follow the numbered steps to replace the magnetic stirrer.

## Troubleshooting

Refer to the following table for common problem messages or symptoms, possible causes and corrective actions.

**Table 1 Calibration warnings and errors**

Error/Warning	Solution
Calibration out of range	Measured value out of range. Calibrate again. Connect a new probe.
Unknown buffer	Calibrate again.
Same buffer / standard. Re-calibrate	Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.  Examine the buffer solution: Make sure that the buffer used matches the buffer specified in configuration; make sure of the temperature specification in configuration; use a new buffer solution.
Unstable measurement Time > 100 s (pH, EC and DO calibrations) Time > 240 s (ISE calibrations)	Calibrate again.  Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.  Make sure that the membrane and the diaphragm are properly immersed in the sample.
Temperature difference > 3.0 °C	Adjust the calibration solutions to the same temperature. Examine the temperature sensor.
Temperature out of range	Examine the temperature sensor. Connect a new probe.

**Table 1 Calibration warnings and errors (continued)**

Error/Warning	Solution
Outside allowable range	Offset or Slope outside the range. Examine the buffer solution: Make sure that the buffer used matches the buffer specified in configuration; make sure of the temperature specification in configuration; use a new buffer solution.  Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Signal too low / high (DO calibrations)	DO probe error. Examine the probe. Connect a different probe to make sure that if problem is with probe or meter.  Use a new standard solution.
Cell constant over limits (EC calibrations)	Insert the probe into the appropriate standard and re-read. Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Cell constant deviation error (EC calibrations)	Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Not calibrated	There are no calibration data stored in the instrument. Perform calibration.

**Table 2 Measurement warnings and errors**

Error/Warning	Solution
-----	Measured value is out of range. Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to find if the problem is with the probe or the meter.
Unstable measurement (pH, EC and DO measurements) Time > 120 s	Make sure that the membrane and the diaphragm are properly immersed in the sample. Examine the temperature. Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to make sure that if problem is with probe or meter.
Unstable measurement (ISE measurements) Time > 240 s	Make sure that the membrane and the diaphragm are properly immersed in the sample. Examine the temperature. Examine the probe: Clean the probe (refer to <a href="#">Clean the probe</a> on page 16 for more information); make sure that there are no air bubbles in the membrane. Shake the probe like a thermometer; connect a different probe to make sure that if problem is with probe or meter.
TC = 0 Impossible to measure salinity	Modify TC
10.389 mg/L >>>>> 00012 00:13	ISE measurement: The measured value is superior to the highest standard used in calibration. Measure again.
0.886 mg/L <<<<<< 00018 00:11	ISE measurement: The measured value is inferior to the lowest standard used in calibration. Measure again.

## Replacement parts and accessories

*Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.*

### Replacement parts

Description	Item no.
sensION+ PH3 Lab pH-meter with accessories, without probe	LPV2000.97.0002
sensION+ PH31 Lab pH-meter, GLP, with accessories, without probe	LPV2100.97.0002
sensION+ MM340 Lab pH & Ion-meter, GLP, 2 channels, with accessories, without probe	LPV2200.97.0002
sensION+ EC7 Lab conductivity meter, with accessories, without probe	LPV3010.97.0002
sensION+ EC71 Lab conductivity meter, GLP, with accessories, without probe	LPV3110.97.0002
sensION+ MM374, 2 channel Lab meter, GLP, accessories, without probes	LPV4110.97.0002
sensION+ MM378, 2 channel Lab meter, GLP, accessories, without probes	LPV4130.97.0002

### Consumables

Description	Item no.
pH buffer solution 4.01, 250 mL	LZW9463.99
pH buffer solution 7.00, 250 mL	LZW9464.97
pH buffer solution 10.01, 250 mL	LZW9471.99
Electrolytic solution (KCl 3M), 125 mL	LZW9510.99
Electrolytic solution (KCl 3M), 250 mL	LZW9500.99
Electrolytic solution (KCl 3M), 50 mL	LZW9509.99

## Consumables (continued)

Description	Item no.
Electrolytic solution (DO), 50 mL	2759126
Conductivity standard 147 µS/cm, 250 mL	LZW9700.99
Conductivity standard 1413 µS/cm, 250 mL	LZW9710.99
Conductivity standard 12.88 mS/cm, 250 mL	LZW9720.99
Electrolytic solution 0.1 M, 125 mL	LZW9901.99
Pepsin Cleaning Solution	C20C370
Renovo N, electrode cleaning solution	S16M001
Renovo X, 0.1 N HCl solution	S16M002
Ethanol, 95% (Grease, oils, fats)	2378900

## Accessories

Description	Item no.
Thermal printer, USB, for sensION+ benchtop instruments	LZW8203.97
Thermal paper for printer LZW8203, bag with 4 rolls	LZM078
Power supply for sensION+ benchtop instruments, 230-115 VAC	LZW9008.99
LabCom Easy PC SW, for sensION+ GLP, CD, cable, USB adapter	LZW8997.99
LabCom PC SW, for sensION+ GLP, CD, cable, USB adapter	LZW8999.99
Magnetic stirrer with sensor holder, for sensION+ multimeters	LZW9319.99
3x50 mL printed beakers for pH calibration	LZW9110.97
3x50 mL printed beaker for conductivity calibration	LZW9111.99
Three-sensor holder, for sensION+ benchtop instruments	LZW9321.99

## Accessories (continued)

Description	Item no.
Radiometer probe holder for sensION+ benchtop instruments	LZW9325.99
Pyrex glass chamber, continuous flow measurements	LZW9118.99
PP protector, electrode storage	LZW9161.99
Teflon coated stirring bar, 20 x 6 mm	LZW9339.99
Protector flask for DO probe	LZW5123.99
Replacement membrane for 5131 DO probe	5197300

## Standard solutions

### Technical buffer solutions (DIN 19267)

Refer to [Table 3](#) pH and ORP (mV) values of specific buffer sets at varying temperatures.

**Table 3 pH, ORP (mV) and temperature values**

Temperature		pH					mV
°C	°F						
0	32	2.01	4.01	7.12	9.52	10.30	—
10	50	2.01	4.00	7.06	9.38	10.17	245
20	68	2.00	4.00	7.02	9.26	10.06	228
<b>25</b>	<b>77</b>	<b>2.00</b>	<b>4.01</b>	<b>7.00</b>	<b>9.21</b>	<b>10.01</b>	<b>220</b>
30	86	2.00	4.01	6.99	9.16	9.96	212
40	104	2.00	4.03	6.97	9.06	9.88	195
50	122	2.00	4.06	6.97	8.99	9.82	178
60	140	2.00	4.10	6.98	8.93	9.76	160
70	158	2.01	4.16	7.00	8.88	—	—

**Table 3 pH, ORP (mV) and temperature values (continued)**

Temperature		pH					mV
°C	°F						
80	176	2.01	4.22	7.04	8.83	—	—
90	194	2.01	4.30	7.09	8.79	—	—

**Buffer solutions (DIN 19266)**

Refer to [Table 4](#) for pH values of specific buffer sets at varying temperatures.

**Table 4 pH and temperature values**

Temperature		pH						
°C	°F							
5	32	1.668	4.004	6.951	7.087	9.395	10.245	13.207
10	50	1.670	4.000	6.923	7.059	9.332	10.179	13.003
20	68	1.675	4.001	6.881	7.016	9.225	10.062	12.627
<b>25</b>	<b>77</b>	<b>1.679</b>	<b>4.006</b>	<b>6.865</b>	<b>7.000</b>	<b>9.180</b>	<b>10.012</b>	<b>12.454</b>
30	86	1.683	4.012	6.853	6.987	9.139	9.966	12.289
40	104	1.694	4.031	6.838	6.970	9.068	9.889	11.984
50	122	1.707	4.057	6.833	6.964	9.011	9.828	11.705
60	140	1.723	4.085	6.836	6.968	8.962	—	11.449
70	158	1.743	4.126	6.845	6.982	8.921	—	—
80	176	1.766	4.164	6.859	7.004	8.885	—	—
90	194	1.792	4.205	6.877	7.034	8.850	—	—

**Conductivity standard solutions**

Refer to [Table 5](#) for conductivity values of standard solutions at varying temperatures.

**Table 5 Conductivity and temperature values**

Temperature		Conductivity (EC)			
°C	°F	µS/cm	µS/cm	mS/cm	mS/cm
15.0	59	119	1147	10.48	92.5
16.0	60.8	122	1173	10.72	94.4
17.0	62.6	125	1199	10.95	96.3
18.0	64.4	127	1225	11.19	98.2
19.0	66.2	130	1251	11.43	100.1
20.0	68	133	1278	11.67	102.1
21.0	69.8	136	1305	11.91	104.0
22.0	71.6	139	1332	12.15	105.4
23.0	73.4	142	1359	12.39	107.9
24.0	75.2	145	1386	12.64	109.8
<b>25.0</b>	<b>77</b>	<b>147</b>	<b>1413</b>	<b>12.88</b>	<b>111.8</b>
26.0	78.8	150	1440	13.13	113.8
27.0	80.6	153	1467	13.37	115.7
28.0	82.4	156	1494	13.62	—
29.0	84.2	159	1522	13.87	—
30.0	86	162	1549	14.12	—
31.0	87.8	165	1581	14.37	—
32.0	89.6	168	1609	14.62	—
33.0	91.4	171	1638	14.88	—
34.0	93.2	174	1667	15.13	—
35.0	95	177	1696	15.39	—

**Demal (D) and NaCl 0.05% EC standards**

Refer to [Table 6](#) for conductivity values according to temperature.

**Table 6 Conductivity and temperature values**

Temperature		KCl 1D (mS/cm)	KCl 0.1D (mS/cm)	KCl 0.01D (µS/cm)	NaCl 0.05% (µS/cm)
°C	°F				
0	32	65.14	7.13	773	540.40
1	33.8	66.85	7.34	796	557.73
2	35.6	68.58	7.56	820	575.20
3	37.4	70.32	7.77	843	592.79
4	39.2	72.07	7.98	867	610.53
5	41	73.84	8.20	891	628.40
6	42.8	75.62	8.42	915	646.40
7	44.6	77.41	8.64	940	664.55
8	46.4	79.21	8.86	965	682.83
9	48.2	81.03	9.08	989	701.26
10	50	82.85	9.31	1014	719.82
11	51.8	84.68	9.54	1039	738.53
12	53.6	86.54	9.76	1065	757.37
13	55.4	88.39	9.99	1090	776.36
14	57.2	90.26	10.22	1116	795.48
15	59	92.13	92.13	1142	814.74
16	60.8	94.02	10.69	1168	834.14
17	62.6	95.91	10.93	1194	853.68
18	64.4	97.81	11.16	1220	873.36
19	66.2	99.72	11.40	1247	893.18
20	68	101.63	11.64	1273	913.13

**Table 6 Conductivity and temperature values (continued)**

Temperature		KCl 1D (mS/cm)	KCl 0.1D (mS/cm)	KCl 0.01D (µS/cm)	NaCl 0.05% (µS/cm)
°C	°F				
21	69.8	103.56	11.88	1300	933.22
22	71.6	105.49	12.12	1327	953.44
23	73.4	107.42	12.36	1354	973.80
24	75.2	109.36	12.61	12.61	994.28
<b>25</b>	<b>77</b>	<b>111.31</b>	<b>12.85</b>	<b>1409</b>	<b>1014.90</b>
26	78.8	113.27	13.10	1436	1035.65
27	80.6	115.22	13.35	1464	1056.53
28	82.4	—	13.59	1491	1077.54
29	84.2	—	13.84	1519	1098.67
30	86	—	14.09	1547	1119.92
31	87.8	—	14.34	1575	1141.30
32	89.6	—	14.59	1603	1162.80
33	91.4	—	14.85	1632	1184.41
34	93.2	—	15.10	1660	1206.15
35	95	—	15.35	1688	1228.00
36	96.8	—	15.61	1717	1249.96
37	98.6	—	15.86	1745	1272.03
38	100.4	—	16.12	1774	1294.96
39	102.2	—	16.37	1803	1316.49
40	104	—	16.63	1832	1338.89
41	105.8	—	16.89	1861	1361.38
42	107.6	—	17.15	1890	1383.97

**Table 6 Conductivity and temperature values (continued)**

Temperature		KCl 1D (mS/cm)	KCl 0.1D (mS/cm)	KCl 0.01D (μS/cm)	NaCl 0.05% (μS/cm)
°C	°F				
43	109.4	—	17.40	1919	1406.66
44	111.2	—	17.66	1948	1429.44
45	113	—	17.92	1977	1452.32
46	114.8	—	18.18	2007	1475.29
47	116.6	—	18.44	2036	1498.34
48	2065	—	118.4	18.70	1521.48
49	120.2	—	18.96	2095	1455.71
50	122	—	19.22	2124	1568.01

**Table 7 Conversion from conductivity to salinity (continued)**

Conductivity (mS/cm)	Salinity (g/L NaCl)
56.3	35.8
62.2	40.0
69.4	45.3
75.7	50.0

## Salinity values for dissolved oxygen meters

Table 7 can be used to convert conductivity values to salinity for dissolved oxygen meters.

**Table 7 Conversion from conductivity to salinity**

Conductivity (mS/cm)	Salinity (g/L NaCl)
1.9	1.0
3.7	2.0
7.3	4.0
10.9	6.0
17.8	10.1
25.8	15.1
33.6	20.2
41.2	25.4
48.9	30.6

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)