

In-Line Methods List for Automated Ion Analyzers

Flow Injection Analysis



10 June 2019

In-line Modules:

Method Number	Analyte	Lamp/Tubing	Heater/tubing
10-107-06-6-A	Ammonia	Not used	850cm Tefzel* (PN 41300L)
10-107-06-6-B	Ammonia	Not used	850cm Tefzel* (PN 41300L)
10-106-06-6-E	Ammonia	Not used	745cm 0.8mm i.d. (PN 50028)
12-140-39-5-A	Carbon	6W germicidal/810 zeus	1200 cm 0.8"i.d. (PN 50028)
10-204-00-2-C	Cyanide (total)	6W, BLB/ 750cm zeus	850cm Tefzel (PN 41300L)
10-204-00-2-D	Cyanide (total)	6W, BLB/ 750cm zeus	850cm Tefzel (PN 41300L)
10-204-00-2-E	Cyanide (total)	6W, BLB/ 810cm zeus	850cm Tefzel (PN 41300L)
10-204-00-2-G	Cyanide (total)	6W, BLB/ 700cm zeus	850cm Tefzel (PN 41300L)
10-204-00-2-H	Cyanide (total)	8W, BLB/550cm zeus	475cm 0.8mm i.d.* (PN 50028)
10-204-00-2-I	Cyanide (total)	6W, BLB/ 700cm zeus	850cm Tefzel (PN 41300L)
10-204-00-4-B	Cyanide (WAD)	Not used	1200 cm 0.8"i.d. (Adjacent channel) (PN 50028)
10-204-00-4-C	Cyanide (WAD)	Not used	1225 total cm 0.8"i.d.
12-204-00-2-A	Cyanide (Total)	6W, BLB/ 700cm zeus	850cm Tefzel (PN 41300L)
10-107-06-6-D	Kjeldahl Nitrogen**	Not used	850cm Tefzel* (PN 41300L)
10-107-04-3-A	Total Nitrogen	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
10-107-04-3-B	Total Nitrogen	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
10-107-04-3-C	Total Nitrogen	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
10-107-04-3-D	Total Nitrogen	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
10-107-04-3-E	Total Nitrogen	6W germicidal/850 zeus	1200 cm 0.8"i.d. (PN 50028)
10-107-04-3-P	Total Nitrogen	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
12-107-04-3-B	Total Nitrogen	6W germicidal/810 zeus	1200 cm 0.8"i.d. (PN 50028)
12-107-04-3-C	Total Nitrogen	6W germicidal/810 zeus	1200 cm 0.8"i.d. (PN 50028)
31-107-04-3-A	Total Nitrogen	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
31-107-04-3-B	Total Nitrogen	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
10-115-01-3-A	Total Phosphorus	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
10-115-01-3-B	Total Phosphorus	6W germicidal/850 zeus	1200 cm 0.8"i.d. (PN 50028)
10-115-01-3-D	Total Phosphorus	6W germicidal/550 zeus	780 cm 0.8"i.d. (PN 50028)
10-115-01-3-F	Total Phosphorus	6W germicidal/680 zeus	1200 cm 0.8"i.d. (PN 50028)
31-115-01-3-D	Total Phosphorus	6W germicidal/650 zeus	1200 cm 0.8"i.d. (PN 50028)
31-115-01-3-F	Total Phosphorus	6W germicidal/680 zeus	1200 cm 0.8"i.d. (PN 50028)

(*or for the QC8500 ONLY, on an adjacent heater) **Post digestion!

6W BLB Lamp: PN50717 8W BLB Lamp PN 50817 6W 254nm lamp: PN 50775 8W 312 nm lamp: PN9334400 Zeus tubing: PN50728

Comparison tables are available for all methods that are equivalent to NPDES methods

In the list of methods that follows:

Designation in the methods list means the method is EPA accepted as equivalent for NPDWR, NPDES, or both (Check the table above)

^ Designation in the methods list means the method is equivalent for NPDES reporting under the MUR

R Designation means the manifold is compliant with RoHS-2. (Verify with your sales person that the method can be sold into the EU.

Flow Injection Analysis

Ammonia

10-107-06-5-B ^R	0.10 – 1.0 1.0-10.0	0.01	mg N/L as NH ₃	Waters	Gas diffusion method; low-flow method; ISO (11732) 590nm	19-Mar-04
10-107-06-5-J ^{R^A}	0.01-1.0 0.1-20	0.002 0.02	mg N/Las NH ₃	Waters	Gas Diffusion method. Salicylate/DCIC. May be used for TKN as well as brackish/saline samples. 660 nm. Requires a standard heater.	16-Jan-15
10-107-06-6-A ^{R ^}	0.25 – 20	0.13	mg N/L as NH ₃	Waters	Sodium salicylate-based method; 660 nm. inline distillation method; Requires a standard heater and in-line module for distillation step. Samples w/ particulates not suitable. NPDES Equivalent (350.1);	24-Jul-08
10-107-06-6-B ^{R^A}	0.25 – 10	0.066	mg N/L as NH ₃	Waters	Alkaline phenol-based method; 630 nm. inline distillation method; low-flow method; Requires a standard heater and an in-line module for the distillation. Samples w/ particulates not suitable. NPDES Equivalent (350.1);	29-Jul-08
10-107-06-6-E ^{R^A}	10-250	5 (pres.) 1 (un-pres.)	µg N/L as NH ₃	Waters	Alkaline phenol-based method; 630 nm inline distillation method; Requires a standard heater and an in-line module for the distillation. Low-flow method; samples w/ particulates not suitable NPDES Equivalent (350.1);	15-Apr-11

Carbon (Total Dissolved)

12-140-39-5-A ^R	5-400	0.7	mg C/L	Soil Extracts	0.5M K2SO4 extracts. Phenol Red Method. 440 nm. Can measure TN from the same digest with method 12-107-04-3-C. Requires an in-line digestion module (One supplies both channels)	19-Dec-11
----------------------------	-------	-----	--------	---------------	--	-----------

Cyanide

10-204-00-2-C	2 – 100	0.21	µg CN ⁻ /L	Waters	Total Cyanide ; inline method; low-flow method; NPDES Equivalent; Pyridine/barbituric acid, 570 nm. Samples w/ particulates not suitable Inline module and Standard heater required.	14-Sep-07
10-204-00-2-D	5 – 500	0.51	µg CN ⁻ /L	Waters	Total Cyanide ; inline method; low-flow method; Pyridine/barbituric acid, 570 nm. Samples w/ particulates not suitable Inline module and Standard heater required.	19-Sep-07
10-204-00-2-E ^R	2 – 100	0.5	µg CN ⁻ /L	Waters	Total Cyanide ; inline method; low-flow method; lower recovery of ferricyanide; Pyridine/barbituric acid, 570 nm. Samples w/ particulates not suitable Inline module, Standard heater required	3-Dec-08
10-204-00-2-G ^R	0.002 – 0.01 0.1 – 5.0	0.00016 0.015	mg CN ⁻ /L	Waters	Total Cyanide ; inline method; pyridine-free reagents; can be used w/ 10-204-00-2-H for free cyanide; 600 nm multi-range method; Samples w/ particulates not suitable. In-line module and standard heater required.	22-Jun-07
10-204-00-4-B ^R	2.0 – 100	0.16	µg CN ⁻ /L	Waters	WAD Cyanide ; inline method; pyridine-free reagents; Isonicotinic/barbituric acid 600 nm. samples w/ particulates not suitable; MANIFOLD ONLY	27-Jul-07
10-204-00-4B52 ^R					DEDICATED 220V CHANNEL FOR QC8500 Isonicotinic/barbituric acid 600 nm. samples w/ particulates not suitable Requires 2 channels, two heaters, one detector and one valve)	
10-204-00-4-C ^R	2.0 – 100	0.17	µg CN ⁻ /L	Waters	WAD Cyanide ; Low Flow Method. In line distillation method; pyridine-free reagents; Isonicotinic/barbituric acid 600 nm. Samples w/ particulates not suitable; Requires an inline module and a standard heater.	11-Mar-13

Kjeldahl Nitrogen (TKN)

10-107-06-2-S ^R	0.2-20	0.01	mg N/L	Waters	Simplified TKN (s-TKN™) . 520 nm, cadmium reduction . <u>Two channel method</u> TN and NO ₂ + NO ₃ . S-TKN by subtraction. Requires an in-line module	14-Jul-10
10-107-06-5-J ^A	0.1-5.0 0.25-20	0.02 0.05	mg N/L	Waters	Kjeldahl Digests, Salicylate/DCIC 660 nm. copper catalyst . Gas diffusion method. Sea/brackish water. Can also be used for Ammonia. NPDES Equivalent	26-Sept-12
10-107-06-6-D ^A	0.5 – 20	0.25	mg N/L	Waters	Kjeldahl digests; copper catalyst ; inline distillation method; NPDES Equivalent (351.2); samples w/ particulates not suitable. 660 nm. <u>Can be used with brackish/seawater digests.</u> Requires an in-line module and a standard heater or two heated channels (with one heater non-standard).	31-Jul-09

Nitrogen - Total Nitrogen

10-107-04-3-A ^{R*}	200 – 2000	5.6	µg N/L	Waters	Sulfanilamide/NED Cadmium Reduction;540 nm. Total N; inline method; alkaline persulfate digestion; samples w/ particulates not suitable. In-line sample prep module required. Nitrate/Nitrite support added.	16-Nov-09
10-107-04-3-B ^{R*}	0.5 – 30.0	0.1	mg N/L	Waters	Sulfanilamide/NED <u>imidazole buffer</u> ; Cadmium Reduction; 540 nm. Total N; inline method; alkaline persulfate digestion method; samples w/ particulates not suitable. In-line sample prep module required. Nitrate/Nitrite support added.	16-Nov-09

10-107-04-3-C ^R	0.5 – 10.0	0.011	mg N/L	Waters	Sulfanilamide/NED Cadmium Reduction; 540 nm. Total N; inline method; alkaline persulfate digestion method; samples w/ particulates not suitable In-line sample prep module required.	29-Jun-07
10-107-04-3-D ^R	0.05 – 5.0	0.003	mg N/L	Waters	Sulfanilamide/NED Cadmium Reduction; 540 nm. Total N; inline method; alkaline persulfate digestion; samples w/ particulates not suitable.	2-Dec-12
	0.2 – 20.0	0.008			In-line sample prep module required. Nitrate/Nitrite support added.	
10-107-04-3-E ^R	0.05 – 10	0.005	mg N/L	Waters	Sulfanilamide/NED; Cadmium Reduction; 540 nm. Total N; inline method; alkaline persulfate digestion; samples w/ particulates not suitable In-line sample prep module required.	12-Nov-10
10-107-04-3-P ^R	0.2 – 10.0	0.05	mg N/L	Waters	Sulfanilamide/NED; Cadmium Reduction; 540 nm. Total N; inline method; alkaline persulfate digestion; follows Standard Methods (4500-N-B); samples w/ particulates not suitable.	29-Jun-07
12-107-04-3-B	0.2 – 30.0	0.04	mg N/L	Soil extracts	Sulfanilamide/NED Cadmium reduction. 540 nm. Total N; alkaline persulfate digestion; 0.5M K₂SO₄ extracts of soils; inline module required ; samples w/ particulates not suitable.	13-Nov-09
12-107-04-3-C	0.375-30	0.05	mg N/L	Soil extracts	Sulfanilamide/NED Cadmium reduction ,. 540 nm. Total N; 0.5M K₂SO₄ extracts of soils; inline module required ; persulfate digestion; samples w/ particulates not suitable. <u>Dissolved Organic Carbon may be measured in the same digest using 12-140-39-5-A</u>	3-Feb-10

31-107-04-3-A ^R	25 – 1000 1.79-71.43	4.90	µg N/L µM N/L	Brackish / Seawaters	Sulfanilamide/NED Cadmium reduction , 540 nm. Total N; alkaline persulfate inline digestion method; samples w/ particulates not suitable.
----------------------------	--------------------------------	------	-------------------------	-------------------------	--

Phenol

10-210-00-3-A ^R	2 – 200	0.28	µg phenol/L	Waters	Volatile phenol; 4-amino antipyrene method; 500 nm. inline method; samples w/ particulates not suitable; This PN manifold only QC8500 115V dedicated channel QC8500 220V dedicated channel	20-Dec-06
10-210-00-3A51 ^R 10-210-00-3A52 ^R						
10-210-00-3-C* [^]	2 – 200	0.61	µg phenol/L	Waters	Volatile phenol; 4-aminoantipyrene method; 500 nm. inline method; NPDES Equivalent (420.4); samples w/ particulates not suitable. This PN for the manifold only QC8500 115V dedicated channel QC8500 220V dedicated channel	15-Oct-08
10-210-00-3C51* [^] 10-210-00-3C52* [^]						

Phosphorus, Total (Acidic Persulfate)

10-115-01-3-A ^{R^}	0.1 – 10.0	0.007	mg P/L	Waters	Total P; acidic persulfate digests; molybdate method, 880 follows Standard Methods (4500-P-I); samples w/ particulates not suitable. Requires an in-line sample prep module. Can also use for orthophosphorus over the same range.	18-Nov-09
10-115-01-3-B ^{R*^}	0.1 – 4.0	0.01	mg P/L	Waters	Total P; acidic persulfate digests; molybdate method; 880 nm; samples w/ particulates not suitable Requires an in-line sample prep module. Can also use for orthophosphorus over the same range.	18-Nov-09
10-115-01-3-C ^{R*^}	0.05 – 1.0	0.0011	mg P/L	Waters	Total P; acidic persulfate digests; molybdate method; 880 nm. Samples w/ particulates not suitable Requires an in-line sample prep module. Can also use for orthophosphorus over the same range.	18-Nov-09

10-115-01-3-E ^{RA}	10 – 500	1.4	µg P/L	Waters	Total P; acidic persulfate digests; molybdate method; 880 nm. Samples w/ particulates not suitable Requires an in-line sample prep module and standard heater. Can also use for orthophosphorus over the same range.	5-Jul-07
10-115-01-3-F ^{RA}	2 – 100	0.42	µg P/L	Waters	Total P; acidic persulfate digests; molybdate method; 880 nm method; samples w/ particulates not suitable Requires an in-line sample prep module and non-standard heater.	13-Nov-06
31-115-01-3-D ^R	0.050 – 1.0 1.63-32.36	0.002	mg P/L µM P/L	Brackish / Seawaters	Total P; molybdate method; 880 nm. inline persulfate digestion; samples w/ particulates not suitable. glass calibration vials. Requires an in-line sample prep module.	5-Jul-07
31-115-01-3-F ^R	2-100 0.065-3.23	0.59	µg P/L µM P/L	Brackish / Seawaters	Total P; molybdate based method; 880nm, inline persulfate digestion; samples w/ particulates not suitable Requires glass standard and sample vials an in-line sample prep module, and non-standard heater.	13-Oct-08

Sulfide

10-116-29-3-A ^R	0.01 – 2.0	0.006	mg S/L	Waters	In line distillation method; 660 nm. Requires two dedicated channels with one standard and one non-standard heater; samples w/ particulates not suitable manifold only	4-Oct-07
10-116-29-3A51 ^R					Dedicated channels; QC8500 115V	
10-116-29-3A52 ^R					Dedicated channels; QC8500 220V	
10-116-29-3B52 ^R	1.0 – 10.0	0.2	mg S/L	Waters	In line distillation method; 660 nm. requires two dedicated channels; Requires two dedicated channels with one standard and one non-standard heater samples w/ particulates not suitable; Dedicated channels; QC8500 220V	5-Jul-07


Sulfite

21-116-11-2-D ^R	0.5 – 30.0 25 -300.0	0.15 1.5	mg SO ₃ ²⁻ /L mg SO ₃ ²⁻ /L	Beer,Wine	Pararosaniline method; 560 nm Requires a standard heater	06-Oct-16
21-116-11-2-E ^R	0.5-20.0	0.05	mg SO ₃ ²⁻ /L	Beverages	Pararosaniline method 560 nm. Determination in coconut water. Requires a standard heater.	30-Jan-15


Surfactants (MBAS)

10-306-00-1-C ^R	0.025 – 2.0 0.010 – 1.0	0.004 0.0056	mg/L as LAS mg/L as SDS	Waters	Methylene blue method; 650 nm. dual extraction method. SDS or LAS. Glass calibration and standard vials must be used.	19-Dec-08
10-306-00-1-D ^{R^A}	0.010 – 1.0	0.0024	mg SDS/L	Waters	Methylene blue method; 650 nm.single extraction method; NPDES Equivalent; follows Standard Methods (5540-C). SDS. Glass calibration and standard vials must be used.	25-Mar-08
10-306-00-1-E ^R	0.1 – 20.0	0.05	mg SDS/L	Waters	Methylene blue method; 650 nm. dual extraction method (SDS only). Glass calibration and standards vials must be used.	29-Sep-05
10-306-00-1-F	0.06-2.4	0.009	mg LAS/L	Waters	Methylene blue method; 650 nm. dual extraction method (LAS only). Glass calibration and standards vials must be used. SM5540C/ASTM2330-02	9-Apr-14


Selected Parts. Please note this is NOT a complete listing




50902 or 50902L
Tee Fitting




50045L
Reagent Line Weight




50046L
50036



50065 Flow cell flare, 13 cm




50067 Flow cell flare, 35 cm




31933 Flow cell, 10mm


Nipples




50015 Large




50015A Small




50014L Reducing




50019 Stainless Steel




Interference Filters
89### ### indicates the wavelength
400 nm



24927 Switching Valve
50962 Flare Kit for Nitrate
50963 Flare Kit for Sulfate
24012 Connector, Amber
24015 Washer




50237A Cadmium Column




50254 Cadmium Column Maintenance kit

Union Fitting




50913




50901

Small Pump Tube Adapter




50906




85287 for QC8500

Large Pump Tube Adapter




85293 for QC8500




50907


Pump Tube Adapter




50905




85269 for QC8500




50009
Tube Connector, Amber




31077 for QC8000 Valve
Valve Connector, White




85245
One-piece fitting
for QC8500 valve




28051L
Rheodyne Nut & Ferrule




28058 Rheodyne Ferrule




28057
Upchurch Nut & Ferrule




50013
O-ring, Viton




50024L
Large Clip, Gray




50023L
Small Clip, Gray




50913 Teflon
Union Fitting
(ultem nuts)




50007 or 85258 for QC8500
Large Collar, Black




50006 or 85257 for QC8500
Small Collar, White



50060
O-ring Probe



50100L
Degassing Tube



50031
Scalpel



28081
Blades for PEEK
Tubing Cutter




50021
Transmission Tubing, PVC, 0.030" ID



50029
Transmission Tubing, PVC, 0.060" ID




50091 Microloop Tubing, Teflon, 16 cm
50092 Microloop Tubing, Teflon, 12.5 cm



28080 PEEK Tubing Cutter



Teflon Tubing
50927 Manifold, 0.022" ID, Green
50928 Manifold, 0.032" ID
41300L Tefzel Tubing, 0.040" ID x 0.060" OD
30928 Zeus Tubing, 0.032" ID



Coils

Size	Coil Support Only	Wrapped, Teflon 0.022" ID	Wrapped, Teflon 0.032" ID
1" or 4.5 cm	50016L	50981	50916
2" or 7cm	50018L	50982	50918
2.5" or 8 cm	50017	50983	50917
4" or 12cm	50020	50984	50920
8" or 22 cm	50022L	N/A	50922



Alternating Coils 12 cm

Coil Support Only	Wrapped, Teflon 0.022" ID	Wrapped, Teflon 0.032" ID
50039	50985	50921



Pump Tubing

534XX	PVC
544XX	Duraprene
494XX	Silicone
654XX	Acidflex

XX is the number that specifies the pump tube color:

05	Orange-Yellow, 0.020" ID	13	Blue-Blue, 0.065" ID
06	Orange-White, 0.025" ID	14	Green-Green, 0.073" ID
07	Black-Black, 0.030" ID	15	Purple-Purple, 0.081" ID
08	Orange-Orange, 0.035" ID	16	Purple-Black, 0.090" ID
09	White-White, 0.040" ID	17	Purple-Orange, 0.100" ID
10	Red-Red, 0.045" ID	18	Purple-White, 0.110" ID
11	Gray-Gray, 0.051" ID	19	Yellow-Blue, 0.060" ID
12	Yellow-Yellow, 0.056" ID		



Bottles

28193L	Glass Bottle, 1000 mL
35102	Glass Bottle, 2000 mL
43915	Glass Bottle, 100 mL (includes cap)



50012
O-ring Remover



28101L
End Plug



50062
Drill Bit

Lachat Instruments Brand Loveland, Colorado USA

United States: 800-247-7613 tel | 970-669-3050 tel | 970-461-3915 fax | sales@lachatinstruments.com | techhelp@hach.com
Outside United States: Contact the Lachat office or distributor serving you. | sales@lachatinstruments.com | intltech@hach.com
www.lachatinstruments.com

Printed in U.S.A. ©Hach Company, 2016. All rights reserved.

LACHAT
INSTRUMENTS

