



**52 61**

**DOC022.98.90550**

**09-2015, Edition 1**

# 52 61

2

# Puntos clave. Punti chiave.

## Points clé. Key points.

1

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*Protective cap.*

2

Cabezal con conector S7.

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*S7 connector head.*

3

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*Orifizio riempimento.*

Orifice de remplissage.

*Filling port cap.*

4

Electrolito de referencia.

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*Reference electrolyte.*

5

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*Storage protector.*

6

Elemento de referencia.

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*Reference element.*

7

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*Diaframma ceramico.*

Diaphragme en céramique.

*Ceramic diaphragm.*

8

Anillo de platino.

*Anello de platino.*

Anneau en platine.

*Platinum ring.*

# Electrodo de Pt 52 61

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# Pt electrode 52 61

## Introduction

The **52 61** is a platinum electrode for oxide-reduction potential measurements (ORP) and redox titrations.

## Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it **⑤**.

Remove the cap from the filling port **③**.

## Measurements. Recommendations

- The sample must cover the electrode's diaphragms **⑦**.
- After measurement, clean the electrode with distilled water.

## Maintenance

- Periodically refill the electrode with electrolyte. Maintain the level as high as possible **④**.
- Diaphragm cleaning: use the diaphragm cleaning solution.
- Platinum ring cleaning: can be ribbed with soft abrasive, for example tooth paste.

## Troubleshooting

... the reading is always the same.  
Short circuit in the electrode or connector.

*Substitute the electrode.*

... unstable reading.  
Broken cable.

*Substitute the cable.*

... very slow response.  
Dirt on the diaphragm or the silver ring.  
Reference system worn out.

*Clean, regenerate or substitute the electrode.*

## Specifications

Measuring range	± 2000 mV
Indicator metal	platinum
Operating temp.	0 ... 80 °C
Reference	encapsulated Ag / AgCl crystals
Diaphragm	ceramic
Electrolyte	KCl 3M
Body material	glass
Min. puncture depth	12 mm 

## Applications

General measurements of redox potential.

Redox titrations.

## Limits

Products with colloids or solids in suspension.

## Warranty

HACH guarantees this electrode against manufacturing defects.

**Validity:** 6 months from shipment date.

## Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

## QUALITY CERTIFICATE

The attached electrode 52 61 has passed the quality test and follows the next specifications:

- Deviation  $< \pm 10$  mV in redox standard.
- Response time in redox standard  $< 20$  s.

**HACH COMPANY** World Headquarters

P.O. Box 389, Loveland, CO 80539-0389 U.S.A.

Tel. (970) 669-3050

(800) 227-4224 (U.S.A. only)

Fax (970) 669-2932

[orders@hach.com](mailto:orders@hach.com)

[www.hach.com](http://www.hach.com)

**HACH LANGE GMBH**

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D-40549 Düsseldorf, Germany

Tel. +49 (0) 2 11 52 88-320

Fax +49 (0) 2 11 52 88-210

[info@de.hach.com](mailto:info@de.hach.com)

[www.de.hach.com](http://www.de.hach.com)

**HACH LANGE Sàrl**

6, route de Compois

1222 Vésenaz

SWITZERLAND

Tel. +41 22 594 6400

Fax +41 22 594 6499

---

© Hach Company/Hach Lange GmbH, 2015.

All rights reserved. Printed in Germany..