


### Step 1 - UltraCheck™ 1 Calibration

Perform one UltraCheck 1 calibration per day or per each set of samples analyzed.




**NOTE:** If  $RLU_{ATP1} \leq 50,000$  using a PhotonMaster or Lumitester C-110, rehydrate a new bottle of Luminase<sup>XL</sup> for maximum sensitivity.

### Step 2 – Total ATP (tATP™) Analysis (1 per sample)

Included in QG21S™ and QG21St test kits.

#### 2.1 – PRE-DILUTION

Pre-dilute sample.



#### 2.4 – ASSAY


Measure ATP concentration.



**NOTE:** If  $RLU_{tATP} \leq 10$  using a PhotonMaster or Lumitester C-110, you are below the low- detection limit.

#### 2.2 – EXTRACTION

Add pre-diluted sample to extract ATP.




**NOTE:** If  $RLU_{tATP} \leq 50$  using a PhotonMaster or Lumitester C-110, consider accounting for background ( $RLU_{bg}$ ). See Test Kit Instructions for guidance.

**Total ATP (tATP) Calculation:**

$$tATP (pg \text{ ATP/mL}) = \frac{RLU_{tATP}}{RLU_{ATP1}} \times 100,000 (pg \text{ ATP / mL})$$

#### 2.3 – DILUTION

Dilute out interferences.




### Step 3 – Dissolved ATP (dATP™) Analysis (1 per sample)

Included in QG21S test kit only.


#### 3.1 – DILUTION

Add pre-diluted sample to recover ATP.



#### 3.2 – ASSAY

Measure ATP concentration.



**NOTE:** If  $RLU_{dATP} \leq 10$  using a PhotonMaster or Lumitester C-110, you are below the low-detection limit.

**NOTE:** If  $RLU_{dATP} \leq 50$  using a PhotonMaster or Lumitester C-110, consider accounting for background ( $RLU_{bg}$ ). See Test Kit Instructions for guidance.

**Dissolved ATP (dATP) Calculation:**

$$dATP (pg \text{ ATP/mL}) = \frac{RLU_{dATP}}{RLU_{ATP1}} \times 100,000 (pg \text{ ATP / mL})$$

### Calculations:

**NOTE:** When using the QG21S tATP – only kit, skip final calculations and interpret tATP results as you would cATP results using the Interpretation Guidelines.

**NOTE:** If the results show for a given sample that  $dATP (ng/mL) > tATP (ng/mL)$ , report  $cATP^{TM} = 0$  and  $BSI^{TM} = 100\%$

**Cellular ATP (cATP) Calculation:**

$$cATP \left( \frac{pg \text{ ATP}}{mL} \right) = tATP (pg \text{ ATP / mL}) - dATP (pg \text{ ATP / mL})$$

**Microbial Equivalent (ME/mL) Calculation:**

$$cATP \left( \frac{ME}{mL} \right) = cATP (pg \text{ ATP / mL}) \times \frac{1 \text{ ME}}{0.001 \text{ pg ATP}}$$

**Biomass Stress Index (BSI) Calculation:**

$$BSI (\%) = \frac{dATP (pg \text{ ATP / mL})}{tATP (pg \text{ ATP / mL})} \times 100\%$$

**NOTE:** 1 ME (Microbial Equivalent) assumes 0.001 pg (1 fg) ATP per cell.

### Interpretations Guidelines

Application	Good Control (pg cATP/mL)	Preventative Action (pg cATP/mL)	Corrective Action (pg cATP/mL)
<b>Product Quality Control</b> (Paint, Coating, Slurries)	< 100	100 to 1,000	> 1,000

For BSI (when applicable), it can generally be interpreted that good control is achieved at levels of 75% or above. Preventive action should be taken at levels between 50% and 75%, and corrective action should be taken at levels below 50%.

**NOTE:** Interpretation Guidelines provided for general guidance. For best results, establish your own baseline and control levels.