



STABLCAL[®] STABILIZED FORMAZIN STANDARDS*
CERTIFICATE OF COMPLIANCE
StablCal[®] Calibration Set # 2659105

TEST PERFORMED: The Standards listed below were measured using a Hach EPA Turbidimeter. This instrument was calibrated on fresh formazin standards that are traceable to the Stock Primary Standard Formazin Suspension as specified in Standard Methods, Section 2130 B, part 3b.

	Cat. No.	Value	Variation
Standard #1	2659801	1.0 NTU	+/- .05 NTU
Standard #2	2659901	10.0 NTU	+/- 0.5 NTU
Standard #3	2660201	100 NTU	+/- 5.0 NTU
Standard #4	2660601	1000 NTU	+/- 50 NTU

*StablCal[®] Stabilized Formazin Standards are cited as primary calibration standards in Hach Method 8195, an acceptable version of USEPA Method 180.1.

Download Certificates of Analysis by entering the StablCal Calibration Set part number and lot number on our web site

Signed: _____

on behalf of Hach Company

Hach Cat. No 2659187