



Ready-to-use Reagents for EZ Series Online Analyzers



Are you tired of mixing your own reagents?

How often do you have to prepare the reagents for your EZ Series analyzer and experienced the situation where you have other things to do instead of spending your time preparing the chemistries for your analyzer?

And do you fear that the preparation of the reagents could be a potential source of error where you lose confidence in results?

Would you prefer to use your online instruments with ready-to-use reagents from the same supplier?

With Hach's EZ reagent kits you can trust your results!

Ready-to-use

Hach EZ reagents are ready-to-use with EZ Series Online Analyzers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Product Table

Parameter	Prod. No.	Description	
Total Hardness	APPC1016KTO	Reagent Kit for EZ1016 Total Hardness Analyzer	
	APPC1016-01	EZ1016 Total Hardness Acid Solution	
	APPC1016-02	EZ1016 Total Hardness Buffer Solution	
	APPC1016-03	EZ1016 Total Hardness Color Solution	
	APPC1016-04	EZ1016 Total Hardness EDTA Solution	
	12153	Calcium Standard Solution, 1000 mg/L	
Aluminium	APPC1001KTO	Reagent Kit for EZ1001 Aluminium Analyzer	
	APPC1001-01	EZ1001 Aluminium Buffer Solution	
	APPC1001-02A	EZ1001 Aluminium Color Indicator	
	APPC1001-03A	EZ1001 Aluminium Mixed Reagent	
	APPC1001-03B	EZ1001 Aluminium Ascorbic Acid	
	27253	Deionized Water (1000 mL)	
	1417442	Aluminium Standard Solution, 100 mg/L	
Total Aluminium	APPC2000KTO	Reagent Kit for EZ2000 Total Aluminium Analyzer	
	APPC1001-01	EZ1001 Aluminium Buffer Solution	
	APPC1001-02A	EZ1001 Aluminium Color Indicator	
	APPC1001-03A	EZ1001 Aluminium Mixed Reagent	
	APPC2000-04	EZ2000 Total Aluminium Nitric Acid Solution	
	27253	Deionized Water (1000 mL)	
	1417442	Aluminium Standard Solution, 100 mg/L	
	2305842	Aluminium Standard Solution, 10 mg/L, 100 mL	
Boron	APPC1004KTO	Reagent Kit for EZ1004 Boron Analyzer	
	APPC100401	EZ1004 Boron Buffer Solution	
	APPC100402	EZ1004 Boron Glycerol Solution	
	APPC100402A	EZ1004 Boron Color Indicator	
Calcium		2305442	Calcium Standard Solution, 10 mg/L
	Hydrogen Peroxide	APPC1022KTO	Reagent Kit for EZ1022 H₂O₂ Analyzer
APPC1022-01		EZ1022 H ₂ O ₂ Buffer Solution	
APPC1022-02		EZ1022 H ₂ O ₂ Color Solution	
APPC1022-DPD		EZ1022 H ₂ O ₂ DPD Powder	
Iron	APPC1024KTO	Reagent Kit for EZ1024 Iron Analyzer	
	APPC1024-01	EZ1024 Iron Buffer Solution	
	APPC1024-02	EZ1024 Iron Color TPTZ Solution	
	APPC1024-03	EZ1024 Iron Reducing Reagent	
	14049	Iron Standard Solution, 10 mg/L	
Manganese	APPC1025KTO	Reagent Kit for EZ1025 Manganese Analyzer	
	APPC1025-01	EZ1025 Manganese Buffer Solution	
	APPC1025-02	EZ1025 Manganese Color Solution	
	APPC1025-03	EZ1025 Manganese Reducing Reagent	
	APPC1025-04	EZ1025 Manganese EDTA Solution	
	1279142	Manganese Standard Solution, 1000 mg/L	
Total Manganese	APPC2003KTO	Reagent Kit for EZ2003 Total Manganese Analyzer	
	1279142	Manganese Standard Solution, 1000 mg/L	
	APPC1025-01	EZ1025 Manganese Buffer Solution	
	APPC1025-02	EZ1025 Manganese Color Solution	
	APPC1025-03	EZ1025 Manganese Reducing Reagent	
	APPC1025-04	EZ1025 Manganese EDTA Solution	
	APPC2003-05	EZ2003 Total Manganese Nitric Solution	
Phosphate	APPC1031KTO	Reagent Kit for EZ1031 Phosphate Analyzer	
	APPC1031-01	EZ1031 Phosphate Color Solution	
	17149	Phosphate Standard Solution, 50 mg/L	
Total Hardness	APPC4041KTOL	Reagent Kit for EZ4041 Total Hardness Analyzer (10-250 mg/L)	
	APPC4041-01	EZ4041 Total Hardness Acid Solution	
	APPC4041-02	EZ4041 Total Hardness Buffer Solution	
	APPC4041-03	EZ4041 Total Hardness Color Solution	
	APPC4041-04	EZ4041 Total Hardness EDTA Solution, 0.005 M	
	12153	Calcium Standard Solution, 1000 mg/L	

Product Table

Parameter	Prod. No.	Description
Total Hardness	APPC4041KTOM	Reagent Kit for EZ4041 Total Hardness Analyzer (50-500 mg/L)
	APPC4041-01	EZ4041 Total Hardness Acid Solution
	APPC4041-02	EZ4041 Total Hardness Buffer Solution
	APPC4041-03	EZ4041 Total Hardness Color Solution
	APPC4041-05	EZ4041 Total Hardness EDTA Solution, 0.01 M
	12153	Calcium Standard Solution, 1000 mg/L
Total Hardness	APPC4041KTOS	Reagent Kit for EZ4041 Total Hardness Analyzer (100-1000 mg/L)
	APPC4041-01	EZ4041 Total Hardness Acid Solution
	APPC4041-02	EZ4041 Total Hardness Buffer Solution
	APPC4041-03	EZ4041 Total Hardness Color Solution
	APPC4041-06	EZ4041 Total Hardness EDTA Solution, 0.02 M
	12153	Calcium Standard Solution, 1000 mg/L
Total Hardness	APPC4043KTO	Reagent Kit for EZ4043 Total Hardness Analyzer (0.25-10 mg/L)
	APPC4041-01	EZ4041 Total Hardness Acid Solution
	APPC4041-02	EZ4041 Total Hardness Buffer Solution
	APPC4041-03	EZ4041 Total Hardness Color Solution
	APPC4043-04	EZ4043 Total Hardness EDTA Solution, 0.005 N
	4457649	Calcium Standard Solution, 100 mg/L
Fluoride	APPC3007KTO	Reagent Kit for EZ3007 Fluoride Analyzer
	APPC3007-01	EZ3007 Fluoride Buffer Solution (TISAB)
	35949	Fluoride Standard Solution, 10 mg/L
Fluoride	APPC3507KTO	Reagent Kit for EZ3507 Fluoride Analyzer
	APPC3507-01	EZ3507 Fluoride Buffer Solution (TISAB)
	23249	Fluoride Standard Solution, 100 mg/L
Total Alkalinity <i>The Alkalinity titrants can be used in other EZ Series Alkalinity Titration Analyzers.</i>	APPC4004KTOL	Reagent Kit for EZ4003 & EZ4004 Total Alkalinity Analyzer (10-500 mg/L)
	APPC4004-01	EZ4003/EZ4004 Total Alkalinity Acid Solution, 0.025 N HCl
	2826253	Alkalinity Standard Solution, 500 mg/L
Total Alkalinity	APPC4004-02	EZ4004 Total Alkalinity Acid Solution, 0.05 N HCl (25-1250 mg/L)
Total Alkalinity	APPC4004-03	EZ4004 Total Alkalinity Acid Solution, 0.1 N HCl (50-2500 mg/L)
Total Alkalinity	APPC4004-04	EZ4004 Total Alkalinity Acid Solution, 0.2 N HCl (100-5000 mg/L)
Ammonia	APPC1102KTO	Reagent Kit for EZ1102 Ammonia Analyzer
	APPC1102-01	EZ1102 Ammonia Chlorine Solution
	APPC1102-02	EZ1102 Ammonia Buffer Solution
	APPC1102-03	EZ1102 Ammonia Indicator Solution
	2354153	Ammonia Standard Solution, 1000 mg/L
Sulfate	APPC1036KTO	Reagent Kit for EZ1036 Sulfate Analyser
	APPC1036-01	EZ1036 Sulfate Buffer Solution
	APPC1036-02	EZ1036 Sulfate Barium Chloride Solution
	APPC1036-03	EZ1036 Sulfate EDTA Solution
	1425249	Sulfate Standard Solution, 2500 mg/L
Total Nitrogen, Total Phosphorus	APPC76NP-01NA	EZ7600 TN/TP Series Color Solution
Total Nitrogen, Total Phosphorus	APPC76NP-02NA	EZ7600 TN/TP Series Persulfate Solution
Total Nitrogen, Total Phosphorus	APPC76NP-03NA	EZ7600 TN/TP Series NaOH Solution
Total Nitrogen, Total Phosphorus	APPC76NP-04NA	EZ7600 TN/TP Series HCl Solution
Total Nitrogen, Total Phosphorus	APPC76NP-05NA	EZ7600 TN/TP Series Reducing Reagent
Total Nitrogen, Total Phosphorus	APPC76NP-06NA	EZ7600 TN/TP Series Reference 1 Solution
Total Nitrogen, Total Phosphorus	APPC76xx-07NA*	EZ76xx TN/TP Analyzer Reference 2 Solution

*Article number depending on EZ analyzer model, 76xx = 7621/7632/7641/7642/7652/7653/7654/7663/7664/7665/7675/7676/7685/7686



World Headquarters: Loveland, Colorado USA | hach.com

United States 800-227-4224 fax: 970-669-2932 email: orders@hach.com
 Outside United States 970-669-3050 fax: 970-461-3939 email: int@hach.com

©Hach Company, 2023. All rights reserved.

In the interest of improving and updating its equipment, Hach Company reserves the right to alter specifications to equipment at any time.

DOC053.53.35315.Feb24