



TL23 Series

Laboratory Turbidimeters

The TL23 is an essential part of your quality control and compliance program. Let Hach® Service help ensure your instrument operates at peak performance to deliver the accurate turbidity measurements you require.

Benefits of Service

- Ensures critical turbidity measurements are accurate and reliable
- Detects and prevents issues that can degrade performance or reduce instrument life
- Provides documentation to support your quality management and compliance programs
- Loaner instruments ensure continuous operations during service
- Technical support hotline provides rapid and direct resolution to many issues

Annual Average Maintenance Time: Up to 1.5 hours
Recommended Maintenance Frequency: 1x / year

Essential Maintenance Tasks Performed

Routine:

- Clean and visually inspect housing, filter assembly, and vial compartment
- Replace lamp (EPA version only)
- Update firmware to latest version
- Calibrate with certified StablCal standards

Advanced:

- Replace lens (EPA version only)
- Perform hardware check to verify system performance and component functionality
- Produce service report and performance certificate to support quality management and compliance

Parts Replaced During Maintenance (EPA Version Only)

Part #	Description	Qty.	Frequency
9659000	Lamp	1	36 months
4705100	Lens	1	12 months

These parts and replacement labour are included with a Hach Service Plan at the recommended frequency as long as the plan is renewed (EPA version only). Coverage is available for repairs at the Hach Service Centre. Loaner instruments are subject to availability.

For more information, visit: www.uk.hach.com/service