

Ready-to-use Reagents for EZ Series Online Analysers



Are you tired of mixing your own reagents?

How often do you have to prepare the reagents for your EZ Series analyser and experienced the situation where you have other things to do instead of spending your time preparing the chemistries for your analyser?

And do you fear that the preparation of the reagents could be a potential source of error where you lose confidence in results?

Would you prefer to use your online instruments with ready-to-use reagents from the same supplier?

With Hach's EZ reagent kits you can trust your results!

Ready-to-use

Hach EZ reagents are ready-to-use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Product Table

Parameter	Part number	Description
Total Hardness	APPC1016KTO-EU	Reagent Kit for EZ1016 Total Hardness Analyser
	APPC1016-01-EU	EZ1016 Total Hardness Acid Solution
	APPC1016-02-EU	EZ1016 Total Hardness Buffer Solution
	APPC1016-03	EZ1016 Total Hardness Colour Solution
	APPC1016-04	EZ1016 Total Hardness EDTA Solution
Iron <i>The TPTZ Iron reagents can be used in other EZ Series Iron Analysers.</i>	APPC1024KTO-EU	Reagent Kit for EZ1024 Iron Analyser
	APPC1024-01-EU	EZ1024 Iron Buffer Solution
	APPC1024-02-EU	EZ1024 Iron Colour TPTZ Solution
	APPC1024-03-EU	EZ1024 Iron Reducing Reagent
Manganese	APPC1025KTO-EU	Reagent Kit for EZ1025 Manganese Analyser
	APPC1025-01-EU	EZ1025 Manganese Buffer Solution
	APPC1025-02-EU	EZ1025 Manganese Colour Solution
	APPC1025-03-EU	EZ1025 Manganese Reducing Reagent
	APPC1025-04-EU	EZ1025 Manganese EDTA Solution
Total Manganese	APPC2003KTO-EU	Reagent Kit for EZ2003 Total Manganese Analyzer
	APPC1025-01-EU	EZ1025 Manganese Buffer Solution
	APPC1025-02-EU	EZ1025 Manganese Colour Solution
	APPC1025-03-EU	EZ1025 Manganese Reducing Reagent
	APPC1025-04-EU	EZ1025 Manganese EDTA Solution
	APPC2003-05-EU 1279142	EZ2003 Total Manganese Nitric Solution Manganese Standard Solution, 1000 mg/L
Total Hardness <i>The Hardness titrants can be used in other EZ Series Hardness Titration Analysers.</i>	APPC4041KTOL-EU	Reagent Kit for EZ4041 Total Hardness Analyser (10-250 mg/L)
	APPC4041-01-EU	EZ4041 Total Hardness Acid Solution
	APPC4041-02-EU	EZ4041 Total Hardness Buffer Solution
	APPC4041-03	EZ4041 Total Hardness Colour Solution
	APPC4041-04	EZ4041 Total Hardness EDTA Solution, 0.005 M
Total Hardness	APPC4041KTOM-EU	Reagent Kit for EZ4041 Total Hardness Analyser (50-500 mg/L)
	APPC4041-01-EU	EZ4041 Total Hardness Acid Solution
	APPC4041-02-EU	EZ4041 Total Hardness Buffer Solution
	APPC4041-03	EZ4041 Total Hardness Colour Solution
	APPC4041-05	EZ4041 Total Hardness EDTA Solution, 0.01 M
Total Hardness	APPC4041KTOS-EU	Reagent Kit for EZ4041 Total Hardness Analyser (100-1000 mg/L)
	APPC4041-01-EU	EZ4041 Total Hardness Acid Solution
	APPC4041-02-EU	EZ4041 Total Hardness Buffer Solution
	APPC4041-03	EZ4041 Total Hardness Colour Solution
	APPC4041-06	EZ4041 Total Hardness EDTA Solution, 0.02 M
Total Hardness	APPC4043KTO-EU	Reagent Kit for EZ4043 Total Hardness Analyser (0.25-10 mg/L)
	APPC4041-01-EU	EZ4041 Total Hardness Acid Solution
	APPC4041-02-EU	EZ4041 Total Hardness Buffer Solution
	APPC4041-03	EZ4041 Total Hardness Colour Solution
	APPC4043-04	EZ4043 Total Hardness EDTA Solution, 0.005 N
	4457649	Calcium Standard Solution, 100 mg/L
Total Alkalinity <i>The Alkalinity titrants can be used in other EZ Series Alkalinity Titration Analysers.</i>	APPC4004-01-EU	EZ4003/EZ4004 Total Alkalinity Acid Solution, 0.025 N HCl (10-500 mg/L)
	APPC4004-02-EU	EZ4004 Total Alkalinity Acid Solution, 0.05 N HCl (25-1250 mg/L)
Total Alkalinity	APPC4004-03-EU	EZ4004 Total Alkalinity Acid Solution, 0.1 N HCl (50-2500 mg/L)
Total Alkalinity	APPC4004-04-EU	EZ4004 Total Alkalinity Acid Solution, 0.2 N HCl (100-5000 mg/L)

Product Table

Parameter	Part number	Description
Phosphate	APPC1031-01-EU	EZ1031 Phosphate Colour Solution
Ammonia	APPC1102KTO	Reagent Kit for EZ1102 Ammonia Analyser
	APPC1102-01	EZ1102 Ammonia Chlorine Solution
	APPC1102-02	EZ1102 Ammonia Buffer Solution
	APPC1102-03	EZ1102 Ammonia Indicator Solution
	2354153	Ammonium Standard Solution, 1000 mg/L
Sulfate	APPC1036KTO-EU	Reagent Kit for EZ1036 Sulfate Analyser
	APPC1036-01	EZ1036 Sulfate Buffer Solution
	APPC1036-02-EU	EZ1036 Sulfate Barium Chloride Solution
	APPC1036-03	EZ1036 Sulfate EDTA Solution
	1425249	Sulfate Standard Solution, 2500 mg/L
Fluoride	APPC3007-01	EZ3007 Fluoride Buffer Solution (TISAB)
Fluoride	APPC3507-01	EZ3507 Fluoride Buffer Solution (TISAB)
Aluminium	APPC1001KTO	Reagent Kit for EZ1001 Aluminium Analyser
	APPC1001-01	EZ1001 Aluminium Buffer Solution
	APPC1001-02A	EZ1001 Aluminium Colour Indicator
	APPC1001-03A	EZ1001 Aluminium Mixed Reagent
	APPC1001-03B	EZ1001 Aluminium Ascorbic Acid
	27253	Deionized Water
	1417442	Aluminium Standard Solution, 100 mg/L
Total Aluminium	APPC2000KTO-EU	Reagent Kit for EZ2000 Total Aluminium Analyser
	APPC1001-01	EZ1001 Aluminium Buffer Solution
	APPC1001-02A	EZ1001 Aluminium Colour Indicator
	APPC1001-03A	EZ1001 Aluminium Mixed Reagent
	APPC1001-03B	EZ1001 Aluminium Ascorbic Acid
	APPC2000-04-EU	EZ2000 Total Aluminium Nitric Acid Solution
	27253	Deionized Water
	1417442	Aluminium Standard Solution, 100 mg/L
	2305842	Aluminium Standard Solution, 10 mg/L
Boron	APPC1004KTO-EU	Reagent Kit for EZ1004 Boron Analyser
	APPC100401-EU	EZ1004 Boron Buffer Solution
	APPC100402	EZ1004 Boron Glycerol Solution
	APPC100402A	EZ1004 Boron Colour Indicator