

# DO6 ⇔ 5120 / 5130 sensor

## COMPATIBILITY

COMPATIBILITÉ / KOMPATIBILITÄT / COMPATIBILIDAD / COMPATIBILITÀ

Instruments	Instrument label **	DO sensor 5120	DO sensor 5130
sensION+ DO6	-	✓	✘
	"LPV...7772" - ?	*	✓
sensION+ DO6 DL	-	*	✓

Kits	Article Number	DO sensor 5120	DO sensor 5130
sensION+ DO6	LPV4551...	✓	✘
	LPV4552...	*	✓
sensION+ DO6 DL	LPV4552DL...	*	✓

\*\* **sensION™+ DO6**  
 LPV4500.97.7772  
 S/N: 221021  
 Batteries 3x AA/LR6, 1.5V IP67  
CRISON INSTRUMENTS, S.A.  
 E-08328 ALELLA - Barcelona  
 Developed and designed in Spain  
 Assembled in China



✓	Compatible
✘	Not compatible
*	ADJUST meter settings



- \* EN *How to adjust DO6 instrument settings?*
- FR *Comment régler les paramètres de l'instrument DO6?*
- DE *Anpassen der DO6 Messgeräteeinstellungen*
- ES *Ajuste de los parámetros del instrumento DO6*
- IT *Come regolare le impostazioni dello strumento DO6*

